NADAL[®] H. pylori Antigen Test (test cassette)

REF 262002

DE Gebrauchsanweisung	2	PL Sposób użycia	17
EN Instruction for use	5	Symbols	20
FR Instructions d'utilisation	8	Our Teams	20
ES Instrucciones de uso	11		
IT Istruzioni per l'uso	14		

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Carl-Zeiss-Strasse 12
47445 Moers
Germany

Tel: +49 (2841) 99820-0
Fax: +49 (2841) 99820-1

Friedenstrasse 32
93053 Regensburg
Germany

Tel: +49 941 29010-0
Fax: +49 941 29010-50

www.nal-vonminden.com
info@nal-vonminden.com

Directors:
Sandra von Minden
Roland Meißner
Thomas Zander

Commercial reg. Kleve
HRB 5679
Steuer-Nr. 244/133/00130
UST-ID-Nr. DE 189 016 086

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2. Soll, AH. Pathogenesis of peptic ulcer and implications for therapy. New England J. Med. (1990), 322: 909-16.
3. Hazell, SL, et al. Campylobacter pyloridis and gastritis I: Detection of urease as a marker of bacterial colonization and gastritis. Amer. J. Gastroenterology. (1987), 82(4): 292-96.
4. Cutler AF. Testing for *Helicobacter pylori* in clinical practice. Am J. Med. 1996; 100:355-415.
5. Anand BS, Raed AK, Malaty HM, et al. Loe point prevalence of peptic ulcer in normal individual with *Helicobacter pylori* infection. Am J Gastroenterol. 1996;91:1112-1115.

1. Intended Use

The NADAL® H. pylori Antigen Rapid Test Device (Feces) is a rapid visual immunoassay for the qualitative presumptive detection of *Helicobacter pylori* antigens in human fecal specimens. This kit is intended to be used as an aid in the diagnosis of H. pylori infection.

2. Introduction and/or Clinical Significances

Helicobacter pylori, also known as *Campylobacter pylori*, is a spiral-shaped Gram negative bacteria with a typical flagellum. The bacteria infects the gastric mucosa by secreting the Vacuolating cytotoxin A (VacA) and causes several gastric enteric diseases such as non-ulcerous dyspepsia, gastric and duodenal ulcer and active gastritis. H. pylori is classified as carcinogen agent type I, infection by which may increase the risk of stomach adenocarcinoma.

Many H. pylori strains have been isolated. The one expressing the cytotoxin associated gene A (CagA) antigen is found to be highly immunogenic and is of utmost clinical importance because of its association with cytotoxic factors. It is widely reported in many literature articles that the risk of gastric cancer in infected patients, who are positive for the antibody against CagA, is up to five times higher than the reference group infected with a CagA negative bacterial strain. The presence of the CagA antibody determines the persistence of H. pylori infection and ulceration.

The CagA antigen is able to bind to other antibodies such as CagII and CagC and lead to a sudden inflammatory response, which may provoke peptic ulceration and allergic episodes, and decrease the efficacy of therapy.

At present several invasive and non-invasive approaches are available to detect H. pylori infection. Invasive methodologies, which require endoscopy of the gastric mucosa and traditional urease investigation, are cost- and time-consuming for an accurate final diagnosis.

Non-invasive methods include classical ELISA, immunoblotting assays and Breath Test. But the latter is extremely complicate and highly unselective.

3. Principle of the Test

The NADAL® H. pylori Antigen Rapid Test Device (Feces) has been designed to detect *Helicobacter pylori* through visual interpretation of color development in the internal strip. The membrane was immobilized with anti-H. pylori antibodies on the test region. During the test, the specimen is allowed to react with colored anti-H. pylori antibodies colloidal gold conjugates, which were precoated on the sample pad of the test. The mixture then moves on the membrane by a capillary action, and interact with reagents on the membrane.

If there are enough H. pylori antigens in specimens, a colored band will form at the test region of the membrane. Presence of this colored band indicates a positive result, while its absence indicates a negative result. Appearance of a colored band at the control region serves as a procedural control. This indicates that proper volume of specimen has been added and membrane wicking has occurred.

4. Reagents and Materials Supplied

- Individually packed test devices
- Specimens dilution tube with buffer

- Package insert
- Specimen collection cards

5. Additional Required Materials

- Centrifuge
- Timer

6. Storage & Stability

The kit should be stored at 2-30°C until the expiry date printed on the sealed pouch.

Do not freeze.

The test must remain in the sealed pouch until use.

Care should be taken to protect the components of the kit from contamination. Do not use if there is evidence of microbial contamination or precipitation. Biological contamination of dispensing equipments, containers or reagents can lead to false results.

7. Warnings and Precautions

- For professional *in vitro* diagnostic use only.
- Do not use after the expiration date indicated on the package. Do not use the test if the foil pouch is damaged. Do not reuse tests.
- This kit contains products of animal origin. Certified knowledge of the origin and/or sanitary state of the animals does not completely guarantee the absence of transmissible pathogenic agents. It is therefore, recommended that these products be treated as potentially infectious, and handled by observing usual safety precautions (e.g., do not ingest or inhale).
- Avoid cross-contamination of specimens by using a new specimen collection container for each specimen obtained.
- Read the entire procedure carefully prior to testing.
- Do not eat, drink or smoke in the area where the specimens and kits are handled. Handle all specimens as if they contain infectious agents. Observe established precautions against microbiological hazards throughout the procedure and follow standard procedures for the proper disposal of specimens. Wear protective clothing such as laboratory coats, disposable gloves and eye protection when specimens are assayed.
- Humidity and temperature can adversely affect results.
- Used testing materials should be discarded according to local regulations.

8. Specimen Collection and Preparation


- The NADAL® H. pylori Antigen Rapid Test Device (Feces) is intended only for use with human fecal specimens.
- Perform the testing immediately after the specimen collection. Do not leave the specimens at room temperature for prolonged periods. Specimens may be stored at 2-8°C for up to 72 hours.
- Bring specimens to room temperature prior to testing.
- Pack the specimens in compliance with applicable regulations for transportation of etiological agents, in case they need to be shipped.

9. Procedure of the Test

Bring tests, specimens, buffer and/or controls to room temperature (15-30°C) before use.

Specimen collection and pre-treatment:

1. Use the specimens collection cards provided in the kit for specimens collection. Follow the operation procedure written on it for instructions. Other clean dry containers could also be used for the same purpose. Best results will be obtained if the assay is performed within 6 hours after collection.
2. Unscrew and remove the dilution tube applicator. Be careful not to spill or spatter solution from the tube. Collect specimens by inserting the applicator stick into at least 3 different sites of the feces to collect approximately 50 mg of feces (equivalent to 1/4 of a pea).
3. Push the applicator back into the tube and screw the cap tightly. Be careful not to break the tip of the dilution tube.
4. Shake the specimen collection tube vigorously to mix the specimen and the extraction buffer.

Testing

1. Remove the test from its sealed pouch, and place it on a clean, level surface. Label the test with patient or control identification. To obtain a best result, the assay should be performed within one hour.
2. Using a piece of tissue paper, break the tip of the dilution tube. Hold the tube vertically and dispense 2 drops of solution into the specimen well (S) of the test device. Avoid trapping air bubbles in the specimen well (S), and do not drop any solution in observation window. As the test begins to work, you will see color move across the membrane.
3. Wait for the colored band(s) to appear. The result should be read at 10 minutes. Do not interpret the result after 20 minutes.

Note: If the specimen does not migrate (presence of particles), centrifuge the extracted specimens contained in the extraction buffer vial. Collect 80 µL of supernatant, dispense into the specimen well (S) of a new test device and start afresh following the instructions mentioned above.

10. Interpretation of the Results

Positive

Two colored bands appear on the membrane. One band appears in the control region (C) and another band appears in the test region (T).

Negative

Only one colored band appears, in the control region (C). Non apparent colored band appears in the test region (T).

Invalid

Control band fails to appear. Results from any test which has not produced a control band at the specified read time must be discarded. Please review the procedure and repeat with a new test. If the problem persists, discontinue using the kit immediately and contact your local distributor.

Note:

The intensity of color in the test region (T) may vary depending on the concentration of analytes present in the specimen. Therefore, any shade of color in the test region should be considered positive. Note that this is a qualitative test only, and cannot determine the concentration of analytes in the specimen.

Insufficient specimen volume, incorrect operating procedure or expired tests are the most likely reasons for control band failure.

11. Quality Control

Internal procedural controls are included in the test. A colored band appearing in the control region (C) is considered an internal positive procedural control, confirming sufficient specimen volume and correct procedural technique.

External controls are not supplied with this kit. It is recommended that positive and negative controls be tested as a good laboratory practice to confirm the test procedure and to verify proper test performance.

12. Limitations

- The NADAL® H. pylori Antigen Rapid Test Device (Feces) is for professional *in vitro* diagnostic use, and should be used for the qualitative detection of *Helicobacter pylori* only.
- Following certain antibiotic treatments, the concentration of H. pylori antigens may decrease to the concentration below the minimum detection level of the test. Therefore, diagnosis should be made with caution during antibiotic treatment.
- As with all diagnostic tests, a definitive clinical diagnosis should not be based on the results of a single test, but should only be made by the physician after all clinical and laboratory findings have been evaluated.

13. Performance Characteristics

The NADAL® H.pylori Antigen Test vs. Biopsy/Histology/RUT

Biopsy/ Histology/ RUT	NADAL® H. pylori Antigen Test		
	+	-	Total
+	132	0	132
-	0	154	154
	132	154	286

Relative Sensitivity: >99.9% (97.3% - 100.0%)*

Relative Specificity: >99.9% (97.6% - 100.0%)*

Overall Agreement: >99.9% (98.7% - 98.8%)*

*95% Confidence Interval

Specificity:

Cross reactivity with following organisms has been studied at 1.0×10^9 organisms/ml. The following organisms were found

negative when tested with the NADAL®H. pylori Antigen Test (Feces).

<i>Staphylococcus aureus</i>	<i>Pseudomonas aeruginosa</i>	<i>Enterococcus faecalis</i>
Group <i>C Streptococcus</i>	<i>Klebsiella pneumoniae</i>	<i>Branhamella catarrhalis</i>
<i>Candida albicans</i>	<i>Proteus mirabilis</i>	<i>Acinetobacter spp</i>
<i>Salmonella choleraesuis</i>	<i>Gardnerella vaginalis</i>	<i>Acinetobacter calcoaceticus</i>
<i>E. coli</i>	<i>Chlamydia trachomatis</i>	<i>Neisseria gonorrhoea</i>
Group <i>B Streptococcus</i>	<i>Proteus vulgaris</i>	<i>Enterococcus faecium</i>
<i>Hemophilus influenzae</i>	<i>Neisseria meningitidis</i>	<i>Rotavirus</i>

14. References

1. Marshall, BJ, McGeachie, DB, Rogers, PAR and Glancy, RG. Pyloric Campylobacter infection and gastroduodenal disease. Med. J. Australia. (1985), 149: 439-44.
2. Soll, AH. Pathogenesis of peptic ulcer and implications for therapy. New England J. Med. (1990), 322: 909-16.
3. Hazeli, SL, et al. Campylobacter pyloridis and gastritis I: Detection of urease as a marker of bacterial colonization and gastritis. Amer. J. Gastroenterology. (1987), 82(4): 292-96.
4. Cutler AF. Testing for Helicobacter pylori in clinical practice. Am J. Med. 1996; 100:355-415.
5. Anand BS, Raed AK, Malaty HM, et al. Loe point prevalence of peptic ulcer in normal individual with Helicobacter pylori infection. Am J Gastroenterol. 1996;91:1112-1115.

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14. Références

1. Marshall, BJ, McGeachie, DB, Rogers, PAR and Glancy, RG. Pyloric Campylobacter infection and gastroduodenal disease. Med. J. Australia. (1985), 149: 439-44.
2. Soll, AH. Pathogenesis of peptic ulcer and implications for therapy. New England J. Med. (1990), 322: 909-16.
3. Hazell, SL, et al. Campylobacter pyloridis and gastritis I: Detection of urease as a marker of bacterial colonization and gastritis. Amer. J. Gastroenterology. (1987), 82(4): 292-96.
4. Cutler AF. Testing for Helicobacter pylori in clinical practice. Am j. Med. 1996; 100:355-415.
5. Anand BS, Raed AK, Malaty HM, et al. Loe point prevalence of peptic ulcer in normal individual with Helicobacter pylori infection. Am J Gastroenterol. 1996;91:1112-1115.

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14. Bibliografia

1. Marshall, BJ, McGeachie, DB, Rogers, PAR and Glancy, RG. Pyloric Campylobacter infection and gastroduodenal disease. Med. J. Australia. (1985), 149: 439-44.
2. Soll, AH. Pathogenesis of peptic ulcer and implications for therapy. New England J. Med. (1990), 322: 909-16.
3. Hazell, SL, et al. Campylobacter pyloridis and gastritis I: Detection of urease as a marker of bacterial colonization and gastritis. Amer. J. Gastroenterology. (1987), 82(4): 292-96.
4. Cutler AF. Testing for Helicobacter pylori in clinical practice. Am J. Med. 1996; 100:355-415.
5. Anand BS, Raed AK, Malaty HM, et al. Loe point prevalence of peptic ulcer in normal individual with Helicobacter pylori infection. Am J Gastroenterol. 1996;91:1112-1115.

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nal von minden GmbH