NADAL® Dengue NS1 Ag+IgG/IgM Test (test cassette)

REF 532016N-25

DE Gebrauchsanweisung	2	CZ Návod k použití	37
EN Instructions for use	7	FI Käyttöohje	42
FR Instructions d'utilisation	12	SE Användarinstruktioner	47
ES Instrucciones de uso	17	Symbols	55
IT Istruzioni per l'uso	22	Our Teams	56
PL Sposób użycia	27		
PT Instruções de Utilização	32		

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15. Referenzen

- Halstead SB, Selective primary health care: strategies for control of disease in the developing world: XI, Dengue. Rev. Infect. Dis. 1984; 6:251-264.
- Halstead SB, Pathogenesis of dengue: challenges to molecular biology. Science 1988; 239:476-481.
- Ruechusatsawat K, et al. Daily observation of antibody levels among dengue patients detected by enzyme-linked immunosorbent assay (ELISA). Japanese J. Trop. Med. Hygiene 1994; 22: 9-12.
- Lam SK. Dengue haemorrhagic fever. Rev. Med. Micro. 1995; 6:39-48.
- Dengue haemorrhagic fever: diagnosis, treatment, prevention and control. 2nd edition. Geneva: World Health Organization.
- Yamada K, et al. Antibody responses determined for Japanese dengue fever patients by neutralization and hemagglutination inhibition assays demonstrate cross-reactivity between dengue and Japanese encephalitis viruses. Clin Diagn Lab Immunol. 2003 Jul; 10(4): 725-8.
- Dobler G, et al. Cross reactions of patients with acute dengue fever to tick-borne encephalitis. Wien Med Wochenschr (in German). 1997; 147(19-20): 463-4.
- Makino Y, et al. Studies on serological cross-reaction in sequential flavivirus. Microbiol Immunol. 1994; 38(12): 951-5.

Rev. 0, 2016-09-06 OM/UJA

1. Intended Use

The NADAL® Dengue NS1 Ag + IgG/IgM Test is a rapid chromatographic immunoassay for the simultaneous qualitative detection and differentiation of the dengue NS1 antigen and anti-dengue virus IgG and IgM in human whole blood, serum or plasma. The test is intended for use as an aid in the diagnosis of dengue infections and designed for professional use only.

2. Introduction and Clinical Significance

Dengue virus is a flavivirus, transmitted by tiger mosquitos *Aedes aegypti* and *Aedes albopictus*. It is widely distributed throughout the tropical and subtropical areas of the world¹ and causes up to 100 million infections annually.² Classic Dengue infection is characterised by a sudden onset of fever, intense headache, myalgia, arthralgia and a rash.

Primary dengue infection causes IgM to increase to a detectable level in 3 to 5 days after the onset of fever. IgM generally persist for 30 to 90 days.³ Most dengue patients in endemic regions have secondary infections,⁴ resulting in high levels of specific IgG prior to or simultaneous with IgM response.⁵ Therefore, the detection of specific anti-dengue IgM and IgG can also help to distinguish between primary and secondary infections.

NS1 is one of 7 Dengue virus non-structural proteins which are thought to be involved in viral replication. NS1 exists as a monomer in its immature form, but is rapidly processed in the endoplasmic reticulum to form a stable dimer. A small amount of NS1 remains associated with intracellular organelles where it is thought to be involved in viral replication. The rest of NS1 is found either to be associated with the plasma membrane or to be secreted as a soluble hexadimer. NS1 is essential for viral viability but its precise biological function is unknown. Antibodies produced in response to NS1 in viral infection can cross-react with cell surface antigens on epithelial cells and platelets. This has implications in the development of dengue hemorrhagic fever.

The NADAL® Dengue NS1 Ag + IgG/IgM Test is a rapid test that utilises a combination of coated dengue antibodies bound to coloured particles for the detection of dengue NS1 antigens and dengue antigens bound to coloured particles for the detection of anti-dengue IgG and anti-dengue IgM in human whole blood, serum or plasma.

3. Test Principle

The NADAL® Dengue NS1 Ag + IgG/IgM Test is a qualitative membrane-based immunoassay for the detection of anti-dengue antibodies and dengue NS1 antigens in human whole blood, serum or plasma.

The NADAL® Dengue NS1 Ag + IgG/IgM Test contains two test strips (left side: Dengue NS1 Ag Test; right side: Dengue IgG/IgM Test).

The Dengue IgG/IgM Test on the right side consists of two components. In the IgG component anti-human-IgG are precoated in the testline region 'IgG'. During testing, the specimen reacts with dengue antigen-coated particles in the test cassette. The mixture then migrates along the membrane chromatographically by capillary action and reacts with the anti-human IgG antibodies in the test line region 'IgG'. If the

specimen contains anti-dengue IgG, a coloured line appears in the test line region 'IgG'.

In the IgM component, anti-human IgM antibodies are precoated in the test line region 'IgM'. If anti-dengue IgM are present in the specimen, they react with the dengue antigen-coated particles in the test cassette. The complex is then captured by the anti-human IgM antibodies, forming a coloured line in the test line region 'IgM'.

If the specimen does not contain dengue antibodies, no coloured line appears in either of the test line regions, indicating a negative result.

The Dengue NS1 Ag Test on the left side:

During testing, the specimen reacts with dengue antibody-conjugates on the membrane of the internal test strip. The gold-antibody conjugates bind to dengue antigens in the specimen, which in turn are bound by anti-dengue NS1 antibodies precoated on the membrane. As the mixture moves along the membrane, the anti-dengue NS1 antibodies bind to the antibody-antigen complex, causing pale or dark pink line to form in the test line region 'NS1' of the membrane. The colour intensity of the lines may vary depending upon the amount of antigens present in the sample. The appearance of a pink line in the test line region 'NS1' should be considered as a positive result.

If the specimen does not contain dengue antigens, no coloured line appears in the test line region, indicating a negative result.

The appearance of both coloured lines in the control line regions (C) serve as a procedural control, indicating that proper volume of specimen has been added and membrane wicking has occurred.

4. Reagents and Materials Supplied

- 25 NADAL® Dengue NS1 Ag + IgG/IgM test cassettes
- 25 disposable pipettes (5 µl)
- 25 disposable pipettes (25 µl)
(the provided 5 µl and 25 µl pipettes dispense approximately 5 µl and 25 µl respectively in one drop even if more blood is aspirated into the pipette)
- 1 Buffer (3 ml)
- 1 package insert

5. Additional Materials Required

- Specimen collection containers
- Centrifuge (for serum/plasma)
- Timer
- Lancet (for fingerstick whole blood only)
- Capillary tube (75 µl, if needed)
- Micropipette (if needed)

6. Storage & Stability

NADAL® Dengue NS1 Ag + IgG/IgM Test kits should be stored at 2-30°C and used before the expiry date printed on the packaging. Test cassettes are stable until the expiration date printed on the foil pouches. Test cassettes should remain in sealed pouches until use. Do not freeze tests. Do not use tests beyond the expiration date.

7. Warnings and Precautions

- For professional *in-vitro* diagnostic use only.

- Carefully read through the test procedure prior to testing.
- Do not use the test beyond the expiration date indicated on the package.
- Do not use the test if the foil pouch is damaged.
- Do not reuse tests.
- Do not add samples to the reaction area (result area).
- In order to avoid contamination, do not touch the reaction area (result area).
- Avoid cross-contamination of specimens by using a new specimen collection container for each specimen obtained.
- Do not eat, drink or smoke in the area where specimens and test kits are handled.
- Wear protective clothing such as laboratory coats, disposable gloves and eye protection when specimens are being assayed.
- Handle all specimens as if they contain infectious agents. Observe established precautions for microbiological risks throughout all procedures and standard guidelines for the appropriate disposal of specimens.
- The test kit contains products of animal origin. Certified knowledge of the origin and/or sanitary state of the animals does not completely guarantee the absence of transmissible pathogenic agents. It is therefore recommended that these products be treated as potentially infectious, and handled in accordance with usual safety precautions (e.g., do not ingest or inhale).
- Humidity and temperature can adversely affect test results.
- Used testing materials should be discarded according to local regulations.

8. Specimen Collection and Preparation

The NADAL® Dengue NS1 Ag + IgG/IgM Test can be performed using whole blood, serum or plasma.

To collect fingerstick whole blood specimens:

- Wash the patient's hand with soap and warm water or clean it with an alcohol swab. Allow it to dry.
- Massage the hand without touching the puncture site by rubbing down the hand towards the fingertip of the middle or ring finger.
- Puncture the skin with a sterile lancet. Wipe away the first drop of blood.
- Gently rub the hand from wrist to palm to finger to form a rounded drop of blood over the puncture site.

Venipuncture whole blood specimens

For serum or plasma testing, separate serum or plasma from blood as soon as possible to avoid hemolysis. Use only clear, non-hemolysed specimens.

Testing should be performed immediately after specimen collection. Do not leave specimens at room temperature for prolonged periods of time. Serum and plasma specimens can be stored at 2-8°C for up to 3 days. For long-term storage, specimens should be kept below -20°C. Whole blood collected by venipuncture should be stored at 2-8°C if the test is to be run within 2 days of collection. Do not freeze whole blood specimens. Fingerstick whole blood should be tested immediately. Bring specimens to room temperature prior to testing. Frozen specimens should be completely thawed and mixed well prior to testing. Specimens should not be frozen and thawed repeatedly. If specimens are to be shipped, they

should be packed in compliance with federal regulations for transportation of etiologic agents.

9. Test Procedure


Bring tests, specimens, buffer and/or controls to room temperature (15-30°C) prior to testing.

1. Remove the test cassette from the sealed foil pouch and place it on a clean and level surface. The test cassette should be used immediately after opening the foil pouch. Label the test cassette with the patient or control identification.

2. a) For serum or plasma specimen:

For Dengue NS1 Ag:

Holding the provided 25 µl dropper vertically, add 3 drops (approximately 75 µl) of serum or plasma to the sample well (S) of the test cassette.

and

For anti-Dengue IgG/IgM:

I. To use a pipette:

Holding the provided 5 µl pipette vertically, draw serum or plasma specimen up to the first widening (approximately 5 µl) and add it to the specimen well (S) of the test cassette. Then add 1 drop (approximately 40 µl) of buffer to the sample well (S). **Avoid trapping air bubbles in the specimen well.**

or

II. To use a micropipette:

Using a micropipette, dispense 5 µl of serum or plasma into the specimen well (S) of the test cassette. Then add 1 drop (approximately 40 µl) of buffer to the sample well (S).


b) For fingerstick/venipuncture whole blood specimen:

For Dengue NS1 Ag:

I. To use a pipette:


Holding the provided 25 µl dropper vertically, add 3 drops (approximately 75 µl) of whole blood to the sample well (S) of the test cassette. Then add 1 drop (approximately 40 µl) of buffer to the sample well (S).

or

II. To use a capillary tube:

Fill the capillary tube and transfer approximately 75 µl of fingerstick whole blood specimen into the sample well (S). Then add 1 drop (approximately 40 µl) of buffer to the sample well (S).

or

III. To use hanging drops:

Allow 3 hanging drops (approximately 75 µl) of fingerstick whole blood to fall into the sample well (S) of the test cassette. Then add 1 drop

(approximately 40 µL) of buffer to the sample well (S).

For anti-Dengue IgG/IgM:

I. To use a pipette:

Holding the provided 5 µL pipette vertically, draw whole blood specimen 1 cm above the first widening (approximately 10 µL) and add it to the specimen well (S) of the test cassette. Then add 1 drop (approximately 40 µL) of buffer to the sample well (S).

or

10 min

3. Start the timer.

4. Read the result after 10 minutes. Do not interpret the result after more than 20 minutes.

10. Result Interpretation

Positive for Dengue NS1 Ag:

One coloured line develops in each control line region 'C' and one coloured line develops in the test line region 'NS1'. The result may be indicative of current infection.

Positive for anti-Dengue IgG:

One coloured line develops in each control line region 'C' and one coloured line develops in the test line region 'IgG'. The result may be indicative of secondary dengue virus infection.

Positive for anti-Dengue IgM:

One coloured line develops in each control line region 'C' and one coloured line develops in the test line region 'IgM'. The result may be indicative of primary dengue virus infection.

Positive for anti-Dengue IgG/IgM:

One coloured line develops in each control line region 'C' and one coloured line develops in the test line region 'IgG' and 'IgM' respectively. The result may be indicative of secondary dengue virus infection. The color intensities of the lines do not have to match.

Note: A combination of positive results for various parameters is possible.

Negative:

One coloured line appears in each control line region 'C'. No coloured lines appear in the test line regions 'NS1', 'IgG' and 'IgM'.


Invalid:

At least one control line 'C' fails to appear. Results from any test which has not produced all control lines at the specified reading time must be discarded.

Insufficient specimen volume, incorrect operating procedure or expired tests are the most likely reasons for control line failure. Please review the procedure and repeat the test with a new test cassette. If the problem persists, discontinue using the test kit immediately and contact your distributor.

Note: The colour intensity in the test line region 'NS1' and/or 'IgG' and/or 'IgM' may vary depending on the concentration of the Dengue NS1 antigens and/or IgG and/or IgM present in the specimen. Therefore, any shade of colour in the test line region should be considered positive.

11. Quality Control

An internal procedural control is included in the test cassette:

A coloured line appearing in each control line region (C) is considered an internal procedural control. It confirms sufficient specimen volume, adequate membrane wicking and correct procedural technique.

Good laboratory practice (GLP) recommends the use of control materials to ensure proper test kit performance.

12. Limitations

- The NADAL® Dengue NS1 Ag + IgG/IgM Test only detects the presence of dengue NS1 antigens and anti-dengue antibodies in the specimen and should not be used as the sole criterion for the diagnosis of dengue.
- The sections 'Test Procedure' and 'Result Interpretation' should be followed closely when testing for the presence of dengue antigen and anti-dengue antibody in whole blood, serum or plasma from individual subjects. Failure to follow the procedure may give inaccurate results.
- A negative test result for dengue NS1 antigens does not preclude the possibility of exposure to or infection with Dengue viruses.
- A negative test result for dengue NS1 antigens can occur if the quantity of dengue antigens present in the specimen is below the detection limit of the test or if the detected dengue antigens were not present during the stage of disease in which the sample was collected.
- Specimens containing unusually high titer of heterophile antibodies or rheumatoid factor may affect expected results.
- If symptoms persist, while the result obtained with the NADAL® Dengue NS1 Ag + IgG/IgM Test is negative or non-reactive, it is recommended to re-sample the patient a few days later or test with an alternative testing method, such as PCR or ELISA.

- The results obtained with the NADAL® Dengue NS1 Ag + IgG/IgM Test should only be interpreted in conjunction with other diagnostic procedures and clinical findings.
- After the early onset of fever, anti-dengue IgM concentrations may be below detectable levels. For primary infection, an IgM antibody-capture enzyme-linked immunosorbent assay (MAC-ELISA) showed that 80% of the tested dengue patients exhibited detectable levels of IgM by the fifth day after infection and 99% of the patients tested IgM positive by day 10. Therefore, it is recommended that patients be tested within this period of time. As for the secondary infection, a low molar fraction of anti-dengue IgM and a high molar fraction of anti-dengue IgG that are broadly reactive to flaviviruses, characterise the antibodies.⁵ The test line 'IgM' may be faint, and a cross-reaction in the test line region 'IgG' may occur.
- Serological cross-reactivity across the flavivirus group (Dengue 1, 2, 3 & 4, St. Louis encephalitis, West Nile virus, Japanese encephalitis and yellow fever viruses) is common.^{6,7,8} Positive results should be confirmed using other methods.
- The continued presence or absence of antibodies cannot be used to determine the success or failure of therapy.
- Results from immunosuppressed patients should be interpreted with caution.

13. Expected Values

The NADAL® Dengue NS1 Ag + IgG/IgM Test has been compared with a leading commercially available Dengue Ag EIA. The correlation between the two systems is 96%.

Primary dengue infection is characterised by the presence of detectable IgM 3-5 days after the onset of infection. Secondary dengue infection is characterised by the elevation of dengue-specific IgG. In the majority of cases, this is accompanied by elevated levels of IgM.

The NADAL® Dengue NS1 Ag + IgG/IgM Test has been compared with a leading commercially available Dengue ELISA, demonstrating sensitivity of 83.3% for IgM in primary infection and >99.9% for IgG in secondary infection.

14. Performance Characteristics

Sensitivity and Specificity

A seroconversion panel comprising clinical samples obtained from a population of symptomatic and asymptomatic individuals was tested for dengue NS1 antigens and anti-Dengue IgG/IgM, using the NADAL® Dengue NS1 Ag + IgG/IgM Test. The results for Dengue NS1 Ag were confirmed by a leading, commercially available Dengue Ag EIA, whereas the results for anti-Dengue IgG/IgM were confirmed by a leading, commercially available Dengue ELISA.

The results for anti-Dengue IgG/IgM show that the overall relative sensitivity for the primary and secondary infection of the NADAL® Dengue NS1 Ag + IgG/IgM Test is 95.7%, the relative specificity is >99.9% and the relative accuracy is 99.3%.

Dengue primary infection for IgG/IgM test results

		Dengue ELISA			
		Positive		Negative	
NADAL® Dengue NS1 Ag + IgG/IgM Test (IgG/IgM)	Results		IgM		IgG
	Positive	IgM	15	0	0
		IgG	3	0	0
Negative		0	0	0	

Relative sensitivity (IgM): 83.3%

Dengue secondary infection for IgG/IgM test results

		Dengue ELISA			
		Positive		Negative	
NADAL® Dengue NS1 Ag + IgG/IgM Test (IgG/IgM)	Results		IgM		IgG
	Positive	IgM	37	0	0
		IgG	15	52	0
Negative		0	0	0	

Relative sensitivity (IgM): 71.2%

Relative sensitivity (IgG): >99.9%

Non-dengue infection for IgG/IgM test results

		Dengue ELISA			
		Positive		Negative	
NADAL® Dengue NS1 Ag + IgG/IgM Test (IgG/IgM)	Results		IgM		IgG
	Positive	IgM	0	0	0
		IgG	0	0	0
Negative		0	0	338	

Relative sensitivity (negative): >99.9%

Relative sensitivity: $(15+52)/(18+52) = 95.7\%$ (95%CI*: 88%~99.1%)

Relative specificity: $338/338 >99.9\%$ (95%CI*: 99.1%~100%)

Accuracy: $(15+52+338)/(18+52+338) = 99.3\%$ = (95%CI*: 97.9%~99.8%)

*Confidence Intervals

The results for Dengue NS1 Ag show that the relative sensitivity of the NADAL® Dengue NS1 Ag + IgG/IgM Test is 95.8% and the relative specificity is 96.1%.

Dengue NS1 Ag test results:

		Dengue Ag ELISA		
		+	-	Total
NADAL® Dengue NS1 Ag + IgG/IgM Test (NS1 antigens)	+	137	8	145
	-	6	200	206
	Total	143	208	351

Relative sensitivity: $137/143*100\% = 95.8\%$ (95%CI*: 91.1%~98.4%)

Relative specificity: $200/208*100\% = 96.1\%$ (95%CI*: 92.6%~98.4%)

Accuracy: $(137+208)/(137+6+8+200)*100\% = 96\%$ (95%CI*: 93.4%~97.8%)

*Confidence Intervals

Precision

Intra-Assay

Within-run precision has been determined by using 15 replicates of four specimens: a negative, a low positive, a middle positive and a high positive. The specimens were correctly identified >99% of the time.

Inter-Assay

Between-run precision has been determined by 15 independent assays on the same four specimens: a negative, a low positive, a middle positive and a high positive. Three different lots of the NADAL® Dengue NS1 Ag + IgG/IgM Test were tested using these specimens. The specimens were correctly identified >99% of the time.

Cross-reactivity

The NADAL® Dengue NS1 Ag + IgG/IgM Test has been tested with HAMA, RF, HBsAg, HBsAb, HBeAg, HBeAb, HBcAb, Syphilis, HIV, HCV, *H. pylori*, MONO, CMV, Rubella and TOXO positive specimens. The results showed no cross-reactivity.

Interfering substances

The following potentially interfering substances were added to dengue negative and positive specimens.

Acetaminophen: 20 mg/dl	Caffeine: 20 mg/dL
Acetylsalicylic acid: 20 mg/dL	Gentisic acid: 20 mg/dL
Ascorbic acid: 2 g/dL	Albumin: 2 g/dL
Creatine: 200 mg/dL	Hemoglobin 1000 mg/dL
Bilirubin: 1 g/dL	Oxalic acid: 60 mg/dL

None of the substances at the concentration tested interfered in the assay.

15. References

1. Halstead SB, Selective primary health care: strategies for control of disease in the developing world: XI, Dengue. *Rev. Infect. Dis.* 1984; 6:251-264.
2. Halstead SB, Pathogenesis of dengue: challenges to molecular biology. *Science* 1988; 239:476-481.
3. Ruechusatsawat K, et al. Daily observation of antibody levels among dengue patients detected by enzyme-linked immunosorbent assay (ELISA). *Japanese J. Trop. Med. Hygiene* 1994; 22: 9-12.
4. Lam SK. Dengue haemorrhagic fever. *Rev. Med. Micro.* 1995; 6:39-48.
5. Dengue haemorrhagic fever: diagnosis, treatment, prevention and control. 2nd edition. Geneva: World Health Organization.
6. Yamada K, et al. Antibody responses determined for Japanese dengue fever patients by neutralization and hemagglutination inhibition assays demonstrate cross-reactivity between dengue and Japanese encephalitis viruses. *Clin Diagn Lab Immunol.* 2003 Jul; 10(4): 725-8.
7. Dobler G, et al. Cross reactions of patients with acute dengue fever to tick-borne encephalitis. *Wien Med Wochenschr* (in German). 1997; 147(19-20): 463-4.
8. Makino Y, et al. Studies on serological cross-reaction in sequential flavivirus. *Microbiol Immunol.* 1994; 38(12): 951-5.

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15. Bibliografia

1. Halstead SB, Selective primary health care: strategies for control of disease in the developing world: XI, Dengue. Rev. Infect. Dis. 1984; 6:251-264.

2. Halstead SB, Pathogenesis of dengue: challenges to molecular biology. Science 1988; 239:476-481.
3. Ruechusatsawat K, et al. Daily observation of antibody levels among dengue patients detected by enzyme-linked immunosorbent assay (ELISA). Japanese J. Trop. Med. Hygiene 1994; 22: 9-12.
4. Lam SK. Dengue haemorrhagic fever. Rev. Med. Micro. 1995; 6:39-48.
5. Dengue haemorrhagic fever: diagnosis, treatment, prevention and control. 2nd edition. Geneva: World Health Organization.
6. Yamada K, et al. Antibody responses determined for Japanese dengue fever patients by neutralization and hemagglutination inhibition assays demonstrate cross-reactivity between dengue and Japanese encephalitis viruses. Clin Diagn Lab Immunol. 2003 Jul; 10(4): 725-8.
7. Dobler G, et al. Cross reactions of patients with acute dengue fever to tick-borne encephalitis. Wien Med Wochenschr (in German). 1997; 147(19-20): 463-4.
8. Makino Y, et al. Studies on serological cross-reaction in sequential flavivirus. Microbiol Immunol. 1994; 38(12): 951-5.

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15. Referências

- Halstead SB. Selective primary health care: strategies for control of disease in the developing world: XI, Dengue. Rev. Infect. Dis. 1984; 6:251-264.
- Halstead SB. Pathogenesis of dengue: challenges to molecular biology. Science 1988; 239:476-481.
- Ruechusatsawat K, et al. Daily observation of antibody levels among dengue patients detected by enzyme-linked immunosorbent assay (ELISA). Japanese J. Trop. Med. Hygiene 1994; 22: 9-12.
- Lam SK. Dengue haemorrhagic fever. Rev. Med. Micro. 1995; 6:39-48.
- Dengue haemorrhagic fever: diagnosis, treatment, prevention and control. 2nd edition. Geneva: World Health Organization.
- Yamada K, et al. Antibody responses determined for Japanese dengue fever patients by neutralization and hemagglutination inhibition assays demonstrate cross-reactivity between dengue and Japanese encephalitis viruses. Clin Diagn Lab Immunol. 2003 Jul; 10(4): 725-8.
- Dobler G, et al. Cross reactions of patients with acute dengue fever to tick-borne encephalitis. Wien Med Wochenschr (in German). 1997; 147(19-20): 463-4.
- Makino Y, et al. Studies on serological cross-reaction in sequential flavivirus Microbiol Immunol. 1994; 38(12): 951-5.

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15. Reference

1. Halstead SB, Selective primary health care: strategies for control of disease in the developing world: XI, Dengue. Rev. Infect. Dis. 1984; 6:251-264.
2. Halstead SB, Pathogenesis of dengue: challenges to molecular biology. Science 1988; 239:476-481.
3. Ruechusatsawat K, et al. Daily observation of antibody levels among dengue patients detected by enzyme-linked immunosorbent assay (ELISA). Japanese J. Trop. Med. Hygiene 1994; 22: 9-12.
4. Lam SK. Dengue haemorrhagic fever. Rev. Med. Micro. 1995; 6:39-48.
5. Dengue haemorrhagic fever: diagnosis, treatment, prevention and control. 2nd edition. Geneva: World Health Organization.
6. Yamada K, et al. Antibody responses determined for Japanese dengue fever patients by neutralization and hemagglutination inhibition assays demonstrate cross-reactivity between dengue and Japanese encephalitis viruses. Clin Diagn Lab Immunol. 2003 Jul; 10(4): 725-8.
7. Dobler G, et al. Cross reactions of patients with acute dengue fever to tick-borne encephalitis. Wien Med Wochenschr (in German). 1997; 147(19-20): 463-4.
8. Makino Y, et al. Studies on serological cross-reaction in sequential flavivirus Microbiol Immunol. 1994; 38(12): 951-5.

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

15. Referenser

1. Halstead SB, Selective primary health care: strategies for control of disease in the developing world: XI, Dengue. Rev. Infect. Dis. 1984; 6:251-264.
2. Halstead SB, Pathogenesis of dengue: challenges to molecular biology. Science 1988; 239:476-481.
3. Ruechusatsawat K, et al. Daily observation of antibody levels among dengue patients detected by enzyme-linked immunosorbent assay (ELISA). Japanese J. Trop. Med. Hygiene 1994; 22: 9-12.
4. Lam SK. Dengue haemorrhagic fever. Rev. Med. Micro. 1995; 6:39-48.
5. Dengue haemorrhagic fever: diagnosis, treatment, prevention and control. 2nd edition. Geneva: World Health Organization.
6. Yamada K, et al. Antibody responses determined for Japanese dengue fever patients by neutralization and hemagglutination inhibition assays demonstrate cross-reactivity between dengue and Japanese encephalitis viruses. Clin Diagn Lab Immunol. 2003 Jul; 10(4): 725-8.
7. Dobler G, et al. Cross reactions of patients with acute dengue fever to tick-borne encephalitis. Wien Med Wochenschr (in German). 1997; 147(19-20): 463-4.
8. Makino Y, et al. Studies on serological cross-reaction in sequential flavivirus. Microbiol Immunol. 1994; 38(12): 951-5.

Rev. 0, 2016-09-06 MM

Symbol	Deutsch	English	Français	Español	Italiano	Polski
	CE Konformitätszeichen	CE marking of conformity	Conforme aux normes européennes	Conformidad europea	Conformità europea	Znak zgodności CE
	Gebrauchsanweisung beachten	Consult instructions for use	Consulter la notice d'utilisation	Consúltense las instrucciones de uso	Consultare le istruzioni per l'uso	Przestrzegać instrukcji obsługi
	<i>in-vitro</i> -Diagnostika	<i>in-vitro</i> diagnostic medical device	Dispositif médical de diagnostic <i>in-vitro</i>	Producto sanitario para diagnóstico <i>in-vitro</i>	Dispositivo medico-diagnostico <i>in-vitro</i>	Tylko do diagnostyki <i>in-vitro</i>
	Temperaturbegrenzung	Temperature limitation	Limites de température	Límite de temperatura	Limiti di temperatura	Temperatura przechowywania
	Chargenbezeichnung	Batch code	Numéro de lot	Código de lote	Codice lotto	Numer serii
	Nicht zur Wiederverwendung	Do not reuse	Ne pas réutiliser	No reutilizar	Non riutilizzare	Tylko do jednorazowego użytku
	Verwendbar bis	Use by	Utiliser jusqu'au	Fecha de caducidad	Utilizzare entro	Data ważności
	Bestellnummer	Catalogue Number	Référence du catalogue	Número de catálogo	Riferimento di Catalogo	Numer katalogowy
	Hersteller	Manufacturer	Fabricant	Fabricante	Fabbricante	Producent
	Ausreichend für <n> Ansätze	Sufficient for <n> tests	Suffisant pour "n" tests	Suficiente para <n> utilizations	Sufficiente per "n" saggi	Wystarczający na <n> Powtórzeń

Symbol	Português	Český	Suomi	Svenskt	Nederlands	Dansk	Norsk
	Conformidade com as normas europeias	CE certifikát	CE-merkitty	CE-märkning	CE-markering	CE-mærkning	CE standardisert
	Consultar as instruções de utilização	Viz návod k použití	Katso käyttöohjetta	Läs bruksanvisningen	Raadpleeg de gebruiksaanwijzing	Se brugsanvisningen	Les bruksanvisning nøye
	Dispositivo médico para diagnóstico <i>in-vitro</i>	Diagnostický zdravotnický prostředek <i>in-vitro</i>	<i>in-vitro</i> - diagnostiikkaan tarkoitettu lääkinnällinen laite	Medicinteknisk produkt avsedd för <i>in-vitro</i> -diagnostik	Medisch hulpmiddel voor <i>in-vitro</i> diagnostiek	Medicinsk udstyr til <i>in-vitro</i> -diagnostik	<i>in-vitro</i> diagnostic medisinsk enhet
	Limites de temperatura	Teplotní omezení	Lämpötilarajat	Temperaturbegränsning	Temperatuurlimiet	Temperaturbegrænsning	Temperaturbegrensning
	Código do lote	Kód šarže	Eräkoodi	Satsnummer	Code van de partij	Batchkode	Merking
	Não reutilizar	Pro jednorázové použití	Kertäkäyttöinen	Får inte återanvändas	Niet opnieuw gebruiken	Må ikke genbruges	Må ikke brukes om igjen
	Prazo de validade	Spotřebuje do	Käytettävä viimeistään	Används före	Houdbaar tot	Udløbsdato	Tidtaking
	Número de catálogo	Katalogov číslo	Luettelonumero	Listnummer	Catalogus nummer	Best il l ingsnummer	Katalog nummer
	Fabricante	Výrobce	Valmistaja	Tillverkare	Fabrikant	Fabrikant	Produsent
	Suficiente para <n> test	Dostačuje pro <n> testů	Lukumäärä <n> test	Räcker till <n> test	Voldoende voor <n> test	Tilstrækkeligt til <n> test	Tilstrækkelig for <n> tester

Our Teams

Germany:

Regensburg

Tel: +49 941 290 10-0
Fax: +49 941 290 10-50

Moers

Tel: +49 2841 99820-0
Fax: +49 2841 99820-1

Austria:

Tel: +49 941 290 10-29
Free Tel: 0800 291 565
Fax: +49 290 10-50
Free Fax: 0800 298 197

UK & Ireland:

Tel: +49 941 290 10-18
Free Tel – UK: 0808 234 1237
Free Tel – IRE: 1800 555 080
Fax: +49 290 10-50

France:

France Tel: 0800 915 240
France Fax: 0800 909 493

Switzerland

Swiss Tel: 0800 564 720
Swiss Fax: 0800 837 476

Belgium

Belgium Tel: 0800 718 82
Belgium Fax: 0800 747 07

Luxembourg

Lux. Tel: 800 211 16
Lux. Fax: 800 261 79

Spain:

Tel: +49 941 290 10-759
Free Tel: 900 938 315
Fax: +49 941 290 10-50
Free Fax: 900 984 992

Italy:

Tel: +49 941 290 10-34
Fax: +49 941 290 10-50

Poland:

Tel: +49 941 290 10-44
Free Tel: 00 800 491 15 95
Fax: +49 941 290 10-50
Free Fax: 00 800 491 15 94

Portugal:

Tel: +49 941 290 10-735
Tel. Verde: 800 849 230
Fax: +49 941 290 10-50
Fax Verde: 800 849 229

Netherlands:

Tel: +31 30 75 600
Free Tel: 0800 0222 890
Fax: +31 70 30 30 775
Free Fax: 0800 024 9519

Nordic countries (Finland, Norway, Sweden, Denmark):

Tel: +31 703075 607
Free Tel: +45 80 88 87 53
Tax: +31 703030 775

Laboratory Diagnostics Team:

Tel: +49 941 290 10-40
Fax: +49 941 290 10-50

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Carl-Zeiss-Strasse 12 • 47445 Moers • Germany

www.nal-vonminden.com • info@nal-vonminden.com

Fon: +49 2841 99820-0 • Fax: +49 2841 99820-1

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nal von minden GmbH