NADAL[®] Tuberculosis IgG/ IgM Test (test cassette)

REF 322003N-30

| | | | | |
|-----------|----------------------------|----|-----------|----|
| DE | Gebrauchsanweisung | 2 | Symbols | 23 |
| EN | Instructions for use | 5 | Our Teams | 24 |
| FR | Instructions d'utilisation | 8 | | |
| ES | Instrucciones de uso | 11 | | |
| PL | Sposób użycia | 14 | | |
| PT | Instruções de Utilização | 17 | | |

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14. Referenzen

- Schaaf, H. S., P. Botha, N. Beyers, R. P. Gie, H. A. et al. 1996. The 5-year outcome of multidrug resistant tuberculosis patients in the Cape Province of South Africa. Trop. Med. Int. Health 1:718-722
- Havir, D. V., and P. F. Barnes. 1999. Tuberculosis in patients with human immunodeficiency virus infection. N. Engl. J. Med. 340:367-373
- Dye L., Scheele S., V. Pathania, et al.: 1999 Global Burden of Tuberculosis Estimated Incidence, Prevalence, and Mortality by Country. WHO Global Surveillance and Monitoring Project. JAMA. 282:677-686.
- Kochi, A. 1991 The global tuberculosis situation and the new control strategy of the World Health Organization. Tubercle 72:1-6
- Merlin TL, Gibson DW, and DH Connor 1994. Tuberculosis, p400-404 in Rubein E and Farber JL (ed)-Pathology, 2nd edition. J.B. Lippincott Company
- Daniel, T. M. 1996. Immunodiagnosis of tuberculosis, p. 223-231. In W. N. Rom, and S. Garzy (ed.), Tuberculosis. Letter, Brown & Co., Boston, Mass.
- Wilkins, E. G. L. 1994. The serodiagnosis of tuberculosis, p. 367-379. In P. D. O. Davies (ed.), Clinical tuberculosis. Chapman & Hall, Ltd., London, England.
- Chan ED., Helfets L., and MD Iseman. 2000 Immunological Diagnosis of tuberculosis: a review. Tuberc. Lung Dis. 80: 131-140.
- Foulds, J., and R. O'Brien. 1998. New tools for the diagnosis of tuberculosis: the perspective of developing countries. Int. J. Tubercle. Lung Dis. 2: 778-783

Rev.2, 2015-11-25 OM/UJ

1. Intended Use

The NADAL® Tuberculosis IgG/IgM test is a lateral flow chromatographic immunoassay for the simultaneous detection and differentiation of anti-*Mycobacterium Tuberculosis (M.Tb)* IgM and IgG in human serum, plasma or whole blood. The test is intended to be used as a screening test and as an aid in the diagnosis of infection with *M.Tb*. Any reactive specimen with the NADAL® Tuberculosis IgG/IgM test should be confirmed with alternative testing method(s) and clinical findings.

2. Introduction and Clinical Significance

Tuberculosis (Tb) is a chronic, communicable disease principally caused by *M.Tb hominis* (Koch's bacillus), occasionally by *M.Tb bovis*. The lungs are the primary target but any organ may be infected.

The risk of Tb infection exponentially declined in the 20th century. However, the recent emergence of drug-resistant strains¹, particularly among patients with AIDS², has rekindled interest in Tb. The incidence of infection was reported to be around 8 million cases per year with a death rate of 3 million annually. Mortality exceeded 50% in some African countries with high HIV rates.^{3,4}

The initial clinical suspicion, radiographic findings with subsequent laboratory confirmation by sputum and by culture examination are traditional methods in the diagnosis of active Tb.^{5,6} However, these methods either lack sensitivity or are time consuming. Moreover, they are not suitable for smear-negative patients, patients who are unable to produce adequate sputum or patients suspected to have extrapulmonary Tb.

The NADAL® Tuberculosis IgG/IgM test is developed to alleviate these obstacles. The test detects anti-*M.Tb* IgM and IgG in serum, plasma or whole blood in 15 minutes. An IgM positive result indicates a fresh *M.Tb* infection, while an IgG positive result suggests a previous or chronic infection.⁷⁻⁹ The NADAL® Tuberculosis IgG/IgM test also detects anti-*M.Tb* IgM in patients vaccinated with BCG using *M.Tb* specific antigens. In addition, the test can be performed by untrained or minimally-skilled personnel without cumbersome laboratory equipment.

3. Test Principle

The NADAL® Tuberculosis IgG/IgM test is a lateral flow chromatographic immunoassay. The test cassette consists of: 1) a coloured conjugate pad containing *M.Tb* antigens conjugated with colloidal gold (*M.Tb* conjugates) and rabbit IgG-gold conjugates, 2) a nitrocellulose membrane strip containing two test lines M and G and a control line C. The M line is pre-coated with monoclonal anti-human IgM for the detection of anti-*M.Tb* IgM, the G line is pre-coated with antibodies for the detection of anti-*M.Tb* IgG and the C line is pre-coated with goat anti-rabbit IgG.

When an adequate volume of test specimen is dispensed into the sample well (S) of the test cassette, the specimen migrates by capillary action along the test cassette. If anti-*M.Tb* IgM are present in the specimen, they will bind to the *M.Tb* conjugates. The immunocomplex is then captured by the pre-coated anti-human IgM on the membrane, forming a coloured M line, indicating an anti-*M.Tb* IgM positive test result.

If anti-*M.Tb* IgG are present in the specimen, they will bind to the *M.Tb* conjugates. The immunocomplex is then captured by the pre-coated antibodies on the membrane, forming a coloured G line, indicating an anti-*M.Tb* IgG positive test result.

The absence of either test line (G and M) suggests a negative test result. The test contains an internal control (C line). A coloured C line of the immunocomplex of goat anti-rabbit IgG/rabbit IgG-gold conjugate must always appear regardless of the colour development of either of the test lines. Otherwise, the test result is invalid and the specimen should be retested with another test cassette.

4. Reagents and Materials Supplied

- 30 test cassettes, each sealed in a foil pouch with three items inside:
 - 1 test cassette
 - 1 plastic dropper
 - 1 desiccant
- Sample diluent (1 bottle, 5 mL)
- 1 package insert

5. Additional Materials Required

- Clock or timer
- Lancing device for whole blood test

6. Warnings and Precautions

- For professional *in-vitro* diagnostic use only.
- Carefully read through the test procedure prior to testing. Failure to follow the insert will lead to inaccurate test results.
- Do not open the sealed foil pouch unless ready to conduct the assay.
- Do not use expired test devices.
- Bring all reagents to room temperature (15-30°C) before testing.
- Do not substitute or mix components from different test kits.
- Do not use hemolised blood specimens for testing.
- Wear protective clothing such as laboratory coats, disposable gloves and eye protection when specimens are being assayed.
- Handle all specimens as if they contain infectious agents. Observe established precautions against microbiological hazards throughout the procedure and follow standard procedures for the proper disposal of specimens.
- Do not eat, drink or smoke in the area where specimens and test kits are being handled.
- Dispose of all specimens and materials used to perform the test as biohazardous waste.
- The test results should be read after 15 minutes after specimen is added to the sample well (S) of the test cassette. Reading test results after more than 15 minutes may give erroneous results. Positive results can be visible after as little as 1 minute.
- Do not perform the test in a room with strong air flow, i.e. an electric fan or strong air-conditioning.

7. Storage & Stability

All reagents are ready to use as supplied. Store unused test devices unopened at 2-30°C. If stored at 2-8°C, ensure that the test device is brought to room temperature before opening. The test device is stable through to the expiration date printed on the sealed pouch. Do not freeze the test kit or expose it to temperatures of over 30°C.

8. Specimen Collection and Preparation

Consider any materials of human origin as infectious and handle them following biosafety standard operating procedures.

Plasma

1. Collect blood specimen by venipuncture in a collection tube containing EDTA, citrate or heparin.
2. Separate plasma by centrifugation.
3. Carefully withdraw the plasma into a new pre-labeled tube.

Serum

1. Collect blood specimen by venipuncture in a collection tube containing no anticoagulants.
 2. Allow the blood to clot.
 3. Separate serum by centrifugation.
 4. Carefully withdraw the serum into a new pre-labeled tube.
- Test specimens as soon as possible after collecting. If specimens are not tested immediately after collection, store them at 2-8°C for up to 5 days. Specimens should be frozen at -20°C for longer storage.

Avoid multiple freeze-thaw cycles. Prior to testing, bring frozen specimens slowly to room temperature and mix them gently. Specimens containing visible particulate matter should be clarified by centrifugation before testing. Do not use samples demonstrating gross lipemia, gross hemolysis or turbidity in order to avoid interference with test result.

Whole blood

Whole blood of which blood can be obtained by either fingerstick or by venipuncture. Do not use any haemolysed blood for testing.


Whole blood specimens should be stored refrigerated (2-8°C) if not tested immediately. Specimens must be tested within 24 hours of collection.

9. Test Procedure

1. Bring the test components and specimen to room temperature if refrigerated or frozen. Once thawed, mix the specimen well prior to testing.
2. When ready to test, open the foil pouch and remove the test cassette. Place the test cassette on a clean, flat surface. Be sure to label the device with the specimen's ID number.

3. a) For whole blood test


- Dispense 1 drop of whole blood (about 40-50 µL) into the sample well (S).
- Then add 1 drop (about 35-50 µL) of sample diluent immediately making sure that there are no air bubbles.

b) For serum or plasma test

- Fill a dropper with serum/plasma specimen.

- Holding the dropper vertically, dispense 1 drop (about 30-45 µL) of specimen into the sample well (S).
 - Then add 1 drop (about 35-50 µL) of sample diluent immediately making sure that there are no air bubbles.
4. Set up the timer.
 5. Results can be read after 15 minutes. Positive results can be visible in as little as 1 minute.

Don't read results after more than 15 minutes. To avoid confusion, discard the test cassette after interpreting the result at the allotted time.

10. Result Interpretation

Positive result for anti-*M. Tb* IgM:

In addition to the line in the control line region (C), a coloured line develops in the test line region (M), indicating the presence of anti-*M. Tb* IgM antibodies.

Positive result for anti-*M. Tb* IgG:

In addition to the line in the control line region (C), a coloured line develops in the test line region (G), indicating the presence of anti-*M. Tb* IgG antibodies.

Positive result for anti-*M. Tb* IgG and anti-*M. Tb* IgM:

In addition to the line in the control line region (C), both G and M lines develop, indicating the presence of anti-*M. Tb* IgG and anti-*M. Tb* IgM.

Samples with positive results should be confirmed with alternative testing method(s) and clinical findings before a positive determination is made.

Negative result:

Only the control line C appears. The absence of both G and M lines indicates that no anti-*M. Tb* antibodies have been detected.

Invalid result:

The absence of the control line C, regardless of any colour development of the G and M lines, indicates a failure in the test procedure. The sample should be retested.

Note: Insufficient specimen volume, incorrect procedural techniques are the most likely reasons for control line failure. Review the procedure and repeat the test with a new test cassette. If the problem persists, discontinue using the test kit immediately and contact your distributor.

11. Quality Control

A procedural control is included in the test cassette. A coloured line developing in the control line region (C) is considered an internal procedural control. It confirms

sufficient specimen volume, adequate membrane wicking and correct procedural technique.

12. Performance Characteristics

1. Clinical Performance of IgM Tests

A total of 200 specimens from non-*Tb* patients and 35 specimens from patients under anti-*Tb* treatment were tested with the NADAL® Tuberculosis IgG/IgM test and a commercial TB IgM ELISA. The results are presented in the following table.

| IgM ELISA | NADAL® Tuberculosis IgG/IgM Test | | Total |
|-----------|----------------------------------|----------------|-------|
| | Positive (IgM) | Negative (IgM) | |
| Positive | 30 | 5 | 35 |
| Negative | 7 | 193 | 200 |
| Total | 37 | 198 | 235 |

Relative Sensitivity: 85.7%

Relative Specificity: 96.5%

Overall Agreement: 94.9%

2. Clinical Performance of IgG Tests

A total of 200 specimens from non-*Tb* patients and 35 specimens from patients under anti-*Tb* treatment were tested with the NADAL® Tuberculosis IgG/IgM test and a commercial TB IgG ELISA. The results are presented in the following table.

| IgG ELISA | NADAL® Tuberculosis IgG/IgM Test | | Total |
|-----------|----------------------------------|----------------|-------|
| | Positive (IgG) | Negative (IgG) | |
| Positive | 31 | 4 | 35 |
| Negative | 7 | 193 | 200 |
| Total | 38 | 197 | 235 |

Relative Sensitivity: 88.6%

Relative Specificity: 96.5%

Overall Agreement: 95.3%

13. Limitations

- The points "Procedure of the Test" and "Interpretation of the Results" should be followed closely when testing for the presence of antibodies to *M. Tb* in serum, plasma or whole blood. Failure to follow the procedure may give inaccurate results.
- The NADAL® Tuberculosis IgG/IgM test is limited to the qualitative detection of anti-*M. Tb* IgG and anti-*M. Tb* IgM in human serum, plasma or whole blood. The intensity of the test lines does not have a linear correlation with the antibody titer in the specimen.
- The test also detects antibodies to *M. bovis* and *M. africanum*.
- An IgG positive response may be detected in BCG vaccinated persons.
- A negative result indicates the absence of detectable antibodies to *M. Tb*. However, a negative test result does not preclude the possibility of exposure to, or infection with *M. Tb*.
- A negative result can occur if the quantity of antibodies to *M. Tb* present in the specimen is below the detection limit of the assay or the detected antibodies were not present in the sample during the stage of disease in which the sample was collected.
- Immunocompromised conditions such as HIV infection can reduce test sensitivity.

- Some specimens containing unusually high titers of heterophile antibodies or rheumatoid factors may affect expected results.
- The results obtained with this test should only be interpreted in conjunction with other diagnostic procedures and clinical findings.

14. References

1. Schaaf, H. S., P. Botha, N. Beyers, R. P. Gie, H. A. et al 1996. The 5-year outcome of multidrug resistant tuberculosis patients in the Cape Province of South Africa. Trop. Med. Int. Health 1:718-722
2. Havir, D. V., and P. F. Barnes. 1999. Tuberculosis in patients with human immunodeficiency virus infection. N. Engl. J. Med. 340:367-373
3. Dye L., Scheele S., V. Pathania, et al: 1999 Global Burden of Tuberculosis Estimated Incidence, Prevalence, and Mortality by Country. WHO Global Surveillance and Monitoring Project. JAMA. 282:677-686.
4. Kochi, A. 1991 The global tuberculosis situation and the new control strategy of the World Health Organization. Tubercle 72:1-6
5. Merlin TL, Gibson DW, and DH Connor 1994. Tuberculosis, p. 400-404 in Ruben E and Farber JL (ed)-Pathology, 2nd edition. J.B. Lippincott Company
6. Daniel, T. M. 1996. Immunodiagnosis of tuberculosis, p. 223-231. In W. N. Rom, and S. Garay (ed.), Tuberculosis. Letter, Brown & Co., Boston, Mass.
7. Wilkens, E. G. L. 1994. The serodiagnosis of tuberculosis, p. 367-379. In P. D. O. Davies (ed.), Clinical tuberculosis. Chapman & Hall, Ltd., London, England.
8. Chan ED., Helfets L, and MD Iseman. 2000 Immunological Diagnosis of tuberculosis: a review. Tuber. Lung Dis. 80: 131-140.
9. Foulds, J., and R. O'Brien. 1998. New tools for the diagnosis of tuberculosis: the perspective of developing countries. Int. J. Tubercle. Lung Dis. 2: 778-783

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14. Bibliographie

- Schaaf, H. S., P. Botha, N. Beyers, R. P. Gie, H. A. et, al 1996. The 5-year outcome of multidrug resistant tuberculosis patients in the Cape Province of South Africa. Trop. Med. Int. Health 1:718-722
- Havir, D. V., and P. F. Barnes. 1999. Tuberculosis in patients with human immunodeficiency virus infection. N. Engl. J. Med. 340:367-373
- Dye L, Scheele S, V. Pathania, et al: 1999 Global Burden of Tuberculosis Estimated Incidence, Prevalence, and Mortality by Country. WHO Global Surveillance and Monitoring Project. JAMA. 282:677-686.
- Kochi, A. 1991 The global tuberculosis situation and the new control strategy of the World Health Organization. Tubercle 72:1-6
- Merlin TL, Gibson DW, and DH Connor 1994. Tuberculosis, p. 400-404 in Rubein E and Farber JL (ed)-Pathology, 2nd edition. J.B. Lippincott Company
- Daniel, T. M. 1996. Immunodiagnosis of tuberculosis, p. 223-231. In W. N. Rom, and S. Garay (ed.), Tuberculosis. Letter, Brown & Co., Boston, Mass.
- Wilkins, E. G. L. 1994. The serodiagnosis of tuberculosis, p. 367-379. In P. D. O. Davies (ed.), Clinical tuberculosis. Chapman & Hall, Ltd., London, England.
- Chan ED., Heifets L., and MD Iseman. 2000 Immunological Diagnosis of tuberculosis: a review. Tuber. Lung Dis. 80: 131-140.
- Foulds, J., and R. O'Brien. 1998. New tools for the diagnosis of tuberculosis: the perspective of developing countries. Int. J. Tubercle. Lung Dis. 2: 778-783

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- Schaaf, H. S., P. Botha, N. Beyers, R. P. Gie, H. A. et, al 1996. The 5-year outcome of multidrug resistant tuberculosis patients in the Cape Province of South Africa. Trop. Med. Int. Health 1:718-722
- Havir, D. V., and P. F. Barnes. 1999. Tuberculosis in patients with human immunodeficiency virus infection. N. Engl. J. Med. 340:367-373
- Dye L., Scheele S., V. Pathania, et al: 1999 Global Burden of Tuberculosis Estimated Incidence, Prevalence, and Mortality by Country. WHO Global Surveillance and Monitoring Project. JAMA. 282:677-686.
- Kochi, A. 1991 The global tuberculosis situation and the new control strategy of the World Health Organization. Tubercle 72:1-6
- Merlin TL, Gibson DW, and DH Connor 1994. Tuberculosis, p400-404 in Rubein E and Farber JL (ed)-Pathology, 2nd edition. J.B. Lippincott Company
- Daniel, T. M. 1996. Immunodiagnosis of tuberculosis, p. 223-231. In W. N. Rom, and S. Garay (ed.), Tuberculosis. Letter, Brown & Co., Boston, Mass.
- Wilkins, E. G. L. 1994. The serodiagnosis of tuberculosis, p. 367-379. In P. D. O. Davies (ed.), Clinical tuberculosis. Chapman & Hall, Ltd., London, England.
- Chan ED., Hefets L, and MD Iseman. 2000 Immunological Diagnosis of tuberculosis: a review. Tuberc. Lung Dis. 80: 131-140.
- Foulds, J., and R. O'Brien. 1998. New tools for the diagnosis of tuberculosis: the perspective of developing countries. Int. J. Tubercle. Lung Dis. 2: 778-783

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14. Bibliografia

- Schaaf, H. S., P. Botha, N. Beyers, R. P. Gie, H. A. et al. 1996. The 5-year outcome of multidrug resistant tuberculosis patients in the Cape Province of South Africa. Trop. Med. Int. Health 1:718-722
- Havir, D. V., and P. F. Barnes. 1999. Tuberculosis in patients with human immunodeficiency virus infection. N. Engl. J. Med. 340:367-373
- Dye L, Scheele S, V. Pathania, et al: 1999 Global Burden of Tuberculosis Estimated Incidence, Prevalence, and Mortality by Country. WHO Global Surveillance and Monitoring Project. JAMA. 282:677-686.
- Kochi, A. 1991 The global tuberculosis situation and the new control strategy of the World Health Organization. Tubercle 72:1-6
- Merlin TL, Gibson DW, and DH Connor 1994. Tuberculosis, p400-404 in Rubein E and Farber JL (ed)-Pathology, 2nd edition. J.B. Lippincott Company
- Daniel, T. M. 1996. Immunodiagnosis of tuberculosis, p. 223-231. In W. N. Rom, and S. Garay (ed.), Tuberculosis. Letter, Brown & Co., Boston, Mass.
- Wilkins, E. G. L. 1994. The serodiagnosis of tuberculosis, p. 367-379. In P. D. O. Davies (ed.), Clinical tuberculosis. Chapman & Hall, Ltd., London, England.
- Chan ED., Helfets L, and MD Iseman. 2000 Immunological Diagnosis of Tuberculosis: a review. Tuber. Lung Dis. 80: 131-140.
- Foulds, J., and R. O'Brien. 1998. New tools for the diagnosis of tuberculosis: the perspective of developing countries. Int. J. Tubercle. Lung Dis. 2: 778-783

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14. Referências

1. Schaaf, H. S., P. Botha, N. Beyers, R. P. Gie, H. A. et, al 1996. The 5-year outcome of multidrug resistant tuberculosis patients in the Cape Province of South Africa. Trop. Med. Int. Health 1:718-722
2. Havlir, D. V., and P. F. Barnes. 1999. Tuberculosis in patients with human immunodeficiency virus infection. N. Engl. J. Med. 340:367-373
3. Dye L., Scheele S., V. Pathania, et al: 1999 Global Burden of Tuberculosis Estimated Incidence, Prevalence, and Mortality by Country. WHO Global Surveillance and Monitoring Project. JAMA. 282:677-686.
4. Kochi, A. 1991 The global tuberculosis situation and the new control strategy of the World Health Organization. Tubercle 72:1-6
5. Merlin TL, Gibson DW, and DH Connor 1994. Tuberculosis, p400-404 in Rubein E and Farber JL (ed)-Pathology, 2nd edition. J.B. Lippincott Company
6. Daniel, T. M. 1996. Immunodiagnosis of tuberculosis, p. 223-231. In W. N. Rom, and S. Garay (ed.), Tuberculosis. Letter, Brown & Co., Boston, Mass.
7. Wilkens, E. G. L 1994. The serodiagnosis of tuberculosis, p. 367-379. In P. D. O. Davies (ed.), Clinical tuberculosis. Chapman & Hall, Ltd., London, England.
8. Chan ED., Heifets L., and MD Iseman. 2000 Immunological Diagnosis of tuberculosis: a review. Tuberc. Lung Dis. 80: 131-140.
9. Foulds, J., and R. O'Brien. 1998. New tools for the diagnosis of tuberculosis: the perspective of developing countries. Int. J. Tubercle. Lung Dis. 2: 778-783

Rev.2, 2015-11-25 MNL

| Symbol | Deutsch | English | Français | Español | Italiano | Polski |
|--|-----------------------------|------------------------------------|---|--|---|---------------------------------|
|  | CE Konformitätszeichen | CE marking of conformity | Conformité aux normes européennes | Conformidad europea | Conformità europea | Znak zgodności CE |
|  | Gebrauchsanweisung beachten | Consult instructions for use | Consulter la notice d'utilisation | Consúltense las instrucciones de uso | Consultare le istruzioni per l'uso | Przestrzegać instrukcji obsługi |
|  | In-vitro-Diagnostika | In-vitro diagnostic medical device | Dispositif médical de diagnostic in vitro | Producto sanitario para diagnóstico in vitro | Dispositivo medico-diagnostico in vitro | Tylko do diagnostyki in vitro |
|  | Temperaturbegrenzung | Temperature limitation | Limites de température | Límite de temperatura | Limiti di temperatura | Temperatura przechowywania |
|  | Chargenbezeichnung | Batch code | Code du lot | Código de lote | Codice lotto | Numer serii |
|  | Nicht zur Wiederverwendung | Do not reuse | Né pas réutiliser | No reutilizar | Non riutilizzare | Tylko do jednorazowego użytku |
|  | Verwendbar bis | Use by | Utiliser jusqu'au | Fecha de caducidad | Utilizzare entro | Data ważności |
|  | Bestellnummer | Catalogue Number | Référence du catalogue | Número de catálogo | Riferimento di Catalogo | Numer katalogowy |
|  | Hersteller | Manufacturer | Fabricant | Fabricante | Fabbricante | Producent |
|  | Ausreichend für <n> Ansätze | Sufficient for <n> tests | Suffisant pour "n" tests | Suficiente para <n> utilizaciones | Sufficiente per "n" saggi | Wystarczający na <n> Powtórzeń |

| Symbol | Português | Český | Suomi | Svenskt | Nederlands | Dansk |
|--|--|---|---|---|---|--|
|  | Conformidade com as normas europeias | CE certifikát | CE-merkitty | CE-märkning | CE-markering | CE-mærkning |
|  | Consultar as instruções de utilização | Viz návod k použití | Katso käyttöohjeita | Läs bruksanvisningen | Raadpleeg de gebruiksaanwijzing | Se brugsanvisningen |
|  | Dispositivo médico para diagnóstico in vitro | Diagnostický zdravotnický prostředek in vitro | In vitro - diagnostiikkaan tarkoitettu lääkinnällinen laite | Medicinteknisk produkt avsedd för in vitro-diagnostik | Medisch hulpmiddel voor in-vitrodiagnostiek | Medicinsk udstyr til in vitro-diagnostik |
|  | Limites de temperatura | Teplotní omezení | Lämpötilarajat | Temperaturbegränsning | Temperatuurlimiet | Temperaturbegrænsning |
|  | Código do lote | Kód šarže | Eräkoodi | Satsnummer | Code van de partij | Batchkode |
|  | Não reutilizar | Pro jednorázové použití | Kertakäyttöinen | Får inte återvändas | Niet opnieuw gebruiken | Må ikke genbruges |
|  | Prazo de validade | Spotřebujte do | Käytettävä viimeistään | Används före | Houdbaar tot | Udløbsdato |
|  | Número de catálogo | Katalogov číslo | Luettelonumero | Listnummer | Catalogus nummer | Bestilingsnummer |
|  | Fabricante | Výrobce | Valmistaja | Tillverkare | Fabrikant | Fabrikant |
|  | Suficiente para <n> test | Dostačuje pro <n> testů | Lukumäärä <n> test | Räcker till <n> test | Volgende voor <n> test | Tilstrækkeligt til <n> test |

Our Teams

Germany:

Regensburg

Tel: +49 941 290 10-0

Fax: +49 941 290 10-50

Moers

Tel: +49 2841 99820-0

Fax: +49 2841 99820-1

Austria:

Tel: +49 941 290 10-29

Free Tel: 0800 291 565

Fax: +49 290 10-50

Free Fax: 0800 298 197

UK & Ireland:

Tel: +49 941 290 10-18

Free Tel –UK: 0808 234 1237

Free Tel – IRE: 1800 555 080

Fax: +49 290 10-50

France:

France Tel: 0800 915 240

France Fax: 0800 909 493

Switzerland

Swiss Tel: 0800 564 720

Swiss Fax: 0800 837 476

Belgium

Belgium Tel: 0800 718 82

Belgium Fax: 0800 747 07

Luxembourg

Lux. Tel: 800 211 16

Lux. Fax: 800 261 79

Spain:

Tel: +49 941 290 10-759

Free Tel: 900 938 315

Fax: +49 941 290 10-50

Free Fax: 900 984 992

Italy:

Tel: +49 941 290 10-34

Fax: +49 941 290 10-50

Poland:

Tel: +49 941 290 10-44

Free Tel: 00 800 491 15 95

Fax: +49 941 290 10-50

Free Fax: 00 800 491 15 94

Portugal:

Tel: +49 941 290 10-735

Tel. Verde: 800 849 230

Fax: +49 941 290 10-50

Fax Verde: 800 849 229

Netherlands:

Tel: +31 30 75 600

Free Tel: 0800 0222 890

Fax: +31 70 30 30 775

Free Fax: 0800 024 9519

Denmark:

Tel: +31 703075 603

Free Tel: +45 80 88 87 53

Tax: +31 703030 775

Laboratory Diagnostics Team:

Tel: +49 941 290 10-40

Fax: +49 941 290 10-50

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Carl-Zeiss-Strasse 12 • 47445 Moers • Germany

www.nal-vonminden.com • info@nal-vonminden.com

Fon: +49 2841 99820-0 • Fax: +49 2841 99820-1

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nal von minden GmbH