NADAL[®] Ebola Test (test cassette)

REF 2210001

DE Gebrauchsanweisung	2	SE Användarinstruktioner	21
EN Instruction for use	6	Symbols	27
FR Instructions d'utilisation	9	Our Teams	28
ES Instrucciones de uso	12		
IT Istruzioni per l'uso	15		
PL Sposób użycia	18		

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Spezifität: 98%

14. Referenzen

1. Lucht A, Grunow R, Möller P, Feldmann H, Becker S: Development, characterization and use of monoclonal antibodies for the detection of Ebola virus. *J Virol Methods* 2003; 111:21-8
2. Grolla A, Lucht A, Dick D, Strong JE, Feldmann H: Laboratory diagnosis of Ebola and Marburg hemorrhagic fever. *Bull Soc Pathol Exot* 2005; 98:205-9
3. Becker S, Feldmann H, Will C, Slenczka W: Evidence of occurrence of filovirus antibodies in humans and imported monkeys: Do subclinical filovirus infections occur worldwide? *Med Microbiol Immunol* 1992; 181:43-55
4. World Health Organization. Outbreak(s) of Ebola hemorrhagic fever in the Republic of Congo, January-April 2003. *Wkly Epidemiol Rec* 2003; 78:285-9
5. Hartmann H, Lübbers B, Cararetto M, Bautsch W, Klos A, Köhl J. Rapid quantification of C3a and C5a using a combination of chromatographic and immunoassay procedures. *J Immunol Methods* 1993; 166:35-44

6. Ksiazek TG, Rollin PE, Jahrling PB, Johnson E, Dalgard DW, Peters CJ: Enzyme immunosorbent assay for Ebola virus antigens in tissues of infected primates. J Clin Microbiol 1992; 30:947-50
7. National Committee for Clinical Standards Clinical Waste Management: Approved Guideline. NCCLS Document GPS-A. Villanova, PA: NCCLS 1993; 13:1-18, 29-42
8. Townner JS, Rollin PE, Bausch DG, Sanchez A, Crary SM, Vincent M, Lee WF, Spiropoulou CF, Ksiazek TG, Lukwiya M, Kaducu F, Downing R, Nichol ST: Rapid diagnosis of ebola hemorrhagic fever by reverse transcription-PCR in an outbreak setting and assessment of patient viral load as a predictor of outcome. J Virol 2004; 78:4330-41

Rev. 0, 2015-05-26 OM/UJ

1. Intended Use

The NADAL® Ebola Test (serum/throat swab) is a rapid chromatographic immunoassay for the qualitative detection of *Ebolavirus* antigen in human serum or throat swabs. The test is intended to be used as an aid to detect Ebola infection in febrile conditions in conjunction with clinical information/methods and other diagnostic *in-vitro* procedures.

2. Introduction and Clinical Significance

Ebola virus disease (EVD) also known as Ebola haemorrhagic fever (EHF) is a severe, often fatal disease in humans and non-human primates (monkeys, gorillas and chimpanzees) which has appeared sporadically. The genus *Ebolavirus* is one of the three members of the *Filoviridae* family of RNA viruses. There are five identified species in the genus *Ebolavirus*: *Zaire ebolavirus* (EBOV), *Sudan ebolavirus* (SUDV), *Reston ebolavirus* (RESTV), *Bundibugyo ebolavirus* (BDBV) and *Tai Forest ebolavirus* (TAFV). Up to 90% of the outbreaks are caused by *Zaire ebolavirus*.

3. Test Principle

The NADAL® Ebola Test is a rapid, qualitative, membrane-based immunoassay for the detection of VP40 antigen of the Ebola virus in human serum or throat swabs. After applying the sample to the sample well, it reacts with anti-Ebola antibodies coated to particles (anti-Ebola antibody-particle complexes). The mixture migrates along the membrane by capillary action.

Upon reaching the test line region, antigens bound to the anti-Ebola antibody-particle complexes react with the capture reagent, which are anti-Ebola antibodies immobilised in the test line region of the membrane, and generate a coloured line.

The appearance of a coloured line in the test line region indicates a positive result (the presence of Ebola antigen), while its absence indicates a negative result.

The appearance of a coloured line in the control line region serves as a procedural control, indicating that the proper volume of specimen has been added and membrane wicking has occurred.

4. Reagents and Materials Supplied

- 20 NADAL® Ebola test cassettes (containing particles coated with anti-Ebola antibodies, anti-Ebola antibodies and anti-mouse IgG immobilised on nitrocellulose membrane).
- 20 transfer pipettes
- 20 extraction tubes
- Provided additional material according to 93/42/EEC:
 - 20 sterile polyester-tipped swabs CE 0086
 - Puritan Medical Products Company LLC
 - 31 School Street
 - Guilford, Maine 04443-0149 USA (authorised EU representative EMERGO EUROPE, The Hague, The Netherlands)
- 1 extraction buffer Ebola (10 ml)
- 1 tube rack
- 1 package insert

5. Additional Materials Required

- Timer
- Sample collection device for serum

6. Storage & Stability

Store the test unopened at 4-30°C and use it until the expiration date printed on the sealed foil pouch. The test cassette must remain in the sealed foil pouch until use. Do not use any component of the test kit beyond the labeled expiration date. Do not freeze the test kit.

Test cassettes are sensitive to light and humidity. Therefore, do not use test cassettes that have become wet or have damaged pouches.

After opening the sealed foil pouch, the test should be performed immediately or maximum within 1 hour.

After opening the buffer vial, the extraction buffer can be stored at 4-30°C and used safely for 3 months.

7. Warnings and Precautions

- For professional *in-vitro* diagnostic use only.
- Carefully read through the test procedure prior to testing.
- Do not use the test after the expiration date indicated on the package.
- Do not use the test if the foil pouch is damaged.
- Do not reuse tests.
- Do not add samples to the reaction area (result area).
- In order to avoid contamination, do not touch the reaction area (result area).
- Avoid cross-contamination of specimens by using a new extraction tube for each specimen obtained.
- Do not substitute or mix components from different test kits.
- Do not eat, drink or smoke in the area where specimens and test kits are handled.
- Wear protective clothing such as laboratory coats, disposable gloves and eye protection when specimens are being assayed.
- Handle all specimens as if they contain infectious agents. Observe established precautions for microbiological risks throughout all procedures and standard guidelines for the appropriate disposal of specimens.
- The test kit contains products of animal origin. Certified knowledge of the origin and/or sanitary state of the animals does not completely guarantee the absence of transmissible pathogenic agents. It is therefore recommended that these products be treated as potentially infectious, and handled in accordance with usual safety precautions (e.g., do not ingest or inhale).
- Do not use swabs from damaged pouches.
- Humidity and temperature can adversely affect test results.
- Used testing materials should be decontaminated and discarded according to local regulations.

8. Specimen Collection and Preparation

Samples should be collected, handled and shipped following appropriate infection control precautions for Ebola, other hemorrhagic fever viruses or other infectious diseases.

Shipping can be performed according to the policies of the sender, customs regulations and the requirements of the

receiving laboratory. If samples are to be shipped, they should be packed in compliance with statutory regulations for transportation of etiologic agents.

Serum

Use only human serum collected and prepared according to the instructions of the manufacturer of the sampling device for sample collection.

Use only clear, non-hemolysed samples.

Serum can be stored at 2-8°C for up to 3 days. For long-term storage, samples should be kept below -20°C.

Bring samples to room temperature prior to testing. Frozen samples must be completely thawed and well mixed prior to testing.

Samples should not be frozen and thawed repeatedly.

Throat swabs

Use samples at room temperature (15-30°C) only.



Use only fresh throat swabs. Do not store or freeze throat swabs.

9. Test Procedure


Bring test components and samples to room temperature (15-30°C) prior to testing.

- Remove the test cassette from the sealed foil pouch and use it as soon as possible. The test should be performed within 1 hour of opening the sealed foil pouch.
- Place the test cassette on a clean and level surface.
- Label the test cassette with the patient or control identification.
- Label the extraction tube with the patient or control identification and place it into the tube rack.

A. Serum

- Holding the extraction buffer vial vertically, add 3 drops (approx. 100 µl) of the extraction buffer to the extraction tube. Use only the buffer provided with the test kit. 
- Holding the provided transfer pipette vertically, transfer 3 drops (approx. 100 µl) of serum sample into the extraction tube and incubate the solution for 1 to 15 minutes. 
- Now the mixture is ready to be tested. Fill the provided transfer pipette with the mixture. Holding the pipette vertically over the sample well (S) of the test cassette, drop the whole mixture into the sample well. Wait until each drop is absorbed before adding the next one. Start the timer. 
- Read the test result after 10 minutes. Do not interpret the test result after more than 15 minutes. 

B. Throat swabs

- Hold the extraction buffer vial vertically and add 9 drops of the extraction buffer to the extraction tube. Use only the buffer provided with this test kit. 
- Insert the provided sterile swab into the patient's oral cavity. Very gently push the swab until resistance is met at the level of the pharynx.
- Gently rotate the swab a few times against the mucous membrane.
- Immediately insert the swab with collected throat swab sample into the extraction tube, compress the bottom of the tube and rotate the swab vigorously to mix the reagents for at least 1 minute. Press the swab against the inside of the tube and squeeze as much liquid as possible from the swab by rotating and pressing its fiber part against the wall of the tube. Discard the swab. 
- Let the mixture incubate for 1 to 15 minutes. 
- Mix the contents of the tube again by swirling it gently.
- Now the mixture is ready to be tested. Fill the provided transfer pipette with the mixture. Holding the provided pipette over the sample well (S) of the test cassette, squeeze 6 drops into the sample well (S). Wait until each drop is absorbed before adding the next one. Start the timer. 
- Read the test result after 10 minutes. Do not interpret the test result after more than 15 minutes. 

10. Result Interpretation

Positive:

Two distinct coloured lines appear. One line appears in the control line region (C) and the other line appears in the test line region (T).

Any visible coloured line, appearing in the test line region (T), should be interpreted as positive, even if the test line is lighter or darker than the line in the control line region (C).

Negative:

One coloured line appears in the control line region (C). No visible line appears in the test line region (T).

Invalid:

The control line fails to appear. Results from any test which has not produced a control line at the specified reading time must be discarded.

Insufficient specimen volume, incorrect operating procedure or expired tests are the most likely reasons for control line

failure. Please review the procedure and repeat the test with a new test cassette. If the problem persists, discontinue using the test kit immediately and contact your distributor.

11. Quality Control

An internal procedural control is included in the test cassette:

A coloured line appearing in the control line region (C) is considered an internal procedural control. It confirms sufficient specimen volume, adequate membrane wicking and correct procedural technique.

12. Limitations

- The NADAL® Ebola Test is designed for qualitative *in-vitro* detection of the VP40 antigen of the Ebola virus in human serum and throat swabs. Other body fluids and pooled samples may give inaccurate results and should not be used.
- The NADAL® Ebola Test should only be used with samples from patients with suspected Ebola infections (febrile diseases) as an aid to detect Ebola haemorrhagic fever. It should not be used as the sole criterion in diagnosing Ebola infections.
- The NADAL® Ebola Test detects the presence of Ebola virions at densities of $\geq 10^6$ PFU/ml in the sample. Therefore, it should not be used to exclude Ebola infections.
- Results should be interpreted by a physician only.
- Results should be used only in conjunction with clinical information/symptoms/methods and other *in-vitro* diagnostic procedures (e.g. PCR techniques).
- The intensity of the coloured line in the test line region does not correlate with the antigen titer in the sample.
- A negative result does not exclude the possibility of Ebola infection at any time.
- A positive result should be confirmed by using a different method and should be evaluated together with the overall clinical evaluation.
- False-positive or false-negative results may be produced by:
 - insufficient sample volume or incorrect procedural techniques
 - the use of inappropriate samples or materials (different buffer as the one provided with the test kit, other volumes of samples and buffer than those indicated)
 - disregard for the incubation time of at least 1 minute of sample and extraction buffer
 - disregard for the waiting time of 10 minutes after applying the mixture to the test cassette
- Low titers of Ebola may result in a faint line appearing in the test line region after prolonged time.
- The "Hook effect" in case of very high antigen levels in the sample can lead to false-negative results.
- Apart from EDTA, no substance has been demonstrated to trigger cross-reactions.

13. Performance Characteristics

The performance of the NADAL® Ebola Test has been compared with that of the PCR system RealStar® Ebola virus RT-PCR Kit 1.0 based on samples from African patients from Conakry in Guinea.

		NADAL® Ebola Test		
		+	--	Total
PCR	+	12	1	13
	--	2	103	105
	Total	14	104	118

Analytic sensitivity: \leq CT 23

Specificity: 98 %

14. References

1. Lucht A, Grunow R, Möller P, Feldmann H, Becker S: Development, characterization and use of monoclonal antibodies for the detection of Ebola virus. *J Virol Methods* 2003; 111:21-8
2. Grolla A, Lucht A, Dick D, Strong JE, Feldmann H: Laboratory diagnosis of Ebola and Marburg hemorrhagic fever. *Bull Soc Pathol Exot* 2005; 98:205-9
3. Becker S, Feldmann H, Will C, Slenczka W: Evidence of occurrence of filovirus antibodies in humans and imported monkeys: Do subclinical filovirus infections occur worldwide? *Med Microbiol Immunol* 1992; 181:43-55
4. World Health Organization. Outbreak(s) of Ebola haemorrhagic fever in the Republic of Congo, January-April 2003. *Wkly Epidemiol Rec* 2003; 78:285-9
5. Hartmann H, Lübbers B, Cararetto M, Bausch W, Klos A, Köhl J. Rapid quantification of C3a and C5a using a combination of chromatographic and immunoassay procedures. *J Immunol Methods* 1993; 166:35-44
6. Ksiazek TG, Rollin PE, Jahrling PB, Johnson E, Dalgard DW, Peters CJ: Enzyme immunosorbent assay for Ebola virus antigens in tissues of infected primates. *J Clin Microbiol* 1992; 30:947-50
7. National Committee for Clinical Standards Clinical Waste Management: Approved Guideline. NCCLS Document GP5-A. Villanova, PA: NCCLS 1993; 13:1-18, 29-42
8. Townner JS, Rollin PE, Bausch DG, Sanchez A, Crary SM, Vincent M, Lee WF, Spiropoulou CF, Ksiazek TG, Lukwya M, Kaducu F, Downing R, Nichol ST: Rapid diagnosis of ebola hemorrhagic fever by reverse transcription-PCR in an outbreak setting and assessment of patient viral load as a predictor of outcome. *J Virol* 2004; 78:4330-41

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Spécificité : 98%

14. Bibliographie

1. Lucht A, Grunow R, Möller P, Feldmann H, Becker S: Development, characterization and use of monoclonal antibodies for the detection of Ebola virus. *J Virol Methods* 2003; 111:21-8
2. Grolla A, Lucht A, Dick D, Strong JE, Feldmann H: Laboratory diagnosis of Ebola and Marburg hemorrhagic fever. *Bull Soc Pathol Exot* 2005; 98:205-9
3. Becker S, Feldmann H, Will C, Slenczka W: Evidence of occurrence of filovirus antibodies in humans and imported monkeys: Do subclinical filovirus infections occur worldwide? *Med Microbiol Immunol* 1992; 181:43-55
4. World Health Organization. Outbreak(s) of Ebola haemorrhagic fever in the Republic of Congo, January-April 2003. *Wkly Epidemiol Rec* 2003; 78:285-9
5. Hartmann H, Lübbers B, Cararetto M, Bausch W, Klos A, Köhl J. Rapid quantification of C3a and C5a using a combination of chromatographic and immunoassay procedures. *J Immunol Methods* 1993; 166:35-44
6. Ksiazek TG, Rollin PE, Jahrling PB, Johnson E, Dalgard DW, Peters CJ: Enzyme immunoassay for Ebola virus antigens in tissues of infected primates. *J Clin Microbiol* 1992; 30:947-50
7. National Committee for Clinical Standards Clinical Waste Management: Approved Guideline. *NCCLS Document GPS-A. Villanova, PA: NCCLS* 1993; 13:1-18, 29-42
8. Townner JS, Rollin PE, Bausch DG, Sanchez A, Crary SM, Vincent M, Lee WF, Spiropoulou CF, Ksiazek TG, Lukwiyi M, Kaducu F, Downing R, Nichol ST: Rapid diagnosis of ebola hemorrhagic fever by reverse transcription-PCR in an outbreak setting and assessment of patient viral load as a predictor of outcome. *J Virol* 2004; 78:4330-41

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Sensibilidad analítica: \leq CT 23

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14. Referencias

1. Lucht A, Grunow R, Möller P, Feldmann H, Becker S: Development, characterization and use of monoclonal antibodies for the detection of Ebola virus. *J Virol Methods* 2003; 111:21-8
2. Grolla A, Lucht A, Dick D, Strong JE, Feldmann H: Laboratory diagnosis of Ebola and Marburg hemorrhagic fever. *Bull Soc Pathol Exot* 2005; 98:205-9
3. Becker S, Feldmann H, Will C, Slenczka W: Evidence of occurrence of filovirus antibodies in humans and imported monkeys: Do subclinical filovirus infections occur worldwide? *Med Microbiol Immunol* 1992; 181:43-55
4. World Health Organization. Outbreak(s) of Ebola haemorrhagic fever in the Republic of Congo, January-April 2003. *Wkly Epidemiol Rec* 2003; 78:285-9
5. Hartmann H Lübbers B, Cararetto M, Bautsch W, Klos A, Köhl J. Rapid quantification of C3a and C5a using a combination of chromatographic and immunoassay procedures. *J Immunol Methods* 1993; 166:35-44
6. Ksiazek TG, Rollin PE, Jahrling PB, Johnson E, Dalgard DW, Peters CJ: Enzyme immunosorbent assay for Ebola virus antigens in tissues of infected primates. *J Clin Microbiol* 1992; 30:947-50
7. National Committee for Clinical Standards Clinical Waste Management: Approved Guideline. NCCLS Document GP5-A. Villanova, PA: NCCLS 1993; 13:1-18, 29-42
8. Townner JS, Rollin PE, Bausch DG, Sanchez A, Crary SM, Vincent M, Lee WF, Spiropoulou CF, Ksiazek TG, Lukwya M, Kaducu F, Downing R, Nichol ST: Rapid diagnosis of ebola hemorrhagic fever by reverse transcription-PCR in an outbreak setting and assessment of patient viral load as a predictor of outcome. *J Virol* 2004; 78:4330-41

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Puritan Medical Products Company LLC
31 School Street

Guilford, Maine 04443-0149 USA (authorised EU representative EMERGO EUROPE, The Hague, The Netherlands)

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Specificità: 98 %

14. Bibliografia

1. Lucht A, Grunow R, Möller P, Feldmann H, Becker S: Development, characterization and use of monoclonal antibodies for the detection of Ebola virus. *J Virol Methods* 2003; 111:21-8
2. Grolla A, Lucht A, Dick D, Strong JE, Feldmann H: Laboratory diagnosis of Ebola and Marburg hemorrhagic fever. *Bull Soc Pathol Exot* 2005; 98:205-9
3. Becker S, Feldmann H, Will C, Slenczka W: Evidence of occurrence of filovirus antibodies in humans and imported monkeys: Do subclinical filovirus infections occur worldwide? *Med Microbiol Immunol* 1992; 181:43-55
4. World Health Organization. Outbreak(s) of Ebola haemorrhagic fever in the Republic of Congo, January-April 2003. *Wkly Epidemiol Rec* 2003; 78:285-9
5. Hartmann H, Lübbers B, Cararetto M, Bautsch W, Klos A, Köhl J. Rapid quantification of C3a and C5a using a combination of chromatographic and immunoassay procedures. *J Immunol Methods* 1993; 166:35-44
6. Kisazek TG, Rollin PE, Jahrling PB, Johnson E, Dalgaard DW, Peters CJ: Enzyme immunoassay for Ebola virus antigens in tissues of infected primates. *J Clin Microbiol* 1992; 30:947-50
7. National Committee for Clinical Standards Clinical Waste Management: Approved Guideline. NCCLS Document GP5-A. Villanova, PA: NCCLS 1993; 13:1-18, 29-42
8. Townser JS, Rollin PE, Bausch DG, Sanchez A, Crary SM, Vincent M, Lee WF, Spiropoulou CF, Kisazek TG, Lukwija M, Kaducu F, Downing R, Nichol ST: Rapid diagnosis of ebola hemorrhagic fever by reverse transcription-PCR in an outbreak setting and assessment of patient viral load as a predictor of outcome. *J Virol* 2004; 78:4330-41

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Puritan Medical Products Company LLC
31 School Street
Guilford, Maine 04443-0149 USA (upoważniony reprezentant w UE EMERGO EUROPE, The Hague, The Netherlands)

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Swoistość: 98%

14. Bibliografia

1. Lucht A, Grunow R, Möller P, Feldmann H, Becker S: Development, characterization and use of monoclonal antibodies for the detection of Ebola virus. *J Virol Methods* 2003; 111:21-8
2. Grolla A, Lucht A, Dick D, Strong JE, Feldmann H: Laboratory diagnosis of Ebola and Marburg hemorrhagic fever. *Bull Soc Pathol Exot* 2005; 98:205-9
3. Becker S, Feldmann H, Will C, Slenczka W: Evidence of occurrence of filovirus antibodies in humans and imported monkeys: Do subclinical filovirus infections occur worldwide? *Med Microbiol Immunol* 1992; 181:43-55
4. World Health Organization. Outbreak(s) of Ebola haemorrhagic fever in the Republic of Congo, January-April 2003. *Wkly Epidemiol Rec* 2003; 78:285-9
5. Hartmann H, Lübbers B, Cararetto M, Bausch W, Klos A, Köhl J: Rapid quantification of C3a and C5a using a combination of chromatographic and immunoassay procedures. *J Immunol Methods* 1993; 166:35-44
6. Ksiazek TG, Rollin PE, Jahrling PB, Johnson E, Dalgaard DW, Peters CJ: Enzyme immunoassay for Ebola virus antigens in tissues of infected primates. *J Clin Microbiol* 1992; 30:947-50
7. National Committee for Clinical Standards Clinical Waste Management: Approved Guideline. NCCLS Document GPS-A. Villanova, PA: NCCLS 1993; 13:1-18, 29-42
8. Townner JS, Rollin PE, Bausch DG, Sanchez A, Crary SM, Vincent M, Lee WF, Spiropoulou CF, Ksiazek TG, Lukwima M, Kaducu F, Downing R, Nichol ST: Rapid diagnosis of Ebola hemorrhagic fever by reverse transcription-PCR in an outbreak setting and assessment of patient viral load as a predictor of outcome. *J Virol* 2004; 78:4330-41

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Specifitet: 98 %

14. Referenser

- Lucht A, Grunow R, Möller P, Feldmann H, Becker S: Development, characterization and use of monoclonal antibodies for the detection of Ebola virus. *J Virol Methods* 2003; 111:21-8
- Grolla A, Lucht A, Dick D, Strong JE, Feldmann H: Laboratory diagnosis of Ebola and Marburg hemorrhagic fever. *Bull Soc Pathol Exot* 2005; 98:205-9
- Becker S, Feldmann H, Will C, Slenczka W: Evidence of occurrence of filovirus antibodies in humans and imported monkeys: Do subclinical filovirus infections occur worldwide? *Med Microbiol Immunol* 1992; 181:43-55
- World Health Organization. Outbreak(s) of Ebola haemorrhagic fever in the Republic of Congo, January-April 2003. *Wkly Epidemiol Rec* 2003; 78:285-9
- Hartmann H, Lübbers B, Cararetto M, Bausch W, Klos A, Köhl J. Rapid quantification of C3a and C5a using a combination of chromatographic and immunoassay procedures. *J Immunol Methods* 1993; 166:35-44
- Ksiazek TG, Rollin PE, Jahrling PB, Johnson E, Dalgard DW, Peters CJ: Enzyme immunoassay for Ebola virus antigens in tissues of infected primates. *J Clin Microbiol* 1992; 30:947-50
- National Committee for Clinical Standards Clinical Waste Management: Approved Guideline. NCCLS Document GPS-A. Villanova, PA: NCCLS 1993; 13:1-18, 29-42
- Towner JS, Rollin PE, Bausch DG, Sanchez A, Crary SM, Vincent M, Lee WF, Spiropoulou CF, Ksiazek TG, Lukwya M, Kaducu F, Downing R, Nichol ST: Rapid diagnosis of ebola hemorrhagic fever by reverse transcription-PCR in an outbreak setting and assessment of patient viral load as a predictor of outcome. *J Virol* 2004; 78:4330-41

Rev. 0, 2015-05-26 DL

Symbol	Deutsch	English	Français	Español	Italiano	Polski
	CE Konformitätszeichen	CE marking of conformity	Conformité aux normes européennes	Conformidad europea	Conformità europea	Znak zgodności CE
	Gebrauchsanweisung beachten	Consult instructions for use	Consulter la notice d'utilisation	Consúltense las instrucciones de uso	Consultare le istruzioni per l'uso	Przestrzegać instrukcji obsługi
	In-vitro-Diagnostika	In-vitro diagnostic medical device	Dispositif médical de diagnosticc in vitro	Producto sanitario para diagnóstico in vitro	Dispositivo medico-diagnostico in vitro	Tylko do diagnostyki in vitro
	Temperaturbegrenzung	Temperature limitation	Limites de température	Limitación de temperatura	Limiti di temperatura	Temperatura przechowywania
	Chargenbezeichnung	Batch code	Code du lot	Código de lote	Codice lotto	Numer serii
	Nicht zur Wiederverwendung	Do not reuse	Ne pas réutiliser	No reutilizar	Non riutilizzare	Tylko do jednorazowego użytku
	Verwendbar bis	Use by	Utiliser jusqu'au	Fecha de caducidad	Utilizzare entro	Data ważności
	Bestellnummer	Catalogue Number	Référence du catalogue	Número de catálogo	Riferimento di Catalogo	Numer katalogowy
	Hersteller	Manufacturer	Fabricant	Fabricante	Fabbricante	Producent
	Ausreichend für <n> Ansätze	Sufficient for <n> tests	Suffisant pour pour "n" tests	Válido para para <n> ensayos	Sufficiente per "n" saggi	Wystarczający na <n> Powtórzeń

Symbol	Português	Český	Suomi	Svenskt	Nederlands	Dansk
	Conformidade com as normas europeias	CE certifikát	CE-merkitty	CE-märkning	CE-markering	CE-mærkning
	Consultar as instruções de utilização	Viz návod k použití	Katso käyttöohjetta	Läs bruksanvisningen	Raadpleeg de gebruiksaanwijzing	Se brugsanvisningen
	Dispositivo médico para diagnóstico in vitro	Diagnostický zdravotnický prostředek in vitro	In vitro - diagnostiikkaan tarkoitettu lääkinällinen laite	Medicinteknisk produkt avsedd för in vitro-diagnostik	Medisch hulpmiddel voor in-vitrodiagnostiek	Medicinsk udstyr til in vitro-diagnostik
	Limites de temperatura	Teplotní omezení	Lämpötilarajat	Temperaturbegränsning	Temperatuurlimiet	Temperaturbegrænsning
	Código do lote	Kód šarže	Eräkoodi	Satsnummer	Code van de partij	Batchkode
	Não reutilizar	Pro jednorázové použití	Kertakäyttöinen	Får inte återvändas	Niet opnieuw gebruiken	Må ikke genbruges
	Prazo de validade	Spotřebuje do	Käytettävä viimeistään	Används före	Houdbaar tot	Udløbsdato
	Número de catálogo	Katalogov číslo	Luettelonumero	Listnummer	Catalogus nummer	Bestilingsnummer
	Fabricante	Výrobce	Valmistaja	Tillverkare	Fabrikant	Fabrikant
	Suficiente para <n> test	Dostačuje pro <n> testů	Lukumäärä <n> test	Räcker till <n> test	Voldoende voor <n> test	Tilstrækkeligt til <n> test

Our Teams

Germany:

Regensburg

Tel: +49 941 290 10-0
Fax: +49 941 290 10-50

Moers

Tel: +49 2841 99820-0
Fax: +49 2841 99820-1

Austria:

Tel: +49 941 290 10-29
Free Tel: 0800 291 565
Fax: +49 290 10-50
Free Fax: 0800 298 197

UK & Ireland:

Tel: +49 941 290 10-18
Free Tel –UK: 0808 234 1237
Free Tel – IRE: 1800 555 080
Fax: +49 290 10-50

France:

France Tel: 0800 915 240
France Fax: 0800 909 493

Switzerland

Swiss Tel: 0800 564 720
Swiss Fax: 0800 837 476

Belgium

Belgium Tel: 0800 718 82
Belgium Fax: 0800 747 07

Luxembourg

Lux. Tel: 800 211 16
Lux. Fax: 800 261 79

Spain:

Tel: +49 941 290 10-759
Free Tel: 900 938 315
Fax: +49 941 290 10-50
Free Fax: 900 984 992

Italy:

Tel: +49 941 290 10-34
Fax: +49 941 290 10-50

Poland:

Tel: +49 941 290 10-44
Free Tel: 00 800 491 15 95
Fax: +49 941 290 10-50
Free Fax: 00 800 491 15 94

Portugal:

Tel: +49 941 290 10-735
Tel. Verde: 800 849 230
Fax: +49 941 290 10-50
Fax Verde: 800 849 229

Netherlands:

Tel: +31 30 75 600
Free Tel: 0800 0222 890
Fax: +31 70 30 30 775
Free Fax 0800 024 9519

Denmark:

Tel: +31 703075 603
Free Tel: +45 80 88 87 53
Tax: +31 703030 775

Laboratory Diagnostics Team:

Tel: +49 941 290 10-40
Fax: +49 941 290 10-50

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Carl-Zeiss-Strasse 12 • 47445 Moers • Germany

www.nal-vonminden.com • info@nal-vonminden.com

Fon: +49 2841 99820-0 • Fax: +49 2841 99820-1

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nal von minden GmbH