**Gebrauchsanweisung**  
**Instructions for use**  
**Instructions d'utilisation**  
**Instrucciones de uso**  
**Istruzioni per l'uso**  
**Instrukcja obsługi**

# **NADAL<sup>®</sup> Legionella**

## **Urine Antigen Test**

**(test cassette)**

# **REF 552005**

Version 1.0

2010-06-04

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### 14. Referenzen

- Edelstein, P.H. and R. D. Meyer. "Legionella Pneumonias." In Respiratory infections: Diagnosis and Management, 3rd Edition, James E. Pennington, editor. New York: Raven Press, Ltd. 1994
- Marston, B.J., H.B. Lipman, R.F. Breiman. Surveillance for legionnaires' Disease: risk factors for morbidity and mortality. Arch.Intern. Med 1994;154:2417:2422
- Horwitz, M.A., B.J. Martson, C.V. Broome, and R.F. Freiman. Prospects for vaccine development. Presented at the 4th International Symposium on Legionella, 1992. In: Barbaree, J. M., R.F. Breiman, and A.P. DuFour, eds. Legionella: Current status and emerging perspectives. Washington, DC American Society for Microbiology, 1993.
- Kohler, R.B., W.C. Winn, Jr. and L.J. Wheat. Onset and duration of urinary antigen excretion in Legionnaires' disease. J. Clin. Microbiol. 1984; 20:605-607.
- Lindsay, D.S.J., W.H. Abraham, P. Christie, F. Johnston and G.F.S. Edwards. Laboratory diagnosis of legionnaires' disease due to *Legionella pneumophila* serogroup 1: comparison of phenotypic and genotypic methods. J. Med. Microbiol. 2004; 53:183-187.
- Cloud, J.L., K.C. Carroll, P. Pixton, M. Erali, and D.R. Hillyard. Detection of *Legionella* species in respiratory specimens using PCR with sequencing confirmation. J. Clin. Microbiol. 38; 2000; 1709-1712.
- Kohler, R.B., Winn, W.C., Wheat, L.J. Onset of the duration of urinary antigen excretion in Legionnaires' disease. J. Clin. Microbiol. 1984; 20: 601-605.
- Jonas, Daniel, A. Rosenbaum, S. Weyrick and S. Bhakdi. Enzyme-Linked Immunoassay for Detection of PCR-Amplified DNA and Legionellae in Bronchoalveolar fluid. J. Clin. Microbiol. 33; 1995; 1247-1252.

Rev.: 2010-02-14 SW; 2010-06-04 SW

## 1. Intended Use

The NADAL<sup>®</sup> Legionella Urine Antigen Test is a visually read, *in vitro* immunochromatographic rapid assay for the presumptive qualitative detection of *Legionella pneumophila* serogroup 1 antigens in human urine. This test is intended to aid in the presumptive diagnosis of Legionnaires' disease in conjunction with culture and other methods for patients with signs and symptoms of pneumonia. This rapid test is for professional use only.

## 2. Introduction and Clinical Significances

*Legionella pneumophila* serogroup 1 is gram-negative bacillus that appears off-white in color and circular in shape. This bacterium is associated with Pontiac fever and Legionnaires' disease. An estimated 25,000 to 100,000 cases involving *Legionella pneumophila* infections occur in the United States annually. About 5% to 30% of people with Legionnaires' disease die.

Infections caused by this bacterium include the pulmonary disease of pneumonia. In addition, extra pulmonary infections may be spread through the bloodstream or the lymph system, infecting the heart, brain, kidney, liver, and/or spleen. Early onset symptoms include fever, chills, muscle aches, and lack of appetite. Respiratory symptoms such as coughing or congestion are usually absent. As the disease progresses, a dry, hacking cough develops and may become productive after a few days. Half of people who develop Legionnaires' disease suffer shortness of breath. In many cases, a fever as high as 40°C (104°F) has been reported.

Common methods for the detection of *Legionella pneumophila* are culture, direct fluorescent antibody (DFA), indirect fluorescent antibody (IFA), enzyme immunoassays (EIA) and polymerase chain reaction (PCR) using respiratory specimens of expectorated sputum, bronchial washing, transtracheal aspirate and lung biopsy. *Legionella pneumophila* serogroup 1 antigen has been detected in urine during the acute phase of the disease and has become the method by choice for laboratory diagnosis. The NADAL<sup>®</sup> Legionella Urine Antigen Test is designed for the detection of *Legionella pneumophila* serogroup 1 soluble antigen in a urine sample.

## 3. Principle of the Test

The NADAL<sup>®</sup> Legionella Urine Antigen Test utilises a combination of polyclonal antibodies against the antigens of *Legionella pneumophila*.

The NADAL<sup>®</sup> Legionella Urine Antigen Test begins with the addition of urine to the test device. The specimen is absorbed by the sample pad and then moves through the conjugate pad which contains dried gold conjugated antibodies which are specific for *Legionella pneumophila* antigens; if the Legionella antigens are present in the urine sample, a "half-sandwich" immunocomplex is formed. This immuno-complex then migrates via capillary action along a nitrocellulose membrane containing immobilised antibodies to *Legionella pneumophila* antigens. The immobilised antibodies bind the "half-sandwich" immuno-complex to form a "whole sandwich" immuno-complex. Thus, when the "whole sandwich" is formed, a visible, pink coloured line develops in the specimen zone on the test device. In the absence of a Legionella antigen, a "sandwich" immuno-complex is not formed and a negative result is indicated. To serve as a procedural control, a pink coloured control line will always appear in the control zone regardless of the presence or absence of Legionella antigen. The test is available in cassette format.

## 4. Reagents and Materials Supplied

- 10 NADAL<sup>®</sup> Legionella Urine Antigen Test Devices incl. disposable specimen dropper  
(A membrane coated with a combination of polyclonal antibodies against the antigens of *Legionella pneumophila* serogroup 1)
- 1 User Manual

## 5. Additional Required Materials

- 1 Timer
- Specimen collection containers

## 6. Storage & Stability

NADAL<sup>®</sup> Legionella Urine Antigen Test devices should be kept at room temperature (15°-30°C) in sealed pouches with a desiccant. Do not freeze the test kit.

## 7. Warnings and Precautions

- For *in vitro* diagnostic use only.
- Use pipette provided with kit only.
- Use separate pipettes for each sample or control.
- The test device should remain in the pouch until ready for use.
- All specimens should be considered potentially hazardous, and should be handled in the same manner as an infectious agent.
- Wear disposable gloves while handling samples and wash hands after the assay is complete. Warning: The user should refer to the relevant section of the CDC-NIH manual "Biosafety in Microbiology and Biomedical Laboratories," Fifth Edition, 2007.
- Avoid any contact with eyes.
- The test device and all materials should be discarded in a proper biohazard container after testing.
- Do not use kits after the stated expiration date.
- Do not reuse test devices or kit materials.

## 8. Specimen Collection and Preparation


The urine specimen must be collected in a clean, sterile container. Urine specimens may be refrigerated (2°-8°C) and stored up to 72 hours prior to assay. Specimens may be stored frozen for up to 2 years. If specimens are refrigerated, they must be equilibrated to room temperature before testing.

## 9. Procedure of the Test

Allow specimen and/or controls to reach room temperature (15° – 30°C) prior to testing.

Remove test device from pouch and lay on a flat surface. Label the device with the patient ID.

1. Dispense **4 drops** (approximately 0.15 ml) of specimen into the round sample well (see illustration below). Wait for coloured lines to appear.


2. Read results at **15 minutes** or less. The presence of the control line is not indicative of the test being completed. Interpretation after 15 minutes may be inaccurate.

## 10. Interpretation of the Results


### Negative Results

The test is negative if a coloured line appears in the C area (control) only.

### Positive Results

The test is positive if one coloured line appears in the S (specimen) area and one coloured line appears in the C (control) area. Any coloured line in the specimen area should be considered positive. Coloured lines may appear lighter or darker than one another. Specimens with Legionella antigen levels near the threshold of the test may develop colour (faint lines) overtime after the 15 minute reading; however, a line that appears after 15 minutes is invalid and should be ignored.

### Invalid Results

The test is invalid if no coloured line appears in the C (control) area even if a coloured line appears in the S (specimen) area. If a coloured line does not appear in the C area, the test is invalid and should be repeated. Coloured lines that appear in the S (specimen) area after 15 minutes are invalid and should be ignored.

## 11. Quality Control

### Internal Controls

The appearance of a control line in the C region of the device is a positive procedural control. Correct procedural technique, specimen flow and device performance is confirmed when a coloured line appears in the C (control) area of the membrane. If the colored line fails to appear in the C (control) area, the test result is invalid.

A clear background is an internal negative procedural control. The background colour should be white to light pink and should not interfere with the reading of the test result. If any result is difficult to interpret, the test should be repeated with the same sample to eliminate the potential for error. Obtain a new sample and retest when the original sample repeatedly produces unreadable results. It is recommended that the internal quality controls be recorded for each sample run.

### External Controls

Quality control should be performed on each new test kit opened. Additional testing of the appropriate controls of *Legionella pneumophila* may be performed in accordance with federal, state and/or local regulations, accrediting requirements, and your laboratory's quality control procedures. Positive and negative external controls are supplied separately.

## 12. Limitations

- A negative antigen result does not rule out the possible infection from *Legionella pneumophila* serogroup 1. The antigen concentration may fall below the detectable limit of the test. Culture is recommended for patients with suspected pneumonia to determine if other causative agents other than *Legionella pneumophila* serogroup 1 are responsible.
- The test provides a presumptive diagnosis for only *Legionella pneumophila* serogroup 1 infection. In order to make an accurate diagnosis, culture results, serology testing and antigen detection methods in conjunction with clinical findings should be used.
- Sensitive immunoassays may demonstrate positive results with specimens containing heterophilic antibodies. If the qualitative interpretation is inconsistent with clinical findings, then further testing by an alternate method should be performed.

- Excretion of Legionella antigen may be seen as early as three (3) days after the onset of symptoms. However, antigen may continue to be excreted up to one (1) year afterwards. A positive result with the NADAL<sup>®</sup> Legionella Urine Antigen Test may occur during a current and past infection; therefore, in order to make an accurate diagnosis, culture results, serology testing and antigen detection methods in conjunction with clinical findings should be used.

### 13. Performance Characteristics

#### Clinical Sensitivity and Specificity (Retrospective Study)

Three clinical sites (USA and Netherlands) tested three hundred twenty four (324) frozen specimens using the NADAL<sup>®</sup> Legionella Urine Antigen Test. These samples were previously tested for Legionella by cell culture. Results are reported below.

		Cell Culture		
		+	-	
NADAL <sup>®</sup> Legionella Urine Antigen Test	+	95	11	106
	-	10	208	218
		105	219	324
Sensitivity: 95/105 x 100 =		90.5 % (95% CI: 83.4-94.8)		
Specificity: 208/219 x 100 =		95.0 % (95% CI: 91.2-97.5)		

#### Analytical Specificity

Forty-Nine (49) fresh patient urines from healthy individuals were collected and assayed at one (1) clinical site. One hundred percent (100%) of these were found to be negative by the NADAL<sup>®</sup> Legionella Urine Antigen Test.

#### Cross-Reactivity with other Respiratory Tract Infections

Ninety-nine (99) urine samples from patients diagnosed for other etiological respiratory tract infections (84 culture confirmed, 15 suspected) were tested using the NADAL<sup>®</sup> Legionella Urine Antigen Test. The results showed a lack of reactivity in 98/99 samples (99.0%).

#### Bacterial Cross-Reactivity and Interference

To confirm the analytical specificity of the NADAL<sup>®</sup> Legionella Urine Antigen Test, bacterial cultures likely to be found in the respiratory tract were tested. Bacterial cultures were tested at  $1 \times 10^5$  cfu/ml. All yielded negative results.

To confirm a lack of interference by other bacterial species in the NADAL<sup>®</sup> Legionella Urine Antigen Test, purified Legionella antigen was added to bacterial cultures likely to be found in the respiratory tract. The concentration of Legionella antigen was  $1 \times 10^5$  cfu/ml. All tests yielded positive results.

<i>Legionella pneumophila</i> SG2	<i>Legionella gormanii</i>	<i>Proteus mirabilis</i>
<i>Legionella pneumophila</i> SG3	<i>Legionella longbeachae</i>	<i>Citrobacter freundii</i>
<i>Legionella pneumophila</i> SG4	<i>Legionella feeleii</i>	<i>Enterobacter cloacae</i>
<i>Legionella pneumophila</i> SG5	<i>Staphylococcus aureus</i>	<i>Serratia marcescens</i>
<i>Legionella pneumophila</i> SG6	STAPHYLOCOCCUS EPIDERMIDIS	<i>Pseudomonas aeruginosa</i>
<i>Legionella bozemanii</i>	<i>Enterococcus faecalis</i>	<i>Candida albicans</i>
<i>Legionella micdadei</i>	<i>Escherichia coli</i>	<i>Streptococcus pneumoniae</i>
<i>Legionella dumoffii</i>	<i>Klebsiella pneumoniae</i>	<i>Haemophilus influenzae</i>

#### Common Urine Components and Drugs

The components and drugs listed below were tested to determine if they interfered or cross-reacted with NADAL<sup>®</sup> Legionella Urine Antigen Test. Extracted *Legionella pneumophila* antigen was added at a concentration of  $4.5 \times 10^4$  cfu/mL. In addition, each compound was tested in the absence of *Legionella pneumophila* antigen. None listed interfered or cross-reacted with the test results.

Glucose	2000mg/dL	Protein - BSA	500mg/dL
Urea	2000mg/dL	Bilirubin	20mg/dL
Ascorbic Acid	100mg/dL	Rifampicin	0.09mg/mL
Erythrocytes	$10^9$ /mL	Erythromycin	0.067mg/mL
Leukocytes	$10^9$ /mL	Ciprofloxacin	0.22mg/mL

#### Limit of Detection

The limit of detection of the NADAL<sup>®</sup> Legionella Urine Antigen Test was determined to be  $5 \times 10^4$ . *Legionella pneumophila* serogroup 1 ATCC 323152 was prepared using BCYE agar. A dilution of the working concentration was performed. The limit of detection of the NADAL<sup>®</sup> Legionella Urine Antigen Test was determined from these concentrations. The test was performed in 2 replicates. Results were recorded at 15 minutes.

#### Reproducibility Study

The reproducibility of the NADAL<sup>®</sup> Legionella Urine Antigen Test was evaluated at three clinical laboratory sites. The NADAL<sup>®</sup> Legionella Urine Antigen Test was tested against a panel of six (6) specimens of which included four levels of positives and two negatives. The low and high positives were from the purified Legionella antigen. Negatives were comprised of either urine or Legionella antigen below the detectable limit. Three (3) different laboratory personnel assayed each specimen at each laboratory facility over 3 days. The overall reproducibility for the NADAL<sup>®</sup> Legionella Urine Antigen Test was 100%.

#### 14. References

- Edelstein, P.H. and R. D. Meyer. "Legionella Pneumonias." In Respiratory infections: Diagnosis and Management, 3<sup>rd</sup> Edition, James E. Pennington, editor. New York: Raven Press, Ltd. 1994
- Marston, B.J., H.B. Lipman, R.F. Breiman. Surveillance for legionnaires' Disease: risk factors for morbidity and mortality. Arch.Intern. Med 1994;154:2417:2422
- Horwitz, M.A., B.J. Martson, C.V. Broome, and R.F. Freiman. Prospects for vaccine development. Presented at the 4<sup>th</sup> International Symposium on Legionella, 1992. In: Barbaree, J. M., R.F. Breiman, and A.P. DuFour, eds. Legionella: Current status and emerging perspectives. Washington, DC American Society for Microbiology, 1993.
- Kohler, R.B., W.C. Winn, Jr. and L.J. Wheat. Onset and duration of urinary antigen excretion in Legionnaires' disease. J. Clin. Microbiol. 1984; 20:605-607.
- Lindsay, D.S.J., W.H. Abraham, P. Christie, F. Johnston and G.F.S. Edwards. Laboratory diagnosis of legionnaires' disease due to *Legionella pneumophila* serogroup 1: comparison of phenotypic and genotypic methods. J. Med. Microbiol. 2004; 53:183-187.
- Cloud, J.L., K.C. Carroll, P. Pixton, M. Erali, and D.R. Hillyard. Detection of *Legionella* species in respiratory specimens using PCR with sequencing confirmation. J. Clin. Microbiol. 38; 2000; 1709-1712.
- Kohler, R.B., Winn, W.C., Wheat, L.J. Onset of the duration of urinary antigen excretion in Legionnaires' disease. J. Clin. Microbiol. 1984; 20: 601-605.
- Jonas, Daniel, A. Rosenbaum, S. Weyrick and S. Bhakdi. Enzyme-Linked Immunoassay for Detection of PCR-Amplified DNA and Legionellae in Bronchoalveolar fluid. J. Clin. Microbiol. 33; 1995; 1247-1252.

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## 14. Références

- Edelstein, P.H. and R. D. Meyer. "Legionella Pneumonias." In Respiratory Infections: Diagnosis and Management, 3rd Edition, James E. Pennington, editor. New York: Raven Press, Ltd. 1994
- Marston, B.J., H.B. Lipman, R.F. Breiman. Surveillance for legionnaires' Disease: risk factors for morbidity and mortality. Arch.Intern. Med 1994;154:2417:2422
- Horwitz, M.A., B.J. Martson, C.V. Broome, and R.F. Freiman. Prospects for vaccine development. Presented at the 4th International Symposium on Legionella, 1992. In: Barbaree, J. M., R.F. Breiman, and A.P. DuFour, eds. Legionella: Current status and emerging perspectives. Washington, DC American Society for Microbiology, 1993.
- Kohler, R.B., W.C. Winn, Jr. and L.J. Wheat. Onset and duration of urinary antigen excretion in Legionnaires' disease. J. Clin. Microbiol. 1984; 20:605-607.
- Lindsay, D.S.J., W.H. Abraham, P. Christie, F. Johnston and G.F.S. Edwards. Laboratory diagnosis of legionnaires' disease due to Legionella pneumophila serogroup 1: comparison of phenotypic and genotypic methods. J. Med. Microbiol. 2004; 53:183-187.
- Cloud, J.L., K.C. Carroll, P. Pixton, M. Erali, and D.R. Hillyard. Detection of Legionella species in respiratory specimens using PCR with sequencing confirmation. J. Clin. Micbiol. 38; 2000; 1709-1712.
- Kohler, RB, Winn, WC, Wheat, LJ. Onset of the duration of urinary antigen excretion in Legionnaires' disease. J. Clin. Microbiol. 1984; 20: 601-605.
- Jonas, Daniel, A. Rosenbaum, S. Weyrick and S. Bhakdi. Enzyme-Linked Immunoassay for Detection of PCR-Amplified DNA and Legionellae in Bronchoalveolar fluid. J. Clin. Micbiol. 33; 1995; 1247-1252.

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### 14. Referencias

1. Edelstein, P.H. and R. D. Meyer. "Legionella Pneumonias." In Respiratory infections: Diagnosis and Management, 3<sup>rd</sup> Edition, James E. Pennington, editor. New York: Raven Press, Ltd. 1994
2. Marston, B.J., H.B. Lipman, R.F. Breiman. Surveillance for legionnaires' Disease: risk factors for morbidity and mortality. Arch.Intern. Med 1994;154:2417:2422
3. Horwitz, M.A., B.J. Martson, C.V. Broome, and R.F. Freiman. Prospects for vaccine development. Presented at the 4<sup>th</sup> International Symposium on Legionella, 1992. In: Barbaree, J. M., R.F. Breiman, and A.P. DuFour, eds. Legionella: Current status and emerging perspectives. Washington, DC American Society for Microbiology, 1993.
4. Kohler, R.B., W.C. Winn, Jr. and L.J. Wheat. Onset and duration of urinary antigen excretion in Legionnaires' disease. J. Clin. Microbiol. 1984; 20:605-607.
5. Lindsay, D.S.J., W.H. Abraham, P. Christie, F. Johnston and G.F.S. Edwards. Laboratory diagnosis of legionnaires' disease due to Legionella pneumophila serogroup 1: comparison of phenotypic and genotypic methods. J. Med. Microbiol. 2004; 53:183-187.
6. Cloud, J.L., K.C. Carroll, P. Pixton, M. Erali, and D.R. Hillyard. Detection of Legionella species in respiratory specimens using PCR with sequencing confirmation. J. Clin. Micbiol. 38; 2000; 1709-1712.
7. Kohler, RB, Winn, WC, Wheat, LJ. Onset of the duration of urinary antigen excretion in Legionnaires' disease. J. Clin. Microbiol. 1984; 20: 601-605.
8. Jonas, Daniel, A. Rosenbaum, S. Weyrick and S. Bhakdi. Enzyme-Linked Immunoassay for Detection of PCR-Amplified DNA and Legionellae in Bronchoalveolar fluid. J. Clin. Micbiol. 33; 1995; 1247-1252.

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

## 14. Bibliografia

- Edelstein, P.H. and R. D. Meyer. "Legionella Pneumonias." In Respiratory infections: Diagnosis and Management, 3<sup>rd</sup> Edition, James E. Pennington, editor. New York: Raven Press, Ltd. 1994
- Marston, B.J., H.B. Lipman, R.F. Breiman. Surveillance for legionnaires' Disease: risk factors for morbidity and mortality. Arch.Intern. Med 1994;154:2417:2422
- Horwitz, M.A., B.J. Martson, C.V. Broome, and R.F. Freiman. Prospects for vaccine development. Presented at the 4<sup>th</sup> International Symposium on Legionella, 1992. In: Barbaree, J. M., R.F. Breiman, and A.P. DuFour, eds. Legionella: Current status and emerging perspectives. Washington, DC American Society for Microbiology, 1993.
- Kohler, R.B., W.C. Winn, Jr. and L.J. Wheat. Onset and duration of urinary antigen excretion in Legionnaires' disease. J. Clin. Microbiol. 1984; 20:605-607.
- Lindsay, D.S.J., W.H. Abraham, P. Christie, F. Johnston and G.F.S. Edwards. Laboratory diagnosis of legionnaires' disease due to Legionella pneumophila serogroup 1: comparison of phenotypic and genotypic methods. J. Med. Microbiol. 2004; 53:183-187.
- Cloud, J.L., K.C. Carroll, P. Pixton, M. Erali, and D.R. Hillyard. Detection of Legionella species in respiratory specimens using PCR with sequencing confirmation. J. Clin. Micbiol. 38; 2000; 1709-1712.
- Kohler, RB, Winn, WC, Wheat, LJ. Onset of the duration of urinary antigen excretion in Legionnaires' disease. J. Clin. Microbiol. 1984; 20: 601-605.
- Jonas, Daniel, A. Rosenbaum, S. Weyrick and S. Bhakdi. Enzyme-Linked Immunoassay for Detection of PCR-Amplified DNA and Legionellae in Bronchoalveolar fluid. J. Clin. Micbiol. 33; 1995; 1247-1252.

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nal von minden GmbH