NADAL[®] Hb/Hp Complex plus Test (test cassette)

REF 272011N-25

DE	Gebrauchsanweisung	2	PL	Sposób użycia	17
EN	Instructions for use	5	CZ	Návod k použití	20
FR	Instructions d'utilisation	8		Symbols	23
ES	Instrucciones de uso	11		Our Teams	24
IT	Istruzioni per l'uso	14			

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Carl-Zeiss-Strasse 12
47445 Moers
Germany

Tel: +49 (2841) 99820-0
Fax: +49 (2841) 99820-1

Friedenstrasse 32
93053 Regensburg
Germany

Tel: +49 941 29010-0
Fax: +49 941 29010-50

www.nal-vonminden.com
info@nal-vonminden.com

Directors:
Sandra von Minden
Roland Meißner
Thomas Zander

Commercial reg. Kleve
HRB 5679
Steuer-Nr. 244/133/00130
UST-ID-Nr. DE 189 016 086

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14. Referenzen

1. Dam, J.V., et. al.; Fecal Blood Screening for Colorectal Cancer; *Archive of Internal Medicine*; (1995) 155:2389-2402
2. Frommer, D.J. et. al.; Improved Screening for Colorectal Cancer by Immunological Detection of Occult Blood; *British Medical Journal*; (1988) 296:1092-1094
3. Lieberman, D.; Screening/Early Detection Model for Colorectal Cancer, Why Screen? *Cancer Supplement*; (1994) 74(7): 2023-2027
4. Miller, A.B.; An Epidemiological Perspective on Cancer Screening; *Clinical Biochemistry* (1995) 28(1) : 41-48
5. Ransohoff, D.F. and Lang, C.A.; Improving the Fecal Occult- Blood Test: The New England Journal of Medicine; (1996) 334 (3) : 189-190
6. Screening for Colorectal Cancer-United States, 1992-1993, and New Guidelines; *Morbidity and Mortality Weekly Report*; (1995) 45 (5) : 107- 110
7. St. John, D.J.B., et. al.; Evaluation of New Occult Blood Test for Detection of Colorectal Neoplasia; *Gastroenterology*; (1993) 104:1661-1668
8. Yamamoto M.; Nakama H.; Cost-effectiveness analysis of immunochemical occult blood screening for colorectal cancer among three fecal sampling methods: *Hepatogastroenterology*; 2000 Mar-Apr; 47(32):396-9.

Rev. 3, 2016-02-26 OM

1. Intended Use

The NADAL® Hb/Hp Complex plus is a visual, immuno-chromatographic dual rapid test for the qualitative detection of human hemoglobin and hemoglobin-haptoglobin complexes in human stool samples.

2. Clinical Significance

The NADAL® Hb/Hp Complex plus is intended for professional *in-vitro* diagnostic use. The test is designed to provide evidence of bowel cancer and/or adenomas. Bowel cancer is one of the most frequently diagnosed types of cancer and one of the most frequent causes of death from cancer. Bowel cancer can be detected early by means of tests for hidden blood in stools and the mortality rate can therefore be reduced. The hemoglobin/haptoglobin (Hb/Hp) complex plays an important role in the retrieval of hemoglobin from lysed erythrocytes. It is relatively stable against acid and proteolytic degradation and can also be detected if it is set free in the anterior intestinal region. Blood from intestinal polyps and higher carcinomas is also detectable. The NADAL® Hb/Hp Complex plus shows a high sensitivity for the detection of colorectal adenomas and carcinomas. Detection of hemoglobin and hemoglobin-haptoglobin complex is therefore an innovative and extremely sensitive method for the early diagnosis of colon carcinomas and larger adenomas.

3. Test Principle

The immunological NADAL® Hb/Hp Complex plus is based on an immunochromatographic method, in which both Hb and Hb/Hp complexes are detected through specific antibody reactions.

For Hb:

Hemoglobin contained in stool sample reacts with specific gold-conjugated monoclonal antibodies. This complex diffuses along the membrane and binds to the anti-hemoglobin antibody which is immobilised on the test line (T).

For Hb/Hp:

An antibody attaches to hemoglobin and another antibody attaches to haptoglobin. These Hb/Hp complexes are measured and there is no cross-reactivity with hemoglobin. In case of a positive result, hemoglobin and haptoglobin from the stool sample bound to gold-conjugated antibodies bind to immobilised antibodies on the test line (T) and a red/pink colouration becomes visible. In case of a negative result, there are no Hb/Hp complexes that can bind to the antibodies of the test line (T) and no colouration of the test line (T) is formed. A red/pink control line (C) implies that the test has been performed correctly. The test is valid.

4. Reagents and Materials Supplied

- 25 individually packed test cassettes
- 25 sample-taking tubes containing extraction solution
- 25 sets of instructions on taking stool samples
- 25 stool containers
- 1 package insert

5. Additional Materials Required

- 1 stopwatch

6. Storage & Stability

The test can be stored at +2° to +30°C in a sealed bag until the expiration date.

7. Warnings and Precautions

- Only for professional *in-vitro* diagnostic use.
- Read the entire package insert carefully prior to testing.
- Do not use the test after the expiration date.
- Do not use the test if the foil pouch is damaged. Do not reuse the test.
- Use a new stool container for each specimen.
- It is recommended to use the enclosed stool container in order to avoid mixing with blood from urine samples or interfering components from the toilet water.
- Do not eat, drink or smoke in any area where specimens and kits are handled. Handle all specimens as if they contain infectious agents. Observe established precautions against microbiological hazards throughout testing and follow standard procedures for the proper disposal of specimens. Wear protective clothing such as laboratory coats, disposable gloves and eye protection when specimens are assayed.
- Only open the test packaging when the test is actually going to be performed.
- The extraction solution contains a small amount of sodium azide which may react with lead or copper plumbing to form potentially explosive metal azides. When disposing of extraction solution or extracted samples, always flush with copious quantities of water to prevent azide buildup.
- Do not interchange or mix reagents from different lots.
- Humidity and temperature can adversely affect results.
- Used testing materials should be discarded according to local regulations.
- This kit contains products of animal origin. Certified knowledge of the origin and/or sanitary state of the animals does not completely guarantee the absence of transmissible pathogenic agents. It is therefore recommended that these products be treated as potentially infectious and handled by observing usual safety precautions (e.g. do not ingest or inhale).
- Patients should follow the instructions for obtaining a sample exactly and should not take any samples during menstruation or 3 days before or after, in case of bleeding hemorrhoids, bleeding caused by constipation or rectal medicine administration.
- Alcohol, aspirin and other medications taken in excess may cause gastrointestinal irritation resulting in occult bleeding. Such substances should be discontinued at least 48 hours prior to testing.
- No dietary restrictions are necessary before testing.

8. Specimen Collection and Preparation

Note: Please give patients the following information. During sampling, it is important that the sample-taking tube (spiral rod) is inserted **three times** in a row into **three different locations of the stool sample**.

The sample-taking tube should NOT be placed back into the tube between the three collections. No fluid should leak from the tube. **The nature of the test sample and the compliance with instructions will affect the test result.**

Bring the sample-taking tube to room temperature before use.
Collect the stool using the stool container enclosed.

Hold the sample-taking tube vertically and open the sealing cap. Remove the spiral rod.

Place the spiral rod at three different locations in the stool sample.

Place the spiral rod with the stool sample into the tube and seal well.

Label the well-sealed sample-taking tube with the name of the patient.


Perform testing immediately after specimen collection. Do not leave specimens at room temperature for prolonged periods of time. Specimens may be stored at 2-8°C for up to 72 hours.

If specimens are to be shipped, pack them in compliance with all applicable regulations for transportation of etiological agents.

If the sample is taken at the patient's home, the sample should be protected from light and stored in a cool place.

9. Test Procedure

- Bring the test cassette and the patient's sample to room temperature (+20°C to +30°C) before testing.
- Take the required number of test cassettes from the foil packaging only immediately before performing the test. Mark the cassette with the name of the patient or with any form of identification.
- Shake the stool sample tube carefully to ensure that the stool sample mixes properly with the extraction solution.
- Unscrew the protective cap. Using a paper towel break the seal of the sample-taking tube with a rotary motion. Hold the sample-taking tube upright and add 3 drops of the solution mixture to each round sample window (S) of the cassette. Avoid trapping air bubbles in the specimen wells (S) and do not add any solution to the result area.
- Read the results after 5 minutes. Strong positive results may be evaluated even sooner. Evaluate the result in a maximum of 15 minutes after dropping the solution mixture into the round sample windows.

10. Result Interpretation

Test result is evaluated as "positive" if in addition to the control lines at least one of the two test lines develops.

Positive for Hb:

One red/pink line develops in each control line region (C) and one red/pink line develops in the test line region (T) for Hb.

Positive for Hb/Hp:

One red/pink line develops in each control line region (C) and one red/pink line develops in the test line region (T) for Hb/Hp.

Positive for Hb and Hb/Hp:

In addition to the control lines (C), a red/pink line develops in each test line region for Hb and Hb/Hp.

Note: If you are testing strongly positive samples, the intensity of the control line may be reduced. It is not recommended to compare the intensity of the lines.


Negative:

One red/pink line develops in each control line region (C). No line develops in each test line region (T).

Invalid:

The absence of the control lines (C) indicates a failure in the test procedure. The sample should be retested.

Note: Insufficient specimen volume, incorrect procedural techniques are the most likely reasons for control line failure. Review the procedure and repeat the test with a new test cassette. If the problem persists, discontinue using the test kit immediately and contact your distributor.

11. Quality Control

The test contains a procedural control. A red/pink line that develops in each control line region (C) indicates that the test has been performed properly.

The background can become slightly yellowish during testing, depending on the colour of the stool sample. This is acceptable as long as the evaluation of the test results is not adversely affected.

12. Limitations

The test is used for the qualitative detection of human hemoglobin and hemoglobin-haptoglobin complexes in stool samples.

A positive test result will be obtained if human hemoglobin or hemoglobin-haptoglobin complexes are present in the sample. In addition to intestinal bleeding, blood in stools may also have other causes, such as hemorrhoids or a mixture of the stool sample with urine containing blood.

Not all intestinal bleedings are the result of benign or malignant polyps. The information obtained from this test should be examined in conjunction with other clinical test methods.

Carcinomas and polyps can lead to bleeding in different extents and intermittently. It is advisable to check various stool samples by means of immunological verification procedures.

Urine and excessive dilution of samples with water from the toilet bowl can lead to false test results. It is therefore recommended to use the provided stool container.

Stool samples should not be taken during menstruation or 3 days before or after, in case of bleeding hemorrhoids bleeding caused by constipation or rectal medicine administration. This could lead to false positive results.

Blood-thinning medications, such as for example, acetylsalicylic acid or coumarin as well as iron preparations can lead to bleeding unrelated to tumors.

Unlike test methods that detect only free hemoglobin, the NADAL® Hb/Hp Complex plus is very sensitive in detecting bleeding in the anterior intestinal passage.

13. Performance Characteristics

Analytical sensitivity

The detection limit of the NADAL® Hb/Hp Complex plus is 25 ng of free or complex hemoglobin per milliliter of buffer solution.

Prozone effect:

The NADAL® Hb/Hp Complex plus works reliably at extremely high hemoglobin concentrations (>500,000 ng/mL).

Specificity

The NADAL® Hb/Hp Complex plus is specific for human hemoglobin or hemoglobin-haptoglobin complex and has no cross-reactivity with hemoglobin from cattle, pigs, rabbits, horses or sheep at a concentration of up to 0.5 mg/mL in the extraction buffer. The NADAL® Hb/Hp Complex plus also shows no cross-reactivity with Bilirubin, L-ascorbate (vitamin C), and horseradish peroxidase.

14. References

1. Dam, J.V., et. al.; Fecal Blood Screening for Colorectal Cancer; Archive of Internal Medicine; (1995) 155:2389-2402
2. Frommer, D.J. et. al.; Improved Screening for Colorectal Cancer by Immunological Detection of Occult Blood; British Medical Journal; (1988) 296:1092-1094
3. Lieberman, D.; Screening/Early Detection Model for Colorectal Cancer, Why Screen? Cancer Supplement; (1994) 74(7) : 2023-2027
4. Miller, A.B.; An Epidemiological Perspective on Cancer Screening; Clinical Biochemistry (1995) 28(1) : 41-48
5. Ransohoff, D.F. and Lang, C.A.; Improving the Fecal Occult- Blood Test: The New England Journal of Medicine; (1996) 334 (3) : 189-190
6. Screening for Colorectal Cancer-United States, 1992-1993, and New Guidelines; Mobility and Mortality Weekly Report; (1995) 45 (5) : 107- 110
7. St. John, D.J.B., et. al.; Evaluation of New Occult Blood Test for Detection of Colorectal Neoplasia; Gastroenterology; (1993) 104:1661-1668
8. Yamamoto M.; Nakama H.; Cost-effectiveness analysis of immunochemical occult blood screening for colorectal cancer among three fecal sampling methods: Hepatogastroenterology; 2000 Mar-Apr; 47(32):396-9.

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14. Bibliographie

1. Dam, J.V., et. al.; Fecal Blood Screening for Colorectal Cancer; *Archive of Internal Medicine*; (1995) 155:2389-2402
2. Frommer, D.J. et. al.; Improved Screening for Colorectal Cancer by Immunological Detection of Occult Blood; *British Medical Journal*, (1988) 296:1092-1094
3. Lieberman, D.; Screening/Early Detection Model for Colorectal Cancer, *Why Screen? Cancer Supplement*; (1994) 74(7) : 2023-2027
4. Miller, A.B.; An Epidemiological Perspective on Cancer Screening; *Clinical Biochemistry* (1995) 28(1) : 41-48
5. Ransohoff, D.F. and Lang, C.A.; Improving the Fecal Occult- Blood Test: The *New England Journal of Medicine*; (1996) 334 (3) : 189-190
6. Screening for Colorectal Cancer-United States, 1992-1993, and New Guidelines; *Morbidity and Mortality Weekly Report*; (1995) 45 (5) : 107- 110
7. St. John, D.J.B., et. al.; Evaluation of New Occult Blood Test for Detection of Colorectal Neoplasia; *Gastroenterology*; (1993) 104:1661-1668
8. Yamamoto M.; Nakama H.; Cost-effectiveness analysis of immunochemical occult blood screening for colorectal cancer among three fecal sampling methods: *Hepato-gastroenterology*; 2000 Mar-Apr; 47(32):396-9.

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14. Referencias

1. Dam, J.V., et. al.; Fecal Blood Screening for Colorectal Cancer; Archive of Internal Medicine; (1995) 155:2389-2402
2. Frommer, D.J. et. al.; Improved Screening for Colorectal Cancer by Immunological Detection of Occult Blood; British Medical Journal; (1988) 296:1092-1094
3. Lieberman, D.; Screening/Early Detection Model for Colorectal Cancer, Why Screen? Cancer Supplement; (1994) 74(7): 2023-2027
4. Miller, A.B.; An Epidemiological Perspective on Cancer Screening; Clinical Biochemistry (1995) 28(1) : 41-48
5. Ransohoff, D.F. and Lang, C.A.; Improving the Fecal Occult- Blood Test: The New England Journal of Medicine; (1996) 334 (3) : 189-190

6. Screening for Colorectal Cancer-United States, 1992-1993, and New Guidelines; Mobility and Mortality Weekly Report; (1995) 45 (5) : 107- 110
7. St. John, D.J.B., et. al.; Evaluation of New Occult Blood Test for Detection of Colorectal Neoplasia; Gastroenterology; (1993) 104:1661-1668
8. Yamamoto M.; Nakama H.; Cost-effectiveness analysis of immunochemical occult blood screening for colorectal cancer among three fecal sampling methods: Hepatogastroenterology; 2000 Mar-Apr; 47(32):396-9.

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neanche con Bilirubina, L-Ascorbato (vitamina C) o Perossidasi del Rafano.

14. Bibliografia

1. Dam, J.V., et. al.; Fecal Blood Screening for Colorectal Cancer; *Archive of Internal Medicine*; (1995) 155:2389-2402
2. Frommer, D.J. et. al.; Improved Screening for Colorectal Cancer by Immunological Detection of Occult Blood; *British Medical Journal*; (1988) 296:1092-1094
3. Lieberman, D.; Screening/Early Detection Model for Colorectal Cancer, Why Screen? *Cancer Supplement*; (1994) 74(7) : 2023-2027
4. Miller, A.B.; An Epidemiological Perspective on Cancer Screening; *Clinical Biochemistry* (1995) 28(1) : 41-48
5. Ransohoff, D.F. and Lang, C.A.; Improving the Fecal Occult- Blood Test: The New England Journal of Medicine; (1996) 334 (3) : 189-190
6. Screening for Colorectal Cancer-United States, 1992-1993, and New Guidelines; *Mobility and Mortality Weekly Report*; (1995) 45 (5) : 107- 110
7. St. John, D.J.B., et. al.; Evaluation of New Occult Blood Test for Detection of Colorectal Neoplasia; *Gastroenterology*; (1993) 104:1661-1668
8. Yamamoto M.; Nakama H.; Cost-effectiveness analysis of immunochemical occult blood screening for colorectal cancer among three fecal sampling methods: *Hepatogastroenterology*; 2000 Mar-Apr; 47(32):396-9.

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14. Bibliografia

1. Dam, J.V., et. al.; Fecal Blood Screening for Colorectal Cancer; Archive of Internal Medicine; (1995) 155:2389-2402
2. Frommer, D.J. et. al.; Improved Screening for Colorectal Cancer by Immunological Detection of Occult Blood; British Medical Journal; (1988) 296:1092-1094

3. Lieberman, D.; Screening/Early Detection Model for Colorectal Cancer, Why Screen? Cancer Supplement; (1994) 74(7) : 2023-2027
4. Miller, A.B.; An Epidemiological Perspective on Cancer Screening; Clinical Biochemistry (1995) 28(1) : 41-48
5. Ransohoff, D.F. and Lang, C.A.; Improving the Fecal Occult- Blood Test: The New England Journal of Medicine; (1996) 334 (3) : 189-190
6. Screening for Colorectal Cancer-United States, 1992-1993, and New Guidelines; Mobility and Mortality Weekly Report; (1995) 45 (5) : 107- 110
7. St. John, D.J.B., et. al.; Evaluation of New Occult Blood Test for Detection of Colorectal Neoplasia; Gastroenterology; (1993) 104:1661-1668
8. Yamamoto M.; Nakama H.; Cost-effectiveness analysis of immunochemical occult blood screening for colorectal cancer among three fecal sampling methods: Hepatogastroenterology; 2000 Mar-Apr; 47(32):396-9.

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14. Literatura

1. Dam, J.V., et. al.; Fecal Blood Screening for Colorectal Cancer; Archive of Internal Medicine; (1995) 155:2389-2402.
2. Frommer, D.J. et. al.; Improved Screening for Colorectal Cancer by Immunological Detection of Occult Blood; British Medical Journal; (1988) 296:1092-1094.
3. Lieberman, D.; Screening/Early Detection Model for Colorectal Cancer, Why Screen? Cancer Supplement; (1994) 74(7) : 2023-2027.
4. Miller, A.B.; An Epidemiological Perspective on Cancer Screening; Clinical Biochemistry (1995) 28(1) : 41-48.
5. Ransohoff, D.F. and Lang, C.A.; Improving the Fecal Occult- Blood Test: The New England Journal of Medicine; (1996) 334 (3) : 189-190.
6. Screening for Colorectal Cancer-United States, 1992-1993, and New Guidelines; Mobility and Mortality Weekly Report; (1995) 45 (5) : 107- 110.
7. St. John, D.J.B., et. al.; Evaluation of New Occult Blood Test for Detection of Colorectal Neoplasia; Gastroenterology; (1993) 104:1661-1668.
8. Yamamoto M.; Nakama H.; Cost-effectiveness analysis of immunochemical occult blood screening for colorectal cancer among three fecal sampling methods: Hepatogastroenterology; 2000 Mar-Apr; 47(32):396-9.

Rev. 3, 2016-02-26 TF

Symbol	Deutsch	English	Français	Español	Italiano	Polski
	CE Konformitätszeichen	CE marking of conformity	Conformité aux normes européennes	Conformidad europea	Conformità europea	Znak zgodności CE
	Gebrauchsanweisung beachten	Consult instructions for use	Consulter la notice d'utilisation	Consúltense las instrucciones de uso	Consultare le istruzioni per l'uso	Przestrzegać instrukcji obsługi
	In-vitro-Diagnostika	In-vitro diagnostic medical device	Dispositif médical de diagnostic in vitro	Producto sanitario para diagnóstico in vitro	Dispositivo medico-diagnostico in vitro	Tylko do diagnostyki in vitro
	Temperaturbegrenzung	Temperature limitation	Limites de température	Límite de temperatura	Limiti di temperatura	Temperatura przechowywania
	Chargenbezeichnung	Batch code	Code du lot	Código de lote	Codice lotto	Numer serii
	Nicht zur Wiederverwendung	Do not reuse	Ne pas réutiliser	No reutilizar	Non riutilizzare	Tylko do jednorazowego użytku
	Verwendbar bis	Use by	Utiliser jusqu'au	Fecha de caducidad	Utilizzare entro	Data ważności
	Bestellnummer	Catalogue Number	Référence du catalogue	Número de catálogo	Riferimento di Catalogo	Numer katalogowy
	Hersteller	Manufacturer	Fabricant	Fabricante	Fabbricante	Producent
	Ausreichend für <n> Ansätze	Sufficient for <n> tests	Suffisant pour "n" tests	Suficiente para <n> utilizaciones	Sufficiente per "n" saggi	Wystarczający na <n> Powtórzeń

Symbol	Português	Český	Suomi	Svenskt	Nederlands	Dansk	Norsk
	Conformidade com as normas europeias	CE certifikát	CE-merkitty	CE-märkning	CE-markering	CE-mærkning	CE standardisert
	Consultar as instruções de utilização	Viz návod k použití	Katso käyttöohjetta	Läs bruksanvisningen	Raadpleeg de gebruiksaanwijzing	Se brugsanvisningen	Les bruksanvisning nøye
	Dispositivo médico para diagnóstico in vitro	Diagnostický zdravotnický prostředek in vitro	In vitro -diagnostiikkaan tarkoitettu lääkinnällinen laite	Medicinteknisk produkt avsedd för in vitro-diagnostik	Medisch hulpmiddel voor in-vitrodiagnostiek	Medicinsk udstyr til in vitro-diagnostik	In-vitro diagnostic medisinsk enhet
	Limites de temperatura	Teplotní omezení	Lämpötilarajat	Temperatur-begränsning	Temperatuurlimiet	Temperatur-begrænsning	Temperatur begrensning
	Código do lote	Kód šarže	Eräkoodi	Satsnummer	Code van de partij	Batchkode	Merking
	Não reutilizar	Pro jednorázové použití	Kertakäyttöinen	Får inte återanvändas	Niet opnieuw gebruiken	Må ikke genbruges	Må ikke brukes om igjen
	Prazo de validade	Spotřebujte do	Käytettävä viimeistään	Används före	Houdbaar tot	Udløbsdato	Tidtaking
	Número de catálogo	Katalogov číslo	Luettelonumero	Listnummer	Catalogus nummer	Best il lingsnummer	Katalog nummer
	Fabricante	Výrobce	Valmistaja	Tilverkare	Fabrikant	Fabrikant	Produsent
	Suficiente para <n> test	Dostačuje pro <n> testů	Lukumäärä <n> test	Räcker till <n> test	Voldoende voor <n> test	Tilstrækkeligt til <n> test	Tilstrækkelig for <n> tester

Our Teams

Germany:

Regensburg

Tel: +49 941 290 10-0

Fax: +49 941 290 10-50

Moers

Tel: +49 2841 99820-0

Fax: +49 2841 99820-1

Austria:

Tel: +49 941 290 10-29

Free Tel: 0800 291 565

Fax: +49 290 10-50

Free Fax: 0800 298 197

UK & Ireland:

Tel: +49 941 290 10-18

Free Tel –UK: 0808 234 1237

Free Tel – IRE: 1800 555 080

Fax: +49 290 10-50

France:

France Tel: 0800 915 240

France Fax: 0800 909 493

Switzerland

Swiss Tel: 0800 564 720

Swiss Fax: 0800 837 476

Belgium

Belgium Tel: 0800 718 82

Belgium Fax: 0800 747 07

Luxembourg

Lux. Tel: 800 211 16

Lux. Fax: 800 261 79

Spain:

Tel: +49 941 290 10-759

Free Tel: 900 938 315

Fax: +49 941 290 10-50

Free Fax: 900 984 992

Italy:

Tel: +49 941 290 10-34

Fax: +49 941 290 10-50

Poland:

Tel: +49 941 290 10-44

Free Tel: 00 800 491 15 95

Fax: +49 941 290 10-50

Free Fax: 00 800 491 15 94

Portugal:

Tel: +49 941 290 10-735

Tel. Verde: 800 849 230

Fax: +49 941 290 10-50

Fax Verde: 800 849 229

Netherlands:

Tel: +31 30 75 600

Free Tel: 0800 0222 890

Fax: +31 70 30 30 775

Free Fax: 0800 024 9519

Nordic countries (Finland, Norway, Sweden, Denmark):

Tel: +31 703075 607

Free Tel: +45 80 88 87 53

Tax: +31 703030 775

Laboratory Diagnostics Team:

Tel: +49 941 290 10-40

Fax: +49 941 290 10-50

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Carl-Zeiss-Strasse 12 • 47445 Moers • Germany

www.nal-vonminden.com • info@nal-vonminden.com

Fon: +49 2841 99820-0 • Fax: +49 2841 99820-1

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nal von minden GmbH