NADAL[®] Syphilis Test (test cassette)

REF 202001

DE Gebrauchsanweisung	2	DK Brugervejledning	23
EN Instruction for use	5	Symbols	27
FR Instructions d'utilisation	8	Our Teams	28
ES Instrucciones de uso	11		
IT Istruzioni per l'uso	14		
PL Sposób użycia	17		
CZ Návod k použití	20		

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14. Referenzen

1. Claire M. Fraser. Complete genome sequence of *Treponema Pallidum*, the syphilis spirochete, *Science* 1998; 281 July: 375-381.
2. Center for Disease Control. Recommendations for diagnosing and treating syphilis in HIV-infected patients, *MMWR Morb. Mortal Wkly Rep.* 1988; 37: 601
3. Aral R. Marx. Crack, sex and STD, *Sexually Transmitted Diseases*, 1991; 18:92-101
4. J.N. Wasserheit. Epidemiological Synergy: Interrelationships between human immunodeficiency virus infection and other sexually transmitted diseases, *Sexually Transmitted Diseases* 1992; 19:61-77
5. Johnson Phillip C. Testing for Syphilis, *Dermatologic Clinic* 1994; 12 Jan: 9-17.

Rev. 1, 2015-07-17 OM

1. Intended Use

The NADAL® Syphilis test is a rapid chromatographic immunoassay for the qualitative detection of antibodies (IgG and IgM) to *Treponema pallidum* (*T. pallidum*) in whole blood, serum or plasma to aid in the diagnosis of Syphilis.

2. Introduction and Clinical Significance

Treponema pallidum (*T. pallidum*) is the causative agent of the venereal disease Syphilis. *T. pallidum* is a spirochete bacterium with an outer cell wall and a cytoplasmic membrane.¹

Relatively little is known about the organism in comparison with other bacterial pathogens. According to the Center for Disease Control (CDC), the number of cases of Syphilis infection has increased significantly since 1985.² Some key factors that have contributed to this rise include the crack cocaine epidemic and the high incidence of prostitution among drug users.³ One study reported a substantial epidemiological correlation between the acquisition and transmission of the HIV virus and Syphilis.⁴

Multiple clinical stages and long periods of latent, asymptomatic infection are characteristic of Syphilis. Primary Syphilis is defined by the presence of a chancre. The antibodies response to the *T. pallidum* bacterium can be detected within 4 to 7 days after the chancre appears. The infection remains detectable until the patient receives adequate treatment.⁵

The NADAL® Syphilis test utilises a double antigen combination of a Syphilis antigen coated particle and Syphilis antigen immobilized on membrane to detect *T. pallidum* antibodies (IgG and IgM) qualitatively and selectively in whole blood, serum or plasma.

3. Test Principle

The NADAL® Syphilis test is a qualitative membrane-based immunoassay for the detection of *T. pallidum* antibodies (IgG and IgM) in whole blood, serum or plasma. In the test procedure, recombinant Syphilis antigen is immobilised in the test line region of the test cassette.

After the specimen is added to the specimen well of the test cassette, it reacts with Syphilis antigen coated particles. The mixture migrates chromatographically along the test cassette and interacts with the immobilised Syphilis antigen. The double antigen test format can detect both IgG and IgM in specimens. If the specimen contains *T. pallidum* antibodies, a coloured line will develop in the test line region, indicating a positive result. If the specimen does not contain *T. pallidum* antibodies, no line will develop in the test line region, indicating a negative result. A coloured line should always develop in the control line region. It serves as an internal test control and as verification that sufficient volume of specimen has been added and proper membrane wicking has occurred.

4. Reagents and Materials Supplied

- 20 NADAL® Syphilis test cassettes
- 20 disposable pipettes
- 1 Buffer (only for whole blood)
- 1 package insert

5. Additional Materials Required

- Specimen collection containers

- Lancets (only for fingerstick whole blood)
- Heparinised capillary tubes and dispensing bulb (only for fingerstick whole blood)
- Centrifuge
- Timer

6. Storage & Stability

Store the test as packaged in the sealed foil pouch either at room temperature or refrigerated (2-30°C). The test cassette is stable up to the expiry date printed on the sealed foil pouch. The test must remain in the sealed foil pouch until use. Do not freeze the test. Do not use the test beyond the expiration date indicated on the package.

7. Warnings and Precautions

- For professional *in-vitro* diagnostic use only.
- Do not use the test beyond the expiration date indicated on the package.
- Do not eat, drink or smoke in the area where specimens or test kits are being handled.
- Do not use the test if the foil pouch is damaged.
- Handle all specimens as if they contain infectious agents. Observe established precautions for microbiological risks throughout all procedures and standard guidelines for the appropriate disposal of specimens.
- Wear protective clothing such as laboratory coats, disposable gloves and eye protection when specimens are being assayed.
- Used testing materials should be discarded according to local regulations.
- Humidity and temperature can adversely affect test results.

8. Specimen Collection, Preparation and Storage

The NADAL® Syphilis test can be performed using whole blood (by venipuncture or fingerstick), serum or plasma.

Collection of fingerstick whole blood:

Wash the patient's hand with soap and warm water or clean with an alcohol swab. Allow to dry.

Massage the hand without touching the puncture site by rubbing down the hand towards the fingertip of the middle or ring finger.

Puncture the skin with a sterile lancet. Wipe away the first sign of blood.

Gently rub the hand from wrist to palm to finger to form a rounded drop of blood over the puncture site.

Add the fingerstick whole blood specimen to the test cassette using a **capillary tube**:

- Place the end of the capillary tube to the blood until it is filled to approximately 50 µL. Avoid air bubbles.
- Attach the bulb onto the top end of the capillary tube, then squeeze the bulb to dispense the whole blood to the specimen well (S) of the test cassette.

Add the fingerstick whole blood specimen to the test using **hanging drops**:

- Position the patient's finger so that the drop of blood is right above the specimen well (S) of the test cassette.
- Allow 2 hanging drops of fingerstick whole blood to fall into the center of the specimen well (S) on the test cassette, or move the patient's finger so that the hanging drops touch

the center of the specimen well (S). Avoid direct finger contact with the specimen well (S).

Whole blood collection by venipuncture

Collect whole blood by venipuncture following the standard procedure.

Collection of serum and plasma

Separate serum or plasma from blood by centrifugation as soon as possible to avoid hemolysis. Use only clear non-hemolysed specimens.

Storage

Perform the test immediately after specimen collection. Do not leave the specimens at room temperature for prolonged periods. Serum and plasma specimens may be stored at 2-8°C for up to 3 days. For long term storage, specimens should be kept below -20°C. Whole blood collected by venipuncture should be stored at 2-8°C if the test is to be run within 2 days of collection. Do not freeze whole blood specimens. Whole blood collected by fingerstick should be tested immediately. Bring specimens to room temperature prior to testing. Frozen specimens must be completely thawed and mixed well prior to testing. Specimens should not be frozen and thawed repeatedly.

If specimens are to be shipped, they should be packed in compliance with local regulations covering the transportation of etiologic agents.

9. Test Procedure

Allow the test cassettes, specimens and buffer to reach room temperature (15-30°C) prior to testing.


Bring the pouch to room temperature before opening it. Remove the test cassette from the sealed pouch and use it as soon as possible.

Place the test cassette on a clean and level surface.

1a) Fingerstick whole blood specimens:


I. **To use a capillary tube:** Fill the capillary tube and transfer approximately 50 µL of fingerstick whole blood to the specimen well (S) of the test cassette, then add 1 drop of buffer (approximately 40 µL) and start the timer.

II. **To use hanging drops:** Allow 2 hanging drops of fingerstick whole blood (approximately 50 µL) to fall into the specimen well (S) of test cassette, then add 1 drop of buffer (approximately 40 µL) and start the timer.

1b) Venipuncture whole blood specimens:

Hold the dropper vertically and transfer 2 drops of whole blood (approximately 50 µL) to the specimen well (S), then add 1 drop of buffer (approximately 40 µL), and start the timer.

1c) Serum or plasma specimens:

Hold the dropper vertically and transfer 3 drops of serum or plasma (approximately 75 µL) to the specimen well (S) and start the timer.

- Wait for the coloured line(s) to develop. Read results at 10 minutes. Do not interpret the result after 30 minutes.

10. Result Interpretation

Positive:

Two lines develop. One coloured line should be in the control line region (C) and the other coloured line should be in the test line region (T).

Note: The intensity of the colour of the test line will vary depending on the concentration of *T. pallidum* antibodies present in the specimen. Therefore, any shade of colour of the test line should be considered positive.

Negative:

One coloured line develops in the control line region (C). No line develops in the test line region (T).

Invalid:

The absence of the control line (C) indicates a failure in the test procedure. The sample should be retested.

Note: Insufficient specimen volume, incorrect procedural techniques are the most likely reasons for control line failure. Review the procedure and repeat the test with a new test cassette. If the problem persists, discontinue using the test kit immediately and contact your distributor.

11. Quality Control

A procedural control is included in the test cassette. A coloured line developing in the control line region (C) is considered an internal procedural control. It confirms sufficient specimen volume, adequate membrane wicking and correct procedural technique.

12. Limitations

- The NADAL® Syphilis test is for in-vitro diagnostic use only.
- The test should only be used for the detection of *T. pallidum* antibodies in whole blood, serum or plasma specimens. Neither the quantitative value nor the rate of increase of *T. pallidum* antibodies can be determined by this qualitative test.
- The NADAL® Syphilis test only detects the presence of *T. pallidum* antibodies in the specimen and should not be used as the sole method for the diagnosis of *T. pallidum* infection.
- As with all diagnostic tests, all results should be interpreted together with other clinical information available to the physician.
- If the test result is negative and clinical symptoms persist, additional testing using other clinical methods is

recommended. A negative result does not at any time preclude the possibility of *T. pallidum* infection.

13. Performance Characteristics

Expected values

The NADAL® Syphilis test has been compared with a leading commercial TPHA Syphilis test, demonstrating an overall accuracy equal to or greater than 99.7%.

Clinical Sensitivity, Specificity and Accuracy

The NADAL® Syphilis test has correctly identified specimens of a seroconversion panel and has been compared to a leading commercial TPHA Syphilis test using clinical specimens. The results show that the relative sensitivity of the NADAL® Syphilis test is 99.7%, and the relative specificity is 99.6%.

The NADAL® Syphilis test vs. TPHA

		TPHA		
		Positive	Negative	Total
NADAL® Syphilis Test	Positive	384	2	386
	Negative	1	493	494
	Total	385	495	880

Relative sensitivity: 99.7% (98.6%-100.0%)*

Relative specificity: 99.6% (98.6%-100.0%)*

Accuracy: 99.7% (99.0%-99.9%)*

* 95% Confidence Intervals

Precision

Intra-Assay

Within-run precision has been determined by using 10 replicates of four specimens: a negative, a low positive, a medium positive and a high positive. The negative, low positive, medium positive and high positive values were correctly identified >99% of the time.

Inter-Assay

Between-run precision has been determined by 10 independent assays on the same four specimens: a negative, a low positive, a medium positive and a high positive. Three different lots of the NADAL® Syphilis test have been tested over a 3-day period using negative, low positive, medium positive and high positive specimens. The specimens were correctly identified >99% of the time.

14. References

1. Claire M. Fraser. Complete genome sequence of *Treponema Pallidum*, the Syphilis spirochete, *Science* 1998; 281 July: 375-381
2. Center for Disease Control. Recommendations for diagnosing and treating Syphilis in HIV-infected patients, *MMWR Morb. Mortal Wkly Rep.* 1988; 37: 601
3. Aral R. Marx. Crack, sex and STD, *Sexually Transmitted Diseases*, 1991; 18:92-101
4. J.N. Wasserheit. Epidemiological Synergy: Interrelationships between human immunodeficiency virus infection and other sexually transmitted diseases, *Sexually Transmitted Diseases* 1992; 19:61-77
5. Johnson Phillip C. Testing for Syphilis, *Dermatologic Clinic* 1994; 12 Jan: 9-17

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2. Center for Disease Control. Recommendations for diagnosing and treating syphilis in HIV-infected patients, MMWR Morb. Mortal Wkly Rep. 1988; 37: 601
3. Aral R. Marx. Crack, sex and STD, Sexually Transmitted Diseases, 1991; 18:92-101
4. J.N. Wasserheit. Epidemiological Synergy: Interrelationships between human immunodeficiency virus infection and other sexually transmitted diseases, Sexually Transmitted Diseases 1992; 19:61-77
5. Johnson Phillip C. Testing for Syphilis, Dermatologic Clinic 1994; 12 Jan: 9-17.

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14. Bibliografia

1. Claire M. Fraser. Complete genome sequence of *Treponema Pallidum*, the Syphilis spirochete, Science 1998; 281 July: 375-381
2. Center for Disease Control. Recommendations for diagnosing and treating Syphilis in HIV-infected patients, MMWR Morb. Mortal Wkly Rep. 1988; 37: 601
3. Aral R. Marx. Crack, sex and STD, Sexually Transmitted Diseases, 1991; 18:92-101
4. J.N. Wasserheit. Epidemiological Synergy: Interrelationships between human immunodeficiency virus infection and other sexually transmitted diseases, Sexually Transmitted Diseases 1992; 19:61-77
5. Johnson Phillip C. Testing for Syphilis, Dermatologic Clinic 1994; 12 Jan: 9-17

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14. Bibliografia

1. Claire M. Fraser. Complete genome sequence of *Treponema Pallidum*, the syphilis spirochete, Science 1998; 281 July: 375-381.
2. Center for Disease Control. Recommendations for diagnosing and treating syphilis in HIV-infected patients, MMWR Morb. Mortal Wkly Rep. 1988; 37: 601
3. Aral R. Marx. Crack, sex and STD, Sexually Transmitted Diseases, 1991; 18:92-101
4. J.N. Wasserheit. Epidemiological Synergy: Interrelationships between human immunodeficiency virus infection and other sexually transmitted diseases, Sexually Transmitted Diseases 1992; 19:61-77
5. Johnson Phillip C. Testing for Syphilis, Dermatologic Clinic 1994; 12 Jan: 9-17.

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nal von minden GmbH