NADAL® Chikungunya IgM Test (test cassette)

REF 682002N-20

DE	Gebrauchsanweisung	2	Symbols	16
EN	Instructions for use	5	Our Teams	16
ES	Instrucciones de uso	8		
IT	Istruzioni per l'uso	11		

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Relative Spezifität: 98,6%

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14. Referenzen

1. Shah KV, Gibbs CJ Jr, Banerjee G. Virological investigation of the epidemic of haemorrhagic fever in Calcutta: isolation of three strains of Chikungunya virus. *Indian J Med Res* 1964; 52:676-83.
2. Powers AM, Brault AC, Tesh RB, Weaver SC. Re-emergence of Chikungunya and O'nyong-nyong viruses: evidence for distinct geographical lineages and distant evolutionary relationships. *J Gen Virol* 2000;81:471-9
3. Myers RM and Carey DE. Concurrent isolation from patient of two arboviruses, Chikungunya and dengue type 2. *Science* 1967;157:1307-8.
4. Thein S, La Linn M, Aaskov J, Aung MM, Aye M, Zaw A, Myint A. Development of a simple indirect enzyme-linked immunosorbent assay for the detection of mmunoglobulin M antibody in serum from patients following an outbreak of chikungunya virus infection in Yangon, Myanmar. *Trans R Soc Trop Med Hyg.* 1992 86:438-42.

Rev. 0, 2015-11-16 OM/UJA

1. Intended Use

The NADAL® Chikungunya IgM Test is a lateral flow chromatographic immunoassay for the qualitative detection of anti-chikungunya virus IgM in human serum, plasma or whole blood. It is intended for use as a screening test and as an aid in the diagnosis of chikungunya infection. This test provides only a preliminary test result. Therefore, further reference tests, such as ELISA or PCR, should be used to confirm Chikungunya infection.

2. Introduction and Clinical Significance

The Chikungunya virus (CHIKV) is an enveloped, positive-strand, RNA virus belonging to the family *Togaviridae* together with the genus *Alphavirus*, first identified in 1953. Chikungunya (CHIK) fever is transmitted to humans by the bite of a variety of mosquitoes including *Ae. aegypti*, *Ae. africanus*, *Ae. luteocephalus*, *Ae. furcifer* and *Ae. Taylori*, but mainly by that of the Asian tiger mosquito (*Ae. albopictus*). CHIKV has caused outbreaks in East Africa (Tanzania and Uganda), in Austral Africa (Zimbabwe and South Africa), in West Africa (Senegal and Nigeria) and in Central Africa (Central African Republic and Democratic Republic of the Congo).

In Asia, CHIK outbreaks have been reported in India, Sri Lanka, Myanmar, Thailand, Indonesia, the Philippines, Cambodia, Vietnam, Hong Kong and Malaysia. The symptoms including sudden onset of fever, chills, headache, nausea, vomiting, joint pain with or without swelling, low back pain and rash are very similar to those of dengue. Both diseases are primarily transmitted by the same species of mosquitos - *Ae. aegypti* and *Ae. albopictus*. Mixed outbreaks of chikungunya with sporadic cases of dengue have been reported in Andhra Pradesh state, India. However, unlike dengue, there is no hemorrhagic or shock syndrome form. Therefore, the possibility of distinguishing between CHIKV infection and dengue virus infection is extremely beneficial, particularly in areas where dengue virus infection is endemic or epidemic. The NADAL® Chikungunya IgM Test utilises recombinant antigens which are found in structure proteins. It detects anti-CHIKV IgM in serum, plasma or whole blood samples after 15 minutes. The test can be performed by untrained or minimally skilled personnel and without cumbersome laboratory equipment.

3. Test Principle

The NADAL® Chikungunya IgM Test is a lateral flow chromatographic immunoassay. The test cassette consists of: 1) a burgundy conjugate pad containing CHIKV antigens conjugated to colloidal gold (CHIKV conjugates) and rabbit IgG-gold conjugates, 2) a nitrocellulose membrane strip containing a test line region (T) and a control line region (C). Anti-human IgM are pre-coated in the test line region (T), and goat anti-rabbit IgG are pre-coated in the control line region (C).

When an adequate volume of specimen is dispensed into the sample well (S) of the test cassette, it migrates by capillary action along the membrane. If IgM against CHIKV are present in the specimen, they bind to the CHIKV conjugates. The immunocomplex is then captured on the membrane by the pre-coated anti-human IgM, forming a burgundy test line (T), indicating a CHIKV IgM positive test result.

The absence of the test line (T) suggests a negative result.

The test contains an internal control (C line) which should exhibit a burgundy line of the immunocomplexes of goat anti-rabbit IgG/rabbit IgG-gold conjugate regardless of the the presence of the test line (T). Otherwise, the test result is invalid and the specimen should be retested with another test cassette.

4. Reagents and Materials Supplied

- 20 NADAL® Chikungunya IgM test cassettes (incl. disposable pipettes)
- 1 Buffer
- 1 package insert

5. Additional Materials Required

- Specimen collection containers
- Lancets (for fingerstick whole blood only)
- Centrifuge (for serum/plasma)
- Heparinised capillary tubes and dispensing bulb (for fingerstick whole blood only)

6. Storage & Stability

All reagents are ready to use as supplied. Store unused test cassettes unopened at 2-30°C. If test cassettes are stored at 2-8°C, ensure that they brought to room temperature before opening. Test cassettes are stable until the expiration date printed on the foil pouch. Do not freeze the test kit or expose it to temperatures over 30°C.

7. Warnings and Precautions

- For professional *in-vitro* diagnostic use only.
- Carefully read through the test procedure prior to testing.
- Do not use the test beyond the expiration date indicated on the package.
- Do not use the test if the foil pouch is damaged.
- Do not reuse tests.
- Do not add samples to the reaction area (result area).
- In order to avoid contamination, do not touch the reaction area (result area).
- Avoid cross-contamination of specimens by using a new specimen collection container for each specimen obtained.
- Do not substitute or mix components from different test kits.
- Do not eat, drink or smoke in the area where specimens and test kits are handled.
- Wear protective clothing such as laboratory coats, disposable gloves and eye protection when specimens are being assayed.
- Handle all specimens as if they contain infectious agents. Observe established precautions for microbiological risks throughout all procedures and standard guidelines for the appropriate disposal of specimens.
- The test kit contains products of animal origin. Certified knowledge of the origin and/or sanitary state of the animals does not completely guarantee the absence of transmissible pathogenic agents. It is therefore recommended that these products be treated as potentially infectious, and handled in accordance with usual safety precautions (e.g., do not ingest or inhale).
- Humidity and temperature can adversely affect test results.

- Used testing materials should be discarded according to local regulations.

8. Specimen Collection and Preparation

The NADAL® Chikungunya IgM Test can be performed using whole blood (from venipuncture or fingerstick), serum or plasma.

To collect fingerstick whole blood specimens:

- Wash the patient's hand with soap and warm water or clean with an alcohol swab. Allow it to dry.
- Massage the hand without touching the puncture site by rubbing down the hand towards the fingertip of the middle or ring finger.
- Puncture the skin with a new sterile lancet for each person. Wipe away the first drop of blood.
- Gently rub the hand from wrist to palm and then to finger to form a rounded drop of blood over the puncture site.

Separate serum or plasma from blood as soon as possible to avoid hemolysis. Use only clear, non-hemolysed specimens.

Testing should be performed immediately after specimen collection. Do not leave the specimens at room temperature for prolonged periods of time. Serum and plasma specimens can be stored at 2-8°C for up to 3 days. For long term storage, specimens should be kept below -20°C. Whole blood collected by venipuncture should be stored at 2-8°C if the test is to be run within 2 days of collection. Do not freeze whole blood specimens. Fingerstick whole blood should be tested immediately. Bring specimens to room temperature prior to testing. Frozen specimens must be completely thawed and mixed well prior to testing. Specimens should not be frozen and thawed repeatedly. If specimens are to be shipped, they should be packed in compliance with local regulations for the transportation of etiologic agents.

9. Test Procedure


Bring tests, specimens, buffer and/or controls to room temperature (15-30°C) prior to testing.

- Remove the test cassette from the sealed foil pouch and use it as soon as possible. Best results will be obtained if the assay is performed within one hour.

- Place the test cassette on a clean and level surface.

a) For serum or plasma specimen:

Holding the dropper vertically, add 1 drop (approximately 30 µL) of serum or plasma to the specimen well (S) of the test cassette.

b) For venipuncture whole blood specimen:

Holding the dropper vertically, add 2 drops (approximately 50 µL) of venipuncture whole blood to the specimen well (S) of the test cassette.

c) For fingerstick whole blood specimen:

I. To use a capillary tube:

Apply the end of the capillary tube to the blood until it is filled with approximately 50 µL. Avoid air

bubbles. Place the bulb onto the top end of the capillary tube, squeeze the bulb to dispense the whole blood specimen into the specimen well (S) of the test cassette.

II. To use hanging drops:

Position the patient's finger so that the drop of blood is exactly above the specimen well (S) of the test cassette. Allow 2 hanging drops (approximately 50 µL) of fingerstick whole blood to fall into the center of the specimen well (S) of the test cassette or move the patient's finger so that the hanging drop touches the center of the specimen well (S). Avoid touching the finger directly to the specimen well (S).

- Add 1 drop (approximately 40 µL) of buffer to the sample well (S). Start the timer.

- Wait for the coloured line(s) to appear. The result should be read after 15 minutes. Do not interpret the result after more than 15 minutes.

10. Result Interpretation

Positive


Two coloured lines appear on the membrane. One line appears in the control line region (C) and the other line appears in the test line region (T). The test result indicates the presence of anti-CHIKV IgM in the specimen.

Note: Samples with positive results should be confirmed with alternative testing method(s) and clinical findings before a positive determination is made.

Negative

One coloured line appears in the control line region (C). No coloured line appears in the test line region (T). The test result indicates that no detectable anti-CHIKV IgM are present in the specimen.

Invalid

The control line fails to appear. Results from any test which has not produced a control line at the specified reading time must be discarded. Please review the procedure and repeat the test with a new test cassette. If the problem persists, discontinue using the test kit immediately and contact your local distributor.

Insufficient specimen volume, incorrect operating procedure or expired tests are the most likely reasons for the control line failure.

11. Quality Control

An internal procedural control is included in the test cassette:

A coloured line appearing in the control line region (C) is considered an internal procedural control. It confirms sufficient specimen volume, adequate membrane wicking and correct procedural technique.

Good laboratory practice (GLP) recommends the use of control materials to ensure proper test kit performance.

12. Limitations

- The NADAL® Chikungunya IgM Test detects the presence of IgM antibodies to Chikungunya virus in the specimen and should not be used as the sole criterion for the diagnosis of Chikungunya.
- As with all diagnostic tests, a definitive clinical diagnosis should not be based on the result of a single test, but should only be made by the physician after all clinical and laboratory findings have been evaluated.
- If the test result is negative but clinical symptoms persist, additional follow-up testing using other clinical methods is recommended. Negative results do not preclude the possibility of Chikungunya infection.

13. Performance Characteristics

127 positive samples and 296 negative samples were tested using the NADAL® Chikungunya IgM Test and a commercially available CHIKV IgM capture ELISA.

The results are summarised in the following table.

		CHIKV IgM capture ELISA		
		+	-	Total
NADAL® Chikungunya IgM Test	+	123	4	127
	-	4	292	296
	Total	127	296	423

Relative sensitivity: 96.9%

Relative specificity: 98.6%

Overall agreement: 98.1%

14. References

1. Shah KV, Gibbs CJ Jr, Banerjee G. Virological investigation of the epidemic of haemorrhagic fever in Calcutta: isolation of three strains of Chikungunya virus. *Indian J Med Res* 1964; 52 :676-83.
2. Powers AM, Brault AC, Tesh RB, Weaver SC. Re-emergence of Chikungunya and O'nyong-nyong viruses: evidence for distinct geographical lineages and distant evolutionary relationships. *J Gen Virol* 2000;81:471-9
3. Myers RM and Carey DE. Concurrent isolation from patient of two arboviruses, Chikungunya and dengue type 2. *Science* 1967;157:1307-8.
4. Thein S, La Linn M, Aaskov J, Aung MM, Aye M, Zaw A, Myint A. Development of a simple indirect enzyme-linked immunosorbent assay for the detection of immunoglobulin M antibody in serum from patients following an outbreak of chikungunya virus infection in Yangon, Myanmar. *Trans R Soc Trop Med Hyg.* 1992 86:438-42.

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3. Myers RM and Carey DE. Concurrent isolation from patient of two arboviruses, Chikungunya and dengue type 2. Science 1967;157:1307-8.
4. Thein S, La Linn M, Aaskov J, Aung MM, Aye M, Zaw A, Myint A. Development of a simple indirect enzyme-linked immunosorbent assay for the detection of mungobulbin M antibody in serum from patients following an outbreak of chikungunya virus infection in Yangon, Myanmar. Trans R Soc Trop Med Hyg. 1992 86:438-42.

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Concordancia general: 98,1%

14. Referencias

1. Shah KV, Gibbs CJ Jr, Banerjee G. Virological investigation of the epidemic of haemorrhagic fever in Calcutta: isolation of three strains of Chikungunya virus. Indian J Med Res 1964; 52 :676-83.
2. Powers AM, Brault AC, Tesh RB, Weaver SC. Re-emergence of Chikungunya and O'nyong-nyong viruses: evidence for distinct geographical lineages and distant evolutionary relationships. J Gen Virol 2000;81:471-9

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Andamento complessivo: 98,1%

14. Bibliografia

1. Shah KV, Gibbs CJ Jr, Banerjee G. Virological investigation of the epidemic of haemorrhagic fever in Calcutta: isolation of three strains of Chikungunya virus. *Indian J Med Res* 1964; 52 :676-83.
2. Powers AM, Brault AC, Tesh RB, Weaver SC. Re-emergence of Chikungunya and O'nyong-nyong viruses: evidence for distinct geographical lineages and distant evolutionary relationships. *J Gen Virol* 2000;81:471-9
3. Myers RM and Carey DE. Concurrent isolation from patient of two arboviruses, Chikungunya and dengue type 2. *Science* 1967;157:1307-8.
4. Thein S, La Linn M, Aaskov J, Aung MM, Aye M, Zaw A, Myint A. Development of a simple indirect enzyme-linked immunosorbent assay for the detection of mungoblobulin M antibody in serum from patients following an outbreak of chikungunya virus infection in Yangon, Myanmar. *Trans R Soc Trop Med Hyg.* 1992 86:438-42.

Rev. 0, 2015-11-16 BN

Symbol	Deutsch	English	Français	Español	Italiano	Polski
	CE Konformitätszeichen	CE marking of conformity	Conformité aux normes européennes	Conformidad europea	Conformità europea	Znak zgodności CE
	Gebrauchsanweisung beachten	Consult instructions for use	Consulter la notice d'utilisation	Consúltense las instrucciones de uso	Consultare le istruzioni per l'uso	Przestrzegać instrukcji obsługi
	In-vitro-Diagnostika	In-vitro diagnostic medical device	Dispositif médical de diagnostic in vitro	Producto sanitario para diagnóstico in vitro	Dispositivo medico-diagnostico in vitro	Tylko do diagnostyki in vitro
	Temperaturbegrenzung	Temperature limitation	Limites de température	Límites de temperatura	Limiti di temperatura	Temperatura przechowywania
	Chargenbezeichnung	Batch code	Code du lot	Código de lote	Codice lotto	Numer serii
	Nicht zur Wiederverwendung	Do not reuse	Ne pas réutiliser	No reutilizar	Non riutilizzare	Tylko do jednorazowego użytku
	Verwendbar bis	Use by	Utiliser jusqu'au	Fecha de caducidad	Utilizzare entro	Data ważności
	Bestellnummer	Catalogue Number	Référence du catalogue	Número de catálogo	Riferimento di Catalogo	Numer katalogowy
	Hersteller	Manufacturer	Fabricant	Fabricante	Fabbricante	Producent
	Ausreichend für <n> Ansätze	Sufficient for <n> tests	Suffisant pour "n" tests	Suficiente para <n> utilizaciones	Sufficiente per "n" saggi	Wystarczający na <n> Powtórzeń

Our Teams

Germany:

Regensburg

Tel: +49 941 290 10-0
Fax: +49 941 290 10-50

Moers

Tel: +49 2841 99820-0
Fax: +49 2841 99820-1

Austria:

Tel: +49 941 290 10-29
Free Tel: 0800 291 565
Fax: +49 290 10-50
Free Fax: 0800 298 197

UK & Ireland:

Tel: +49 941 290 10-18
Free Tel –UK: 0808 234 1237
Free Tel –IRE: 1800 555 080
Fax: +49 290 10-50

France:

France Tel: 0800 915 240
France Fax: 0800 909 493

Switzerland

Swiss Tel: 0800 564 720
Swiss Fax: 0800 837 476

Belgium

Belgium Tel: 0800 718 82
Belgium Fax: 0800 747 07

Luxembourg

Lux. Tel: 800 211 16
Lux. Fax: 800 261 79

Spain:

Tel: +49 941 290 10-759
Free Tel: 900 938 315
Fax: +49 941 290 10-50
Free Fax: 900 984 992

Italy:

Tel: +49 941 290 10-34
Fax: +49 941 290 10-50

Poland:

Tel: +49 941 290 10-44
Free Tel: 00 800 491 15 95
Fax: +49 941 290 10-50
Free Fax: 00 800 491 15 94

Portugal:

Tel: +49 941 290 10-735
Tel. Verde: 800 849 230
Fax: +49 941 290 10-50
Fax Verde: 800 849 229

Netherlands:

Tel: +31 30 75 600
Free Tel: 0800 0222 890
Fax: +31 70 30 30 775
Free Fax: 0800 024 9519

Denmark:

Tel: +31 703075 603
Free Tel: +45 80 88 87 53
Tax: +31 703030 775

Laboratory Diagnostics Team:

Tel: +49 941 290 10-40
Fax: +49 941 290 10-50

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Carl-Zeiss-Strasse 12 • 47445 Moers • Germany

www.nal-vonminden.com • info@nal-vonminden.com

Fon: +49 2841 99820-0 • Fax: +49 2841 99820-1

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nal von minden GmbH