# NADAL<sup>®</sup> Entamoeba Test (test cassette)

REF 2090001

<b>DE</b> Gebrauchsanweisung	2	<b>PL</b> Sposób użycia	15
<b>EN</b> Instruction for use	5	Symbols	20
<b>FR</b> Instructions d'utilisation	7	Our Teams	20
<b>ES</b> Instrucciones de uso	10		
<b>IT</b> Istruzioni per l'uso	12		

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#### 15. Referenzen

1. REED S., et al. "Cloning of a Virulence Factor of *Entamoeba histolytica*". *Journal of Clinical Investigation*, Volume 91, April 1993, 1532-1540.
2. HAQUE R., et al. "Diagnosis of Amebic Liver Abscess and Intestinal Infection with the TechLab *Entamoeba histolytica* II Antigen Detection and Antibody Tests". *Journal of clinical microbiology*, Sept. 2000, p. 3235–3239.
3. ESPINOSA-CANTELLANO M., et al. "Pathogenesis of Intestinal Amebiasis: From Molecules to Disease". *Clinical Microbiology Reviews*, Apr. 2000, p. 318–331.

Rev. 0, 2014-08-14 OM

### 1. Intended Use

The NADAL® Entamoeba test is a rapid chromatographic immunoassay for the qualitative detection of *Entamoeba spp.* antigens in human faecal specimens to aid in the diagnosis of amoebiasis.

### 2. Introduction and Clinical Significances

Amoebiasis is the infection of human gastrointestinal tract by *Entamoeba histolytica*, a protozoan parasite which is capable of invading intestinal mucosa and may spread to other organs, mainly to liver. *Entamoeba dispar*, an amoeba morphologically similar to *E. histolytica* that also colonises in human gut, has been recently recognised as a separate species with no invasive potential. The acceptance of *E. dispar* as a distinct but closely related protozoan species has had profound implications for epidemiology of amoebiasis, since most asymptomatic infections found worldwide are now attributed to this non-invasive amoeba.

Invasive amoebiasis due to *E. histolytica* is more common in developing countries. In areas of endemic infection, a variety of conditions including ignorance, poverty, overcrowding, inadequate and contaminated water supplies and poor sanitation favour direct faecal-oral transmission of amoebas from one person to another. Being responsible for approximately 70 thousand deaths annually, amoebiasis is the fourth leading cause of death due to a protozoan infection after malaria, Chagas disease and leishmaniasis, and the third cause of morbidity in this organism group after malaria and trichomoniasis, according to the recent World Health Organisation estimate.

Patients with dysentery have an average of three to five mucosanguineous evacuations per day, with moderate colic pain preceding discharge, and rectal tenesmus. In patients with bloody diarrhoea evacuations are also few but stools are composed of liquid faecal material stained with blood. While there is moderate colic pain, there is no rectal tenesmus. Fever and systemic manifestations are generally absent.

### 3. Principle of the Test

The NADAL® Entamoeba test is a qualitative lateral flow immunoassay for the detection of *Entamoeba spp.* antigen in human faecal samples. Monoclonal and polyclonal antibodies against *Entamoeba* antigens are pre-coated in the test line region of the membrane. During testing, the sample reacts with particles coated with antibodies against *Entamoeba* which were pre-dried on the test strip. The mixture migrates along the membrane by capillary action. In the case of a positive result specific antibodies present on the membrane react with mixture conjugate and generate a red coloured line. A green coloured line should always appear in the control line region. It serves as verification that sufficient volume of sample has been added, that membrane wicking has occurred and as an internal control for reagents.

### 4. Reagents and Materials Supplied

- 10 NADAL® Entamoeba test cassettes
- 10 specimen collection tubes containing buffer
- 1 package insert

### 5. Additional Required Materials

- Specimen collection container

- Disposable gloves
- Timer

### 6. Storage & Stability

Store the test as packaged in the sealed foil pouch either refrigerated or at room temperature (2-30°C). The test is stable through the expiration date printed on the sealed foil pouch. The test must remain in the sealed foil pouch until use. Do not freeze.

### 7. Warnings and Precautions

- For professional *in-vitro* diagnostic use only.
- Do not use after expiration date.
- The test should remain in the sealed foil pouch until use.
- Do not use the test if the foil pouch is damaged.
- Wear protective clothing such as laboratory coats, disposable gloves and eye protection when specimens are being tested.
- Do not eat, drink or smoke in the area where the specimens or test kits are handled.
- Handle all specimens as if they contain infectious agents. Observe established precautions against microbiological hazards throughout testing and follow standard procedures for proper disposal of specimens.
- The test should be discarded in a proper biohazard container after testing.
- The test must be carried out within 2 hours of opening the sealed foil pouch.


### 8. Specimen Collection and Preparation

Collect sufficient quantity of feces (1-2 g or 1-2 mL for liquid sample). Stool samples should be collected in clean and dry containers (without preservatives or transport media). Samples can be stored in the refrigerator (2-4°C) for 1-2 days prior to testing. For longer storage samples must be kept frozen at -20°C. In this case, samples should be completely thawed and brought to room temperature before testing.

### 9. Procedure of the Test

#### To process the collected stool samples:


Use a separate specimen collection tube for each sample. Unscrew the tube cap and introduce the applicator stick three times into three different sites of faecal specimen to collect 125 mg of sample. Close the tube containing buffer and stool sample. Shake the tube in order to assure good sample dispersion. For liquid stool samples, aspirate faecal specimen with a dropper and add 125 µL of sample into the specimen collection tube containing buffer.

#### Test Procedure:

Bring tests, stool samples and buffer to room temperature (15-30°C) prior to testing. Do not open foil pouches until ready to perform the assay.

1. Remove the NADAL® Entamoeba test cassette from the sealed foil pouch and use it as soon as possible.
2. Shake the specimen collection tube to assure good sample dispersion. Break off the tip of the tube cap.
3. Use a separate test cassette for each sample. Dispense exactly 4 drops of sample into the specimen well (S) of the test cassette. Start the timer.
4. Read the result at 10 minutes after dispensing the sample into the specimen well (S).

### 10. Interpretation of the Results

#### Positive:

One green line develops in the control line region (C) and one red line develops in the test line region (T).

#### Note:

The intensity of the red line in the test line region (T) may vary depending on the concentration of antigens present in specimen. However, neither the quantitative value nor the rate of increase in antigens can be determined by this qualitative test.

#### Negative:

One green line develops in the control line region (C). No line develops in the test line region (T).

#### Invalid:

The absence of the green control line (C) regardless of the appearance or absence of the test line (T) indicates a failure in the test procedure. The sample should be retested.

**Note:** Insufficient specimen volume, incorrect procedural techniques or deterioration of reagents are the most likely reasons for control line failure. Review the procedure and repeat the test with a new test cassette. If the problem persists, discontinue using the test kit immediately and contact your distributor.

### 11. Quality Control

An internal procedural control is included in the test cassette:

A green line appearing in the control line region (C) is considered an internal procedural control. It confirms sufficient specimen volume, adequate membrane wicking and correct procedural technique.

### 12. Limitations

- The NADAL® Entamoeba test only indicates the presence of parasites in specimen and should be used for the qualitative detection of *Entamoeba* antigens in faecal specimens only.

Neither the quantitative value nor the rate of increase in antigen concentration can be determined by this test.

- An excess of sample can cause wrong results (brown lines appear). Dilute the sample with buffer and repeat the test.
- If the test result is negative and clinical symptoms persist, additional testing using other clinical methods is recommended. A negative result does not at any time preclude the possibility of amoebiasis.
- After one week of infection, the number of parasites in faeces is decreasing, making the sample less reactive. Stool samples should be collected within one week of the onset of symptoms.
- This test provides a presumptive diagnosis of infection caused by *Entamoeba*. All results must be interpreted in conjunction with other clinical information and laboratory findings available to the physician.

### 13. Expected Values

*Entamoeba histolytica* infects more than 500 million people worldwide. It has been observed almost from the time of its discovery that although *E. histolytica* most often causes mild or asymptomatic infections, 10% of patients develop severe dysentery and life-threatening invasive and extraintestinal disease.

100,000 people are estimated to die each year from amoebic colitis and amoebic liver abscess.

Dysenteric and diarrhetic syndromes account for 90% of cases of invasive intestinal amoebiasis.

### 14. Performance Characteristics

#### Sensitivity, Specificity

Some stool samples from patients of different hospitals were studied using the NADAL® Entamoeba test in comparison with another immunochromatographic test. The results show that the NADAL® Entamoeba test has sensitivity of >99% and specificity of >99%.

#### Cross-Reactivity

An evaluation was performed to determine cross-reactivity of NADAL® Entamoeba test. No cross-reactivity with common intestinal pathogens, other organisms and substances occasionally present in feces were observed with NADAL® Entamoeba test: *Escherichia coli*, *Campylobacter*, *Clostridium difficile*, *Giardia lamblia*, *H. pylori*, *Listeria monocytogenes*, *Salmonella*, *Shigella*.

### 15. References

1. REED S., et al. "Cloning of a Virulence Factor of *Entamoeba histolytica*". *Journal of Clinical Investigation*, Volume 91, April 1993, 1532-1540.
2. HAQUE R., et al. "Diagnosis of Amebic Liver Abscess and Intestinal Infection with the TechLab *Entamoeba histolytica* II Antigen Detection and Antibody Tests". *Journal of clinical microbiology*, Sept. 2000, p. 3235-3239.
3. ESPINOSA-CANTELLANO M., et al. "Pathogenesis of Intestinal Amebiasis: From Molecules to Disease". *Clinical Microbiology Reviews*, Apr. 2000, p. 318-331.

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### 15. Referencias

1. REED S, et al. "Cloning of a Virulence Factor of *Entamoeba histolytica*". *Journal of Clinical Investigation*, Volume 91, April 1993, 1532-1540.
2. HAQUE R, et al. "Diagnosis of Amebic Liver Abscess and Intestinal Infection with the TechLab *Entamoeba histolytica* II Antigen Detection and Antibody Tests". *Journal of clinical microbiology*, Sept. 2000, p. 3235-3239.
3. ESPINOSA-CANTELLANO M., et al. "Pathogenesis of Intestinal Amebiasis: From Molecules to Disease". *Clinical Microbiology Reviews*, Apr. 2000, p. 318-331.

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## 15. Bibliografia

1. REED S., et al. "Cloning of a Virulence Factor of Entamoeba histolytica". *Journal of Clinical Investigation*, Volume 91, April 1993, 1532-1540.
2. HAQUE R., et al. "Diagnosis of Amebic Liver Abscess and Intestinal Infection with the TechLab Entamoeba histolytica II Antigen Detection and Antibody Tests". *Journal of clinical microbiology*, Sept. 2000, p. 3235-3239.
3. ESPINOSA-CANTELLANO M., et al. "Pathogenesis of Intestinal Amebiasis: From Molecules to Disease". *Clinical Microbiology Reviews*, Apr. 2000, p. 318-331.

Rev. 0, 2014-08-14 MW

Symbol	Deutsch	English	Français	Español	Italiano	Polski
	CE Konformitätszeichen	CE marking of conformity	Conformité aux normes européennes	Conformidad europea	Conformità europea	Znak zgodności CE
	Gebrauchsanweisung beachten	Consult instructions for use	Consulter la notice d'utilisation	Consultense las instrucciones de uso	Consultare le istruzioni per l'uso	Przestrzegać instrukcji obsługi
	In-vitro-Diagnostika	In-vitro diagnostic medical device	Dispositif médical de diagnostic in vitro	Producto sanitario para diagnóstico in vitro	Dispositivo medico-diagnostico in vitro	Tylko do diagnostyki in vitro
	Temperaturbegrenzung	Temperature limitation	Limites de température	Limitación de temperatura	Limiti di temperatura	Temperatura przechowywania
	Chargenbezeichnung	Batch code	Code du lot	Código de lote	Codice lotto	Numer serii
	Nicht zur Wiederverwendung	Do not reuse	Ne pas réutiliser	No reutilizar	Non riutilizzare	Tylko do jednorazowego użytku
	Verwendbar bis	Use by	Utiliser jusqu'au	Fecha de caducidad	Utilizzare entro	Data ważności
	Bestellnummer	Catalogue Number	Référence du catalogue	Número de catálogo	Riferimento di Catalogo	Numer katalogowy
	Hersteller	Manufacturer	Fabricant	Fabricante	Fabbricante	Producent
	Ausreichend für <n> Ansätze	Sufficient for <n> tests	Suffisant pour pour "n" tests	Válido para para <n> ensayos	Sufficiente per "n" saggi	Wystarczający na <n> Powtórzeń

### Our Teams

#### Germany:

##### Regensburg

Tel: +49 941 290 10-0  
 Fax: +49 941 290 10-50

##### Moers

Tel: +49 2841 99820-0  
 Fax: +49 2841 99820-1

##### Austria:

Tel: +49 941 290 10-29  
 Free Tel: 0800 291 565  
 Fax: +49 290 10-50  
 Free Fax: 0800 298 197

##### UK & Ireland:

Tel: +49 941 290 10-18  
 Free Tel –UK: 0808 234 1237  
 Free Tel – IRE: 1800 555 080  
 Fax: +49 290 10-50

##### France:

France Tel: 0800 915 240  
 France Fax: 0800 909 493

#### Switzerland

Swiss Tel: 0800 564 720  
 Swiss Fax: 0800 837 476

#### Belgium

Belgium Tel: 0800 718 82  
 Belgium Fax: 0800 747 07

#### Luxembourg

Lux. Tel: 800 211 16  
 Lux. Fax: 800 261 79

#### Spain:

Tel: +49 941 290 10-759  
 Free Tel: 900 938 315  
 Fax: +49 941 290 10-50  
 Free Fax: 900 984 992

#### Italy:

Tel: +49 941 290 10-34  
 Fax: +49 941 290 10-50

#### Poland:

Tel: +49 941 290 10-44  
 Free Tel: 00 800 491 15 95  
 Fax: +49 941 290 10-50  
 Free Fax: 00 800 491 15 94

#### Portugal:

Tel: +49 941 290 10-735  
 Tel. Verde: 800 849 230  
 Fax: +49 941 290 10-50  
 Fax Verde: 800 849 229

#### Netherlands:

Tel: +31 30 75 600  
 Free Tel: 0800 0222 890  
 Fax: +31 70 30 30 775  
 Free Fax: 0800 024 9519

#### Denmark:

Tel: +31 703075 603  
 Free Tel: +45 80 88 87 53  
 Tax: +31 703030 775

#### Laboratory Diagnostics Team:

Tel: +49 941 290 10-40  
 Fax: +49 941 290 10-50

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Carl-Zeiss-Strasse 12 • 47445 Moers • Germany  
 www.nal-vonminden.com • info@nal-vonminden.com  
 Fon: +49 2841 99820-0 • Fax: +49 2841 99820-1

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nal von minden GmbH