NADAL® C. difficile GDH Ag Test (test cassette)

REF 582004

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15. Referenzen

1. Wren, M.W.D, et al. "Laboratory diagnosis of Clostridium difficile infection. An evaluation of tests for faecal toxin, glutamate dehydrogenase, lactoferrin and toxigenic culture in the diagnostic laboratory". British Journal of Biomedical Science, 66 (1), 2009
2. Vaishnavi, Ch., "Clinical spectrum & pathogenesis of Clostridium difficile associated diseases". Indian J. Med. Res. 131, April 2010, pp 487-499
3. POUTANEN, S. M. et al. "Clostridium difficile-associated diarrhoea in adults", CMAJ, 171(1) July 2004, pp. 51-58.

Rev. 2, 2015-11-10 OM/UJa

1. Intended Use

The NADAL® C. difficile GDH Ag Test is a rapid chromatographic immunoassay for the qualitative detection of *Clostridium difficile* glutamate dehydrogenase (*C. difficile* GDH) antigens in human faecal specimens. The test is intended for use as an aid in the diagnosis of *Clostridium difficile* infection and is designed for professional use only.

2. Introduction and Clinical Significance

Clostridium difficile (*C. difficile*) is an anaerobic gram-positive, spore-forming bacterium. The ability of *C. difficile* to form spores is the key feature which enables this bacterium to persist in patients and the physical environment for long periods and thereby facilitates its transmission. *C. difficile* is transmitted by the fecal-oral route.

C. difficile is the principal pathogen related to antibiotic associated diarrhea and/or pseudomembranous colitis in hospitalised patients. Mature colonic bacterial flora in a healthy adult is generally resistant to *C. difficile* colonisation. However, if the normal enteric flora is altered, resistance to colonisation is lost. Any factor associated with alteration of the normal enteric flora increases the risk of *C. difficile* colonisation after exposure to antibiotics, especially those with broad-spectrum activity such as penicillins, cephalosporins and clindamycin.

C. difficile can release two high-molecular-weight toxins: toxins A and B, which are responsible for the clinical manifestations ranging from mild, self-limited watery diarrhea to fulminant pseudomembranous colitis, toxic megacolon and death.

Clostridium difficile glutamate dehydrogenase (*C. difficile* GDH) is an enzyme produced in large quantities by all toxigenic and non-toxigenic strains, making it an excellent biomarker.

3. Test Principle

The NADAL® C. difficile GDH Ag Test is a qualitative lateral flow immunoassay for the detection of *C. difficile* GDH antigens in human faecal samples. Monoclonal antibodies against *C. difficile* GDH antigens are pre-coated in the test line region of the membrane. During testing, the sample reacts with anti-*C. difficile* GDH antibodies coated onto coloured particles which were pre-dried on the internal test strip. The mixture migrates along the membrane by capillary action. In the case of a positive result, specific antibodies present on the membrane react with the mixture conjugate and generate a coloured line.

A green line should always appear in the control line region (C) and serves as a procedural control, indicating that the proper volume of specimen has been added and membrane wicking has occurred.

4. Reagents and Materials Supplied

- 10 NADAL® C. difficile GDH Ag test cassettes
- 10 specimen collection tubes with buffer
- 1 package insert

5. Additional Materials Required

- Specimen collection container
- Disposable gloves
- Timer

6. Storage & Stability

Store the test as packaged in the sealed foil pouch at 2-30°C. The test is stable until the expiration date printed on the foil pouch. The test cassette should remain in the sealed foil pouch until use. Do not freeze tests.

7. Warnings and Precautions

- For professional *in-vitro* diagnostic use only.
- Carefully read through the test procedure prior to testing.
- Do not use the test beyond the expiration date indicated on the package.
- Do not use the test if the foil pouch is damaged.
- Do not reuse tests.
- Do not add samples to the reaction area (result area).
- In order to avoid contamination, do not touch the reaction area (result area).
- Avoid cross-contamination of specimens by using a new specimen collection tube for each specimen obtained.
- Do not substitute or mix components from different test kits.
- Do not eat, drink or smoke in the area where specimens and test kits are handled.
- Wear protective clothing such as laboratory coats, disposable gloves and eye protection when specimens are being assayed.
- Handle all specimens as if they contain infectious agents. Observe established precautions for microbiological risks throughout all procedures and standard guidelines for the appropriate disposal of specimens.
- The test kit contains products of animal origin. Certified knowledge of the origin and/or sanitary state of the animals does not completely guarantee the absence of transmissible pathogenic agents. It is therefore recommended that these products be treated as potentially infectious, and handled in accordance with usual safety precautions (e.g., do not ingest or inhale).
- Humidity and temperature can adversely affect test results.
- Used testing materials should be discarded according to local regulations.

8. Specimen Collection and Preparation

Specimen Collection and Storage:


A sufficient quantity of feces should be collected (min. 1-2 g or min. 1-2 mL for liquid sample). Stool samples should be collected without preservatives or transport media in clean and dry containers. Samples can be stored in the refrigerator (2-8°C) for up to 7 days prior to testing. For longer storage (maximum 1 year) the specimen must be kept frozen at -20°C/-4°F. In this case, the sample should be completely thawed and brought to room temperature before testing. Repeated freezing and thawing cycles are not recommended.

Specimen preparation:

Use a separate specimen collection tube with buffer for each sample. Unscrew the cap of the tube and insert the applicator stick into the faecal specimen four times to collect approx. 125 mg of the sample.



Close the tube with the buffer and stool sample. Shake the tube to ensure good sample dispersion. For liquid stool samples, collect the faecal specimen with a dropper and add 125 µL to the specimen collection tube with buffer.

9. Test Procedure

Bring the tests and samples with buffer to room temperature (15-30°C) prior to testing. Do not open foil pouches until you are ready to perform the assay.

1. Remove the NADAL® C. difficile GDH Ag test cassette from the sealed foil pouch and use it as soon as possible.
2. Shake the specimen collection tube to ensure good sample dispersion. Break off the tip of the tube.
3. Use a separate test cassette for each sample. Dispense exactly 4 drops into the specimen well (S) of the test cassette. Start the timer.
4. Read the result after exactly 10 minutes after dispensing the sample.

10. Result Interpretation

Positive:

Two lines appear in the result window. A green line appears in the control line region (C) and a red line appears in the test line region (T).

Negative:

Only one green line appears in the control line region (C).

Invalid:

If no green control line (C) develops, the assay is invalid, regardless of the presence or absence of the red test line.

Insufficient specimen volume, incorrect operating procedure or expired tests are the most likely reasons for the control line failure. In this case, review the procedure and repeat the test with a new test cassette. If the problem persists, discontinue using the test kit and contact your distributor.

Note: The colour intensity of the red line in the test line region (T) may vary depending on the concentration of *C. difficile* GDH antigens in the specimen. However, neither the quantitative value nor the rate of increase in *C. difficile* GDH antigens can be determined with this qualitative test.

11. Quality Control

An internal procedural control is included in the test cassette:

A green line appearing in the control line region (C) is considered an internal procedural control. It confirms sufficient specimen volume, adequate membrane wicking and correct procedural technique.

Good laboratory practice (GLP) recommends the use of control materials to ensure proper test kit performance.

12. Limitations

- The NADAL® C. difficile GDH Ag Test only indicates the presence of *C. difficile* GDH antigens in the specimen (qualitative detection) and should be used for the detection of *C. difficile* GDH antigens in faecal specimens only. Neither the quantitative value nor the rate of increase in *C. difficile* GDH antigen concentration can be determined with this test.
- An excess of sample may cause an inaccurate result (brown lines appear). Dilute the sample with buffer and repeat the test.
- Some stool samples can decrease the intensity of the control green line.
- The test provides a presumptive diagnosis of *C. difficile* infection. A confirmed diagnosis of infection should only be made by a physician after all clinical and laboratory findings have been evaluated in conjunction with the results of further clinical observations.
- Positive results confirm the presence of *C. difficile* GDH antigens in faecal samples. However, they may be also caused due to toxigenic or non-toxicogenic strains of *Clostridium difficile*. A positive result should be followed up with additional laboratory techniques to determine the strain.

13. Expected Values

Clostridium difficile is associated with 95-100% of cases of pseudomembranous colitis, 60-75% of cases of antibiotic-associated colitis and 35% of cases of antibiotic-associated diarrhea.

14. Performance Characteristics

Sensitivity and specificity

Stool samples from patients were tested using NADAL® C. difficile GDH Ag Test and two commercially available immunoassays (*C. DIFF QUIK CHEK COMPLETE*®, Techlab and Wampole® *C. Diff Chek*TM-60, Techlab). The following results were obtained:

Sensitivity >99% and specificity >99%
Sensitivity >95% and specificity >99%

Cross-reactivity

An evaluation was performed to determine the cross-reactivity of the NADAL® C. difficile GDH Ag Test. There was no cross-reactivity found with common gastrointestinal microorganisms occasionally present in faeces:

<i>Campylobacter</i>	<i>Listeria</i>	<i>Staphylococcus aureus</i>
<i>E. coli</i>	<i>Salmonella</i>	<i>Yersinia</i>
<i>H. pylori</i>	<i>Shigella</i>	

15. References

1. Wren, M.W.D, et al. "Laboratory diagnosis of Clostridium difficile infection. An evaluation of tests for faecal toxin, glutamate dehydrogenase, lactoferrin and toxigenic culture in the diagnostic laboratory". British Journal of Biomedical Science, 66 (1), 2009
2. Vaishnavi, Ch., "Clinical spectrum & pathogenesis of Clostridium difficile associated diseases". Indian J. Med. Res. 131, April 2010, pp 487-499
3. POUTANEN, S. M. et al. "Clostridium difficile-associated diarrhoea in adults", CMAJ, 171(1) July 2004, pp. 51-58.

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15. Referencias

1. Wren, M.W.D, et al. "Laboratory diagnosis of Clostridium difficile infection. An evaluation of tests for faecal toxin, glutamate dehydrogenase, lactoferrin and toxigenic culture in the diagnostic laboratory". British Journal of Biomedical Science, 66 (1), 2009
2. Vaishnavi, Ch., "Clinical spectrum & pathogenesis of Clostridium difficile associated diseases". Indican J. Med. Res. 131, April 2010, pp 487-499
3. POUTANEN, S. M. et al. "Clostridium difficile-associated diarrhoea in adults", CMAJ, 171(1) July 2004, pp. 51-58.

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nal von minden GmbH