# TOXOPLASMA ICT IgG-IgM

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## TOXOPLASMA ICT IgG-IgM

Immuno-chromatographic test  
for *in vitro* diagnostic use

# TOXO Ab ICT20 : 20 tests

### INSTRUCTION FOR USE

#### Intended use

TOXOPLASMA ICT IgG-IgM is a rapid test based on the immuno-chromatography technology (lateral flow), allowing the *simultaneous detection* of both IgG and IgM class anti-*Toxoplasma* antibodies in human sera.

#### Principle of the test

TOXOPLASMA ICT IgG-IgM is a unitary test for a qualitative diagnostic. It is based on the principle of the homogeneous sandwich (immunological reaction of two same antigen epitopes with the two binding sites of a bivalent antibody).

Inside the cassette, the device is composed of:

- a nitrocellulose strip on which are spread two reactive bands: the antigens (*Toxoplasma gondii*) of the “test” band (T band) and the rabbit gamma globulins of the “control” band (C band),
- a fiberglass support (conjugate pad) which is impregnated of red latex particles coupled with *Toxoplasma* antigens (“test” latex = T latex) and blue latex particles coupled with goat anti-rabbit IgG (“control” latex = C latex).

The test is run by successively dispensing the serum sample and an eluting solution (called the eluent) in the “sample well” of the cassette. Adding the eluent starts the concomitant migration (chromatography) of the serum and the latex particles. This migration is completed in 20-30 minutes.

If anti-*Toxoplasma* antibodies (IgG and/or IgM) are present in the sample, a complex is formed between the T latex and the patient’s antibodies which is then captured by the T band. It results in the appearance of a red line: the test is positive.

The direct capture of the C latex by the C band results in the appearance of a blue line, meaning that the chromatography performed well. The appearance of this blue line is systematic and independent of the serological status of the patient.

Both letters “T” and “C” are printed on the plastic frame of the cassette, materializing the position of the corresponding reading area.

## Kit components (20 TESTS)

- 20 ready to use tests: 2 bags (sealed + zip closure) of 10 devices including a desiccant packet,
- 1 x 1.4 ml of the eluent (dropper bottle),
- 1 instruction for use manual.

## Storage and stability

- **Store the original sealed bag between 2 and 8°C.** Cassettes can be used until the expiry date written on the bag label. Do not freeze. Do not use after the expiry date.
- **The first opening of a bag of 10 tests must occur at least 15 minutes after** putting the bag at room temperature in order to avoid condensation.
- **After the first opening** of a bag, keep it **at room temperature (18-30°C)**, carefully closed (zip closure), the desiccant packet inside. After opening, the cassettes can be used **for up to 2 month.**

## Precautions for use

### Safety

- For *in vitro* use only. Handle according to Good Laboratory Practices and consider any reagent and any sample as potentially toxic and/or infectious.
- For professional use only.
- All serum samples must be considered as potentially infectious and handled with care.
- Wear a lab coat, gloves and glasses; do not drink, eat or smoke in the laboratory. Do not mouth the pipettes.
- Dispose of waste (samples, tips, tubes, cassettes, used reagent...) according to good practices used in the industry and current regulations in the country.

### Precautions

- Do not use eluent from another lot number.
- Do not use cassettes from two different lot numbers in the same run.
- Close the vials after use; do not use if a substance was accidentally introduced in the reagents. Do not use reagent from a vial that presents signs of leakage. Do not use cloudy or precipitated solution.
- Use only disposable pipette tips. Avoid any inter-cassette contamination.
- Do not use reagents after their expiration date.
- The omission of a sample or the distribution of an inadequate volume may render the test result positive or negative, regardless of its actual serological status.
- Only use cassettes carefully stored in their closed bag, the desiccant packet inside.

## Serum collection and preparation

The test can be done with either serum or heparinised plasma.

Sample collection must be sterile and can be done either on dry tube or with heparin, citrate or EDTA.

On tubes with gel, do not collect gel: it could lead to false positives.

Avoid hemolysis as much as possible.

Keep the samples at 2-8 °C until they are processed. If they need to be stored, freeze the samples at  $-20 \pm 5$  °C. Do not use a contaminated sample. Avoid freezing and thawing the samples repeatedly.

## Test procedure

Additional material required: micropipette and disposable tips for dispensing volumes of 30µL, timer.

**If the bags of 10 tests are kept at 2-8°C, let them at least 15 minutes at room temperature** before opening: the temperature of the bag must reach the room temperature to avoid condensation.

1. Take out the desired number of cassettes, then close carefully the bag with the zip closure (the desiccant packet inside) while pressing out as much air as possible. **Store the bag at room temperature for up to 2 months.** The eluent can be stored at room temperature until expiration date.
2. **Identify each cassette** with the reference of each sample to be tested. Do not work with runs of more than 10 cassettes. Two successive runs of 10 cassettes must be separated with a few minutes in order to be able to make the reading at recommended times (use 2 timers).
3. Use a micropipette with a disposable tip to dispense 30µL of serum or plasma in the sample well. Do so for all cassettes before moving to the next step.
4. Dispense **3 drops** of the eluent of the kit. **Do not use the eluent of another lot number.** Keep the dropper vertically while dispensing. Close the dropper after use.
5. Start the timer when the eluent is dispensed in all the cassette of the run.

## Reading and interpretation

The reading must be done near a window or under direct light (i.e. a desk lamp). Avoid shadows on the reading area.

The reading must be done between 20 and 30 minutes after starting the timer.

***Do not take into account the results from readings after 30 minutes.***

- **Positive test:** 2 lines, a red “T” and a blue “C” appear in the corresponding areas. Every “T” line must be considered positive, even of very weak intensity. For very weak lines, make the reading with the eye vertically above the reading area.
- **Equivocal test:** The apparition of a grey shadow in the test area, most often between 20 and 30 minutes, is an equivocal result. The blue line must appear normally. The sample is considered negative. In this case, however, it is better to control the result with another technique or on a further sample.
- **Negative test:** No red line appears. Only the blue “C” line is visible.
- **Invalid test:** The “C” line does not appear. Read once again the instruction and repeat the test. If the problem persists, contact the manufacturer or your distributor.

**Notes:** This is a qualitative test. Intensity of the red line does not reflect the quantity of anti-*Toxoplasma* antibodies in the sample.

However, the rate of appearance of the red line seems correlated to the quantity of antibodies and could give an idea of the antibody titer.

Positivity of the test is proof of the contact of the patient with *Toxoplasma* but doesn't prejudge the contact date or the clinical status of the patient.

## Quality control

The blue “C” line allows the validation of the good running of the test. However it is recommended to incorporate from times to times a known weak positive sample in a run.

## Test limitations

- Do not use too old serum sample with this technique. It is recommended to use samples frozen for less than 2 years.
- The positivity can be caused by the presence of IgG and/or IgM anti *Toxoplasma*, the test being based on agglutination principle.
- The use of other body fluids (urine, CSF, saliva, whole blood...) has not been validated.
- The use of hemolytic, icteric or lipidic samples is not recommended. However, no interference in the reaction was recorded. Still, very hemolytic samples can hide a weak positive test, due to important red background from hemolysis.
- Do not dispense 2 or 4 drops of eluent.
- **Warning:** Forgetting to dispense the sample, and only use the three drops of eluent will lead to false positive results.
- Do not use eluent outside the one given with the cassettes (same lot number).
- **The Toxoplasma ICT IgG-IgM results must be interpreted in light of other clinical, serological, parasitological, epidemiological and medical imaging information in order to establish the diagnosis of toxoplasmosis.**

## Performances

### ❖ Sensitivity, Specificity:

#### 1. Material et methods:

The evaluation was done in a Reference Laboratory specialising in the diagnosis of toxoplasmosis and member of the French national reference centres network (CNR) for toxoplasmosis.

The serological status (positive or negative) was the final result of the routine diagnostic in the laboratory from the patient's biological and epidemiological data. The discordant biological results were solved by the use of confirmatory techniques LDBIO-TOXO II IgG, ISAGA IgM and Sabin and Feldman's DYE TEST).

The principle of the evaluation consisted of comparing the results of 486 samples (including 67 cord blood samples) obtained with the LDBIO ICT and the results obtained with 3 commercial serological tests (ELISA IgG, ELISA IgM and Indirect Hemagglutination Assay - IHA).

The results of ELISA IgG and ELISA IgM were compiled (ELISA G+M) in order to have a better comparison as both LDBIO ICT and IHA detect simultaneously IgG and IgM class antibodies.

Two studies were conducted in parallel:

- ✓ A prospective study on 356 samples from the routine activity of the laboratory.
- ✓ A retrospective study on 130 samples selected for their specific profile and spread into 5 groups:
  - Group 1 – seroconversion: retrospective analysis of 9 serum sequence (30 samples) from patients having presented a toxoplasmosis seroconversion during their pregnancy.
  - Group 2 – monitoring of uninfected children: This is a retrospective analysis of 15 samples corresponding to 5 sequences of post-natal monitoring of children born to mothers who presented a toxoplasmosis seroconversion during pregnancy.
  - Group 3 – false positive with ELISA: 4 false positive with ELISA IgG and 11 false positive ELISA IgM including 1 false positive and 6 equivocal results by ISAGA IgM.
  - Group 4 – particular populations :
    - residual IgM (n=5),
    - weak IgG level (n=20), selected in ELISA IgG with antibodies level between 1.1 and 3.3 UI/mL (note : of these 20 samples, 2 were therefore considered negative and 16 equivocal according to the criteria of the ELISA used in the study)
    - strong IgG level (n=10),
    - serological reactivation (n=10),
    - sera from multi-organ transplant donor (n=10, 9 positives, 1 negative).
  - Group 5 – potential cross reactions: sera negative for toxoplasmosis but with anti-CMV and/or anti-EBV (n=10) and sera with rheumatoid factor (RF, n=5).

## 2. Results:

- ✓ Whole population (prospective + retrospective studies): n=486

	ELISA G+M	HAI	LDBIO ICT		ELISA G+M	HAI	LDBIO ICT
Positive	151	175	189	Positive	16	2	7
Negative	9	14	0	Negative	276	295	290
Equivocal	29	0	0	Equivocal	5	0	0
Se	94.6%	92.9%	100%	Sp	94.6%	99.3%	97.7%

**Table 5:** results of the different techniques compared to the serological status (n=486) for the entire study. Equivocal results in ELISA were excluded from sensitivity and specificity calculations.

- ✓ Prospective study: n=356

	ELISA G+M	HAI	LDBIO ICT		ELISA G+M	HAI	LDBIO ICT
Positive	88	96	105	Positive	1	0	6
Negative	4	9	0	Negative	249	251	245
Equivocal	13	0	0	Equivocal	1	0	0
Se	95.7%	91.4%	100%	Sp	99.6%	100%	97.6%

**Table 6:** results of the different techniques compared to the serological status (n=486) for the prospective study. Equivocal results in ELISA were excluded from sensitivity and specificity calculations.

- ✓ Retrospective study n=130

- Group 1 - seroconversion:  
Of the 9 sequences, 1 is detected earlier by the LDBIO ICT. All other results were concordant for all three techniques.
- Group 2 – monitoring of uninfected children:  
LDBIO ICT allowed the detection of a serological positivity on a serum already negative with the two other techniques. All other results were concordant for all three techniques.
- Group 3 – false positive ELISA (NEGATIVE serological status) n=15:

	IHA		ISAGA IgM			LDBIO ICT	
	Pos	Neg	Pos	Eq	Neg	Pos	Neg
false positive ELISA IgG	0	4	NA	NA	NA	0	4
false positive ELISA IgM	2	9	1	6	4	0	11

**Table 7:** Compared performances of the different techniques for the false positive serum in ELISA

- Group 4 – particular populations:

With the exception of the weak IgG level samples, the results were concordant between all the techniques. No prozone effect was observed in any of the three techniques, even with 2000UI/mL IgG level.

For the weak IgG level (POSITIVE serological status), the results are as follow: n=20

	ELISA G+M			IHA		LDBIO ICT	
	Pos	Equivocal	Neg	Pos	Neg	Pos	Neg
IgG between 1,1 et 3,3 UI/ml	2	16	2	18	2	20	0

**Table 8:** Performances of the different techniques for the weak IgG level group

- Group 5 – potential cross reactions:

Of the 10 sera positive for EBV and/or CMV, none were found positive with ICT and IHA, but one was found equivocal with ELISA. Of the 5 sera positive for rheumatoid factor, one was falsy positive in ICT and equivocal in ELISA but correctly identified by IHA. All others samples were correctly identified.

There may be a potential cross reaction with this population, but a study on a greater number of serums is needed in order to conclude.

### 3. Conclusions :

Correlation between the serological status and LDBIO ICT is excellent:

- sensibility Se=100%,
- specificity Sp=97,7%,
- predictive negative value PNV = 100%.

LDBIO ICT has a low detection limit (<2UI/mL) without equivocal results.

The analysis of the 67 cord blood samples did not gave any discordant results between ICT and the conclusion, showing that LDBIO ICT could be use on such samples.

A positive result in LDBIO ICT cannot discriminate between IgM, IgG or the simultaneous presence of IgG and IgM.

#### ❖ Reproducibility:

Inter-series and inter-lot reproducibility were tested. In both cases, the serum to serum correlation is excellent.

#### ❖ Interferences:

Even though no particular cross-reaction has been observed with haemolysed, icteric or lipidic sera, it is recommended to interpret the results from the use of such samples with care.

## Bibliographie

- Franck, Jacqueline, Yves Jean-François Garin, et Henri Dumon. 2008. « LDBio-Toxo II immunoglobulin G Western blot confirmatory test for anti-toxoplasma antibody detection ». *Journal of clinical microbiology* 46 (7): 2334-38.  
doi:10.1128/JCM.00182-08.
- Jost, C, F Touafek, A Fekkar, R Courtin, M Ribeiro, D Mazier, et L Paris. 2011. « Utility of immunoblotting for early diagnosis of toxoplasmosis seroconversion in pregnant women ». *Clinical and vaccine immunology: CVI* 18 (11): 1908-12.  
doi:10.1128/CVI.05303-11.
- Khammari, Imen, Fatma Saghrouni, Sami Lakhall, Aida Bouratbine, Moncef Ben Said, et Jalel Boukadida. 2014. « A New IgG Immunoblot Kit for Diagnosis of Toxoplasmosis in Pregnant Women ». *The Korean Journal of Parasitology* 52 (5): 493-99.  
doi:10.3347/kjp.2014.52.5.493.
- Khammari, Imen, Fatma Saghrouni, Alia Yaacoub, Sondoss Gaied Meksi, Hinda Ach, Lamia Garma, Akila Fathallah, et Moncef Ben Saïd. 2013. « IgG Western Blot for Confirmatory Diagnosis of Equivocal Cases of Toxoplasmosis by EIA-IgG and Fluorescent Antibody Test ». *The Korean Journal of Parasitology* 51 (4): 485-88.  
doi:10.3347/kjp.2013.51.4.485.
- Leslé, F, F Touafek, A Fekkar, D Mazier, et L Paris. 2011. « Discrepancies between a new highly sensitive *Toxoplasma gondii* ELISA assay and other reagents: interest of Toxo IgG Western blot ». *European journal of clinical microbiology & infectious diseases: official publication of the European Society of Clinical Microbiology* 30 (10): 1207-12. doi:10.1007/s10096-011-1214-1.
- Maudry, A., G. Chene, R. Chatelain, B. Bellete, H. Patural, J. Hafid, H. Raberin, R. Tran Manh Sung, et P. Flori. 2009. « Expertise du nouveau test Access® TOXO-IgGII et comparaison avec trois autres techniques automatisées et la technique Western Blot LDBIO TOXO II IgG® ». *Immuno-analyse & Biologie Spécialisée* 24 (1): 42-49.  
doi:10.1016/j.immbio.2008.11.004.
- Maudry, Arnaud, Gautier Chene, Rémi Chatelain, Hugues Patural, Bahrie Bellete, Bernard Tisseur, Jamal Hafid, et al. 2009. « Bicentric evaluation of six anti-toxoplasma immunoglobulin G (IgG) automated immunoassays and comparison to the Toxo II IgG Western blot ». *Clinical and vaccine immunology: CVI* 16 (9): 1322-26.  
doi:10.1128/CVI.00128-09.
- Robert-Gangneux, F., et M.-L. Darde. 2012. « Epidemiology of and Diagnostic Strategies for Toxoplasmosis ». *Clinical Microbiology Reviews* 25 (2): 264-96.  
doi:10.1128/CMR.05013-11.

NF EN ISO 13485 - ISO 9001

19A rue Louis Loucheur - 69009 LYON - FRANCE  
 TEL : +33(0)4 7883 3487 - FAX : +33(0)4 7883 3430  
 www.ldbiodiagnostics.com - email : contact@ldbiodiagnostics.com